Manufactured By: HYPHEN BioMed

LIAPHEN Fibrinogen ON BCS

Latex Immunoassay for the quantitative
determination of Fibrinogen on BCS

1. LIAPHEN Fibrinogen (# A120102) reagents

	NAME	Reconstitution	Stability*	T° Stabilization
R1	Latex reagent	Ready to use	Refer insert	** 30 mn on board before any use

(*) Provided any contamination or evaporation is avoided. Stability can be adjusted according to the exact use conditions.

Stabilization of reagents: (**) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate board before any use.

Storage of reagents: Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C, and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

If the reagents are kept on the automate board, take care and use chimneys to limit as much as possible any evaporation of the reagents.

Homogenize the reagents before each use.

Any reagent of biological origin must be handled with all the required cautions, as being potentially infectious.

Do not interchange the reagents from different lots.

Reagents required but not provided:

- Plasma Calibrator titrated for Fibrinogen (eg: Biophen Plasma Calibrator #A222101)
- Normal and Abnormal quality control plasmas titrated for Fibrinogen (ex: BIOPHEN Normal Control Plasma -#A223201 and BIOPHEN Abnormal Control Plasma #A223301).
- Tris-NaCl-BSA pH 7.50 Buffer (eg: TBSA #AAR005A)

2. Preparation of the calibration curve and controls/samples.

- Calibration curve:

-Calibration is performed with a commercially available plasma calibrator, with a known Fibrinogen Concentration "C" in g/L (eg Biophen Plasma Calibrator #A222101). The calibrator must be loaded **1:10 manually prediluted in TBSA**. The test program for the calibration integrates dilutions of the calibrator (managed by the BCS) at various concentrations (to be defined by the user), that correspond approximately to: "C", 2.5, 2.0, 1.5, 1.0, 0.5, 0.3 g/L Fibrinogen.

- **Tested plasma and controls:**

For normal plasmas and controls, plasmas should be loaded **prediluted manually 1:10 in TBSA** and then assayed at the 1:20 dilution managed by the BCS (ie the test dilution is 1:200 final).

The control is performed with commercially available control plasmas, titrated for Fibrinogen. Various control plasmas are available: **Biophen Normal Control Plasma (#A223201)**, and **Biophen Abnormal Control Plasma (#A223301)** or **Abnormal Fibrinogen Control Low (#ASC070K)**.

Nota: For lyophilized calibrators and controls, following reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature. It is recommended to run the calibration curve with a freshly reconstituted calibrator. It is necessary to let the reagent temperature to stabilize for at least 30 minutes onto the automate before any use. Take care avoiding any contamination or evaporation of the reagents. Stability can be adjusted according to the exact use conditions.

Homogenize before each use.

Do not freeze calibrators and quality control plasmas.

Quality controls must be run regularly, and for each new batch of reagents, after an important maintenance of the instrument, or if measured values are not in compliance with the one expected for the method.

3. Results:

- The calibration curve (working range) is of the **Akima lin-lin type**.
- The values obtained for patients and controls are directly calculated from the calibration curve.
- The results are expressed as **g/l** of fibrinogen.

The calibration curve is validated when reactivity, as well as measured control values, are in compliance.

A new calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).

Performances may present variations according to the instrument used. Validate the expected values in the laboratory working conditions. Performances, as well as values for each new lot of quality controls used, must then be confirmed (and adjusted if necessary) in the laboratory working conditions.

NOTE:

- **The Calibrator, controls and samples are loaded, already manually prediluted 1:10 in TBSA, then a 1:20 (or adequate for calibrator) dilution is performed on BCS-XP.**

PROGRAMMATION OF THE ANALYZER:

Creation of reagents files and diluent:

Definition of test. Reagents without lot data. Create the reagent files (and diluent)

<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="2">Name: R1 LX FNG</td> </tr> <tr> <td>Abbreviation:</td> <td>R1 LX FNG</td> </tr> <tr> <td>clean when the reagent is modified</td> <td>Intensive</td> </tr> <tr> <td colspan="2"><input checked="" type="checkbox"/> New calibration with any new lot of reagents</td> </tr> <tr> <td colspan="2">Authorized lines</td> </tr> <tr> <td>from</td> <td>1 to 14</td> </tr> <tr> <td colspan="2">Stability (*)</td> </tr> <tr> <td>< 15C</td> <td>* > = 15C</td> </tr> <tr> <td colspan="2">Mix by jet</td> </tr> <tr> <td>Intensity</td> <td>1</td> </tr> <tr> <td>Frequency</td> <td>No agitation</td> </tr> <tr> <td colspan="2">Reference numbers</td> </tr> <tr> <td></td> <td>Number</td> </tr> <tr> <td></td> <td>add remove replace</td> </tr> </table>	Name: R1 LX FNG		Abbreviation:	R1 LX FNG	clean when the reagent is modified	Intensive	<input checked="" type="checkbox"/> New calibration with any new lot of reagents		Authorized lines		from	1 to 14	Stability (*)		< 15C	* > = 15C	Mix by jet		Intensity	1	Frequency	No agitation	Reference numbers			Number		add remove replace	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="2">Name: TBSA buffer</td> </tr> <tr> <td>Abbreviation:</td> <td>TBSA</td> </tr> <tr> <td>clean when the reagent is modified</td> <td>Normal/Low</td> </tr> <tr> <td colspan="2"><input checked="" type="checkbox"/> New calibration with any new lot of reagents</td> </tr> <tr> <td colspan="2">Authorized lines</td> </tr> <tr> <td>from</td> <td>1 to 14</td> </tr> <tr> <td colspan="2">Stability (*)</td> </tr> <tr> <td>< 15C</td> <td>* > = 15C</td> </tr> <tr> <td colspan="2">Mix by jet</td> </tr> <tr> <td>Intensity</td> <td>1</td> </tr> <tr> <td>Frequency</td> <td>No agitation</td> </tr> <tr> <td colspan="2">Reference numbers</td> </tr> <tr> <td></td> <td>Number</td> </tr> <tr> <td></td> <td>add remove replace</td> </tr> </table>	Name: TBSA buffer		Abbreviation:	TBSA	clean when the reagent is modified	Normal/Low	<input checked="" type="checkbox"/> New calibration with any new lot of reagents		Authorized lines		from	1 to 14	Stability (*)		< 15C	* > = 15C	Mix by jet		Intensity	1	Frequency	No agitation	Reference numbers			Number		add remove replace
Name: R1 LX FNG																																																									
Abbreviation:	R1 LX FNG																																																								
clean when the reagent is modified	Intensive																																																								
<input checked="" type="checkbox"/> New calibration with any new lot of reagents																																																									
Authorized lines																																																									
from	1 to 14																																																								
Stability (*)																																																									
< 15C	* > = 15C																																																								
Mix by jet																																																									
Intensity	1																																																								
Frequency	No agitation																																																								
Reference numbers																																																									
	Number																																																								
	add remove replace																																																								
Name: TBSA buffer																																																									
Abbreviation:	TBSA																																																								
clean when the reagent is modified	Normal/Low																																																								
<input checked="" type="checkbox"/> New calibration with any new lot of reagents																																																									
Authorized lines																																																									
from	1 to 14																																																								
Stability (*)																																																									
< 15C	* > = 15C																																																								
Mix by jet																																																									
Intensity	1																																																								
Frequency	No agitation																																																								
Reference numbers																																																									
	Number																																																								
	add remove replace																																																								

*To be defined by the user, according to the insert instructions and the exact use conditions.

Creation of calibration and controls files

Name: Cal	
Abbreviation:	Cal
clean when the reagent is modified	Low
<input checked="" type="checkbox"/> New calibration with any new lot of reagents	
Authorized lines	
from	1 to 14
Stability (*)	
< 15C	* > = 15C
Mix by jet	
Intensity	1
Frequency	No agitation
Reference numbers	
	Number
	add remove replace

DRAFT PROTOCOL / NOT VALIDATED

<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="2">Name: NormalCont</td> </tr> <tr> <td>Abbreviation:</td> <td>NormalCont</td> </tr> <tr> <td>clean when the reagent is modified</td> <td>Low</td> </tr> <tr> <td colspan="2"><input checked="" type="checkbox"/> New calibration with any new lot of reagents</td> </tr> <tr> <td colspan="2">Authorized lines</td> </tr> <tr> <td>from</td> <td>1 to 14</td> </tr> <tr> <td colspan="2">Stability (*)</td> </tr> <tr> <td>< 15C</td> <td>* > = 15C</td> </tr> <tr> <td colspan="2">Mix by jet</td> </tr> <tr> <td>Intensity</td> <td>1</td> </tr> <tr> <td>Frequency</td> <td>No agitation</td> </tr> <tr> <td colspan="2">Reference numbers</td> </tr> <tr> <td></td> <td>Number</td> </tr> <tr> <td></td> <td>add remove replace</td> </tr> </table>	Name: NormalCont		Abbreviation:	NormalCont	clean when the reagent is modified	Low	<input checked="" type="checkbox"/> New calibration with any new lot of reagents		Authorized lines		from	1 to 14	Stability (*)		< 15C	* > = 15C	Mix by jet		Intensity	1	Frequency	No agitation	Reference numbers			Number		add remove replace	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="2">Name: LowCont</td> </tr> <tr> <td>Abbreviation:</td> <td>LowCont</td> </tr> <tr> <td>clean when the reagent is modified</td> <td>Low</td> </tr> <tr> <td colspan="2"><input checked="" type="checkbox"/> New calibration with any new lot of reagents</td> </tr> <tr> <td colspan="2">Authorized lines</td> </tr> <tr> <td>from</td> <td>1 to 14</td> </tr> <tr> <td colspan="2">Stability (*)</td> </tr> <tr> <td>< 15C</td> <td>* > = 15C</td> </tr> <tr> <td colspan="2">Mix by jet</td> </tr> <tr> <td>Intensity</td> <td>1</td> </tr> <tr> <td>Frequency</td> <td>No agitation</td> </tr> <tr> <td colspan="2">Reference numbers</td> </tr> <tr> <td></td> <td>Number</td> </tr> <tr> <td></td> <td>add remove replace</td> </tr> </table>	Name: LowCont		Abbreviation:	LowCont	clean when the reagent is modified	Low	<input checked="" type="checkbox"/> New calibration with any new lot of reagents		Authorized lines		from	1 to 14	Stability (*)		< 15C	* > = 15C	Mix by jet		Intensity	1	Frequency	No agitation	Reference numbers			Number		add remove replace
Name: NormalCont																																																									
Abbreviation:	NormalCont																																																								
clean when the reagent is modified	Low																																																								
<input checked="" type="checkbox"/> New calibration with any new lot of reagents																																																									
Authorized lines																																																									
from	1 to 14																																																								
Stability (*)																																																									
< 15C	* > = 15C																																																								
Mix by jet																																																									
Intensity	1																																																								
Frequency	No agitation																																																								
Reference numbers																																																									
	Number																																																								
	add remove replace																																																								
Name: LowCont																																																									
Abbreviation:	LowCont																																																								
clean when the reagent is modified	Low																																																								
<input checked="" type="checkbox"/> New calibration with any new lot of reagents																																																									
Authorized lines																																																									
from	1 to 14																																																								
Stability (*)																																																									
< 15C	* > = 15C																																																								
Mix by jet																																																									
Intensity	1																																																								
Frequency	No agitation																																																								
Reference numbers																																																									
	Number																																																								
	add remove replace																																																								

(*) To be filled by the user according to the insert instructions and the exact working conditions. Homogenize before each use.

Creation of the procedure

Once the various reagents, calibration and controls files are created, the procedure file must be created. Go to the "Definition tests, procedure" and enter the following parameters.

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