**Adaptation of BIOPHEN Factor XIa on STAR**

**I. Reconstitution of the BIOPHEN Factor XIa reagents.**

Chromogenic determination of Factor XIa with Biophen Factor XIa.

**BIOPHEN Factor XIa (Ref A220412) reagents:**

	NAME	Reconstitution	Stability*	Stabilization in T°
R1A	F.X-VIIIc	3 ml of distilled water	24 hours at 2-8°C *** 8 hours at room T°	** 30 mn before any use
R1B	F.IX	3 ml of distilled water	24 hours at 2-8°C *** 8 hours at room T°	** 30 mn before any use
R2	Activator	3 ml of distilled water	24 hours at 2-8°C *** 8 hours at room T°	** 30 mn before any use
R3	SXa-11	3 ml of distilled water	1 months at 2-8°C *** 7 days at room T°	** 30 mn before any use
R4	Buffer TRIS-BSA	Ready to use	Stable until the expiration date at 2-8°C	** 30 mn before any use
Cal	Calibrator FXIa	2 ml of distilled water	24 hours at 2-8°C *** 8 hours at room T°	** 30 mn before any use

**\*The stability data claimed on the insert were obtained on reconstituted vials, kept closed, protected from and provided any contamination or evaporation is avoided. Stability must be controlled, and can be adjusted and validated if required, according to the exact use conditions for each laboratory.**

**Reconstitution:**

(\*) After reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature.

**Stabilization of reagents:**

(\*\*) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate before any use.

**Conservation of reagents:**

Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

**Foot-note: Do not interchange the reagents from different lots.**

**NOTE:** Use microcups when running the assay, to improve management of the reagents volumes by the automate.

## **II. Preparation of the calibration curve and controls/samples**

### **Calibration curve:**

Calibration is performed with the FXIa calibrator provided in the kit, with a known factor XIa concentration “C” (expressed in ng/ml).

The calibrator is loaded “undiluted”, and the calibration curve is directly realised by the automate.

**Controls:** The control is performed with internally or commercially available controls, titrated for Factor XIa.

### **Tested samples and controls:**

In these conditions, the tested samples and controls are loaded “undiluted”.

**Nota:** For lyophilized calibrators and controls, following reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature. It is recommended to run the calibration curve with a freshly reconstituted calibrator. It is necessary to let the reagent temperature to stabilize for at least 30 minutes onto the automate before any use. Take care avoiding any contamination or evaporation of the reagents. Stability can be adjusted according to the exact use conditions.

**Homogenize before each use.**

**Do not freeze calibrators and quality control plasmas.**

**Quality controls must be run regularly, and for each new batch of reagents, after an important maintenance of the instrument, or if measured values are not in compliance with the one expected for the method.**

## **III. Results:**

- The calibration curve (working range) is of the Lin (absorbance) – Lin (concentration) type.
- The values obtained for the samples and controls are directly calculated from the calibration curve.
- The results are expressed as ng/ml factor XIa activity.

**The calibration curve is validated when linearity ( $r^2$ ), as well as measured control values, are in compliance.**

**A new calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).**

**Note: Performances may present slight variations according to the instrument used. Validate the expected values in the laboratory working conditions. Performances, as well as values for each new lot of quality controls used, must then be confirmed (and adjusted if necessary) in the laboratory working conditions.**

#### IV. Programming the STAR analyser :

Click on the icon **set up software** for the **manager program** and create the program according to :

This calibration allows measuring concentrations in the range 0 -5ng/ml of Factor X1a

CONFIGURATION – Page 1 / 3

Identification									
Abbreviation		Name		Method			Date		
B F.X1a		BIOPHEN F.X1a		Two points Chromo					
Sample				Diluent					
Volume	Incu	Dil.	Id.	Name		Vial.	stab	Min.Vol.	
80µl	0 sec	1/1	R4 X1a	R4 X1a		25 ml	8 h	2 ml	
reagents									
Id.	Name	Incu sec	Vol. µl	Vial ml	Stab h	Min.Vol. ml	Washing		
							before	after	
Ra R1A X1a	R1A X1a	0*	80	3	8	0.2	no	Normal	
Rb R1B X1a	R1B X1a	120*	80	3	8	0.2	no	Intensive	
Rc R2 X1a	R2 X1a	120*	80	3	8	0.2	no	Special	
Rd R3 X1a	R3 X1a		80	3	8	0.2	no	Intensive	
Analysis		Result			Validation (ng/ml)				
Rd Heating	yes	Primary Units		ng/ml	Standard	min	max		
Stirring	yes	Corrector		1.00		0	10		
First point	20	Single/Duplicate		single					
Last point	200	Precision		5%					
		Redil. Condition							
*Precision : NO		<							
		>							

Enter data concerning washing.

Washing				
	Ra	Rb	Rc	Rd
Front	No	No	No	No
After	Normal	Intensive	Special	Intensive
Wash			12227	
Name			STA-DESORB U	
Stab. H			120	
Vial ml			15	
Vol in. ml			0.9	

Show the second page and write calibration data.

**CONFIGURATION – Page 2 / 3**

Calibration						
Mode graph – linear regression					Scale	
Calibrators				Vial* ml	Stab.* H	Min.Vol.* ml
Id.	Key	Name	Dil.			
XIaCAL		XIaCAL	1/1	2	24	0.2
XIaCAL		XIaCAL	1/2	2	24	0.2
XIaCAL		XIaCAL	1/4	2	24	0.2
XIaCAL		XIaCAL	1/8	2	24	0.2
XIaCAL		XIaCAL	1/16	2	24	0.2
<b>Control</b>						
	Id.	Key	Name	Vial. ml	Stab. H	Min.Vol. ml
<b>Level 1</b>						
<b>Level 2</b>						

(\*) To be filled by the user . Homogenize before each use

Show the last page and enter data concerning the quality control.

The STAR can manage 3 levels of quality control. For example: the following configuration is possible.

**CONFIGURATION – Page 3 / 3**

Printout / Transmission							
Units	Factor convers.	Print	Transmission Test number	Usual Values	ng/ml		
Main	ng/ml	Yes	0	STANDARD	Min. Max.		
Aux.1	DDO	Yes	0				
Aux.2				0	10		
Aux.3							
Printout limits		min	0				
		max	10				
<b>Controls</b>							
	Id.	Key	Name	Period H	Vial. ml	Stab. H	Min.Vol. ml
<b>Level 1</b>	(*)	(*)	(*)	(*)	(*)	(*)	(*)
<b>Level 2</b>	(*)	(*)	(*)	(*)	(*)	(*)	(*)
<b>Level 3</b>	(*)	(*)	(*)	(*)	(*)	(*)	(*)

(\*) To be filled by the user . Homogenize before each use

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