Manufactured By: HYPHEN BioMed

BIOPHEN HEPARIN ON COAGAMATE MTX

Determination of Anti-Xa activity for HBPM, HNF and Orgaran®

Adaptation of BIOPHEN HEPARIN ON COAGAMATE MTX

DRAFT NOT VALIDATED

1. Reconstitution of the BIOPHEN Heparin 3 reagent, (Reference A221003)

Chromogenic determination of the Anti-Xa activity of LMWH, UFH, and Orgaran.

	NAME	Reconstitution	Stability	Stabilization in T°
R1	SXa-11 Substrate	3 ml of distilled water (*)	3 months at 2-8°C * * 7 days at room T° Do not freeze	30 mn <i>on board</i> before any use(***)
R2	Factor Xa	3. ml of distilled water (*)	3 months at 2-8°C * * 7 days at room T° Do not freeze	30 mn <i>on board</i> before any use (***)
Diluent	Physiological Saline		24 hours on board	

2. Reconstitution of the BIOPHEN Heparin 6 reagent, (Reference A221006)

Chromogenic determination of the Anti-Xa activity of LMWH, UFH, and Orgaran.

	NAME	Reconstitution	Stability	Stabilization in T°
R1	SXa-11 Substrate	6 ml of distilled water (*)	3 months at 2-8°C * * 7 days at room T° Do not freeze	30 mn <i>on board</i> before any use(**)
R2	Factor Xa	6 ml of distilled water (*)	3 months at 2-8°C * * 7 days at room T° Do not freeze	** 30 mn <i>on board</i> before any use (**)

Recommendations for reconstitution:

Reconstitution (*): Reconstitute each vial with distilled water. Shake thoroughly (vortex). Let to homogenize for 30 minutes at room temperature (18-25 °C), while shaking the vial from time to time (vortex), until complete dissolution of the content. Check the absence of any solid at the bottom of the vial.

In current practice, in order to allow a good standardization, reconstitute these two reagents the evening before and put them at 2°-8°C following the 30 minutes at room temperature.

Note: In all cases, before use, check the absence of solids at the bottom of the vial, which confirms that dissolution is complete. If necessary, incubate for 1 hour at RT or better at 37°C, while shaking (vortex) from time to time. If required, then additionally incubate overnight at RT.

Storage of reagents: ()** Stability indicated are validated provided any evaporation or contamination of reconstituted reagents avoided. They must be adjusted according to the actual laboratory working conditions

Stability on STAR board under our laboratory working conditions and in full compliance with the recommendations indicated on this application, and provided that all cautions are strictly respected is: 5 DAYS.

Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C, and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

If the reagents are kept on the automate board 24H/24H, take care to limit as much as possible any evaporation of the reagents by using a chimney or mini-reducer.

Stabilization of reagents: (*)** It is necessary to let the substrate and the Factor Xa temperature to stabilize for at least 30 minutes on the automate before any use. A too low temperature of the reagents can induce an over-estimation. Conversely a too high temperature leads to an under evaluation of heparin.

Homogenize the reagents before each use.

Any reagent of biological origin must be handled with all the required cautions, as being potentially infectious.

Do not interchange the reagents from different lots.

If necessary use microcups when running the assay, to improve management of the reagents volumes by the automate

3. Determination of LMWH Heparins

The determination of Low molecular weight Heparin (LMWH) requires its own configuration of the MTX.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Heparin Calibrator (ref A222001)	1 ml of distilled water (*)	7 days at 2-8°C 48 hours at room T°	30 minutes on board before any use (**)
Quality controls Biophen LMWH Control (ref A223001)	1 ml of distilled water (*)	7 days at 2-8°C 48 hours at room T°	30 minutes on board before any use (**)

4. Determination of UFH Heparins

The determination of Unfractionated Heparin (UFH) requires its own configuration of the MTX This one differs from the configuration used for the *determination* of LMWH only by the denomination of calibrators and controls.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Heparin Calibrator (ref A222001)	1 ml of distilled water (*)	7 days at 2-8°C 48 hours at room T°	30 minutes on board before any use (**)
Quality controls Biophen UFH Control (ref A223101)	1 ml of distilled water (*)	7 days at 2-8°C 48 hours at room T°	30 minutes on board before any use (**)

5. Determination of Orgaran:

The determination of Orgaran requires its own configuration of the MTX. This one differs from the configuration used for the *determination* of UFH/LMWH only by the denomination of calibrators and controls.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Orgaran® Calibrator (ref A222201)	1 ml of distilled water (*)	7 days at room T° 48 hours at room T°	30 minutes on board of the instrument before any use (**)
Quality controls Biophen Orgaran® Control (ref A223501)	1 ml of distilled water (*)	7 days at room T° 48 hours at room T°	30 minutes on board of the instrument before any use (**)

Reconstitution:

(*) After reconstitution of calibrators or controls with distilled water, let them to stabilize for 30 minutes at room temperature. It is better to reconstitute calibrators the very day of calibration.

Conservation of reagents:

(**) Take care of strictly respecting the 30 minutes temperature stabilization time for *calibrators* and *controls* at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C. Homogenize before each use.

Nota: For lyophilized calibrators and controls, following reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature. It is recommended to run the calibration curve with a freshly reconstituted calibrator. It is necessary to let the reagent temperature to stabilize for at least 30 minutes onto the automate before any use. Take care avoiding any contamination or evaporation of the reagents. Stability can be adjusted according to the exact use conditions.

Homogenize before each use.

Do not freeze calibrators and quality control plasmas.

A new calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).

Performances may present slight variations according to the instrument used. Validate the expected values in the laboratory working conditions.

Check the validity of the series by including quality control plasmas at different levels in each one.

Quality controls must be run regularly, and for each new batch of reagents, after an important maintenance of the instrument, or if measured values are not in compliance with the one expected for the method.

6. Results:

- The calibration curve is of the Log (absorbance) - Lin (concentration) type. (*or alternatively lin-lin, to be verified*)
- The values obtained for the patients and controls are directly calculated from the calibration curve.
- The results are expressed in IU/ml.

- When Heparin or Orgaran concentrations are out of the working range, assayed plasma must be diluted in normal plasma, appropriately prepared and platelet poor, in order to keep a sufficient concentration of AT III.
- In presence of low AT III concentrations, as it can be the case in young children, an exogenous source of AT III is necessary, in order to correctly measure the heparin concentration.

Nota:

Biophen Heparin reagents are developed for measuring homogeneously Unfractionated heparin (UFH) and Low Molecular Weight Heparin (LMWH), using the same calibration curve. The superimposition conditions are susceptible to slightly vary according to the combination of reagents and lots used, and the technical characteristics and specific adjustment of each apparatus.

Therefore, this superimposition has to be verified and validated in the exact laboratory working conditions and for each machine, by establishing the calibration curve using Biophen Heparin calibrators (LMWH, ref A222001), and by assaying the homogeneous system of associated Biophen UFH calibrators (ref A222301), UFH control plasmas (ref A223101), and LMWH control plasmas at high and low levels (ref A223001, A223701).

The linearity (r^2) is expected to be ≥ 0.98 for the calibration curve.

Controls must be measured in compliance.

Measured values for calibrators are acceptable when:

Measured value = target value ± 0.10 IU/ml for levels ≤ 0.50 IU/ml (better ± 0.05 IU/ml)

Measured value = target value ± 0.15 IU/ml for levels > 0.5 IU/ml (better 90-110% of the target value, or ± 0.10 IU/ml for levels between 0.5 and 1.0 IU/ml)

Should the superimposition not be obtained, a separate calibration curve must be used for each type of heparin, and validated by assaying the homogeneous quality control system (controls measured in compliance and close to the assigned target value).

The calibration curve is validated when linearity (r^2), as well as measured control values, are in compliance.

A new calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).

Performances may present slight variations according to the instrument used. Validate the expected values in the laboratory working conditions. Performances, as well as values for each new lot of quality controls used, must then be confirmed (and adjusted if necessary) in the laboratory working conditions.

Nota: The 3 months stability at 2-8°C and the 7 days stability at room temperature (18-25°C) claimed are validated on reconstituted vials kept closed, protected from any evaporation or contamination.

For your information, stability data on STAR board tested in our laboratory working conditions are:

-7 days for reconstituted vials loaded on STAR during working time, using mini reducers or chimneys then closed and stored at 2-8 °C overnight.

-5 days for reconstituted vials left loaded continuously on STAR.

Stabilities of reagents on each automate must be controlled, adjusted and validated under each exact laboratory working conditions.

I. PROGRAMMING ANALYZER

From the menu, select appropriate password in order to access programming panels.

Click on the different points of the board and do things in order:

- 1) Reagent Definition
- 2) Analyzer Creation
- 3) Profile Creation
- 4) Reagent Panel
- 5) Calibration

1. Reagent Definition

Choose <<Definition new reagent>>

Click and fill the board for the three reagents: Dilution Buffer (physiological saline NaCl 9g/L), FXa and Substrate.

Definition of reagents:

FXa Biophen Heparin

(FXa BH)

Information	
Lot Number	<input type="text"/>
Stability	<input type="text" value="72 h"/>
<input type="text" value="Buffer"/>	

Washing	
<input type="text" value="Buffer"/>	Washing 2

Mixing	
<input type="text" value="No"/>	

Definition of reagents:

Sub Biophen Heparin

(Sub BH)

Information	
Lot Number	<input type="text"/>
Stability	72 h

Washing	
<input type="text" value="Buffer"/>	Rinse cycle time 1.5 sec

Mixing	
<input type="text" value="No"/>	

Definition of Reagents:

NaCl Buffer

(NaCl)

Information	
Lot Number	<input type="text"/>
Stability	24h

Mixing	
<input type="text" value="No"/>	

2. Analyze creation

Analyze Creation:

Biophen Heparin LMWH

Procedure	
Total Sample volume,	60 µl
Buffer	<input type="text" value="NaCl"/> Dil 1+6
Reagent	<input type="text" value="FXa BH"/> 100µl
	<input type="text"/>
Activation Time	60 sec
Starting Reagent.	<input type="text" value="SUB BH"/> 100 µl
Measured time	50 sec

Normal Values	
Mini	0.1
Maxi	2.0

Unit Chose

-	UI/ml
-	
-	
-	
-	
-	
-	
-	
-	

Determination ou

Max variance 10%

Click on Parameters and fill the following board:

Technical parameter	
<input type="text" value="Kinetic"/>	<input type="button" value="Configuration"/>
Resolution photometer	0
Connexion Hote	

Click on Configuration and fill the following

Parameter

Latence 30 Sec
Post phase 5 Sec

3. Profile Creation

Click on new profile.

Fill profile creation = Name of analysis

Click on add analysis → Pull-Dow menu

Choose the test

4. Reagent Panel : do Nothing

5. Calibration

Name	Lot number	Value
Cal 5	<input type="text"/>	<input type="text"/>
Cal 4	<input type="text"/>	<input type="text"/>
Cal 3	<input type="text"/>	<input type="text"/>
Cal 2	<input type="text"/>	<input type="text"/>
Cal 1	<input type="text"/>	<input type="text"/>

- Start by the standard N°5
- Finish by the standard N°1 (put 0.01 for 0 IU/ml)

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