Manufactured by Hyphen BioMed.

# BIOPHEN HEPARIN

## ON OLYMPUS AU400-AU640-AU2700

Determination of Anti-Xa activity for HBPM, HNF and Orgaran®

### Adaptation of BIOPHEN HEPARIN ON OLYMPUS

**DRAFT NOT VALIDATED**

#### 1. Reconstitution of the BIOPHEN Heparin 6 reagent, (Reference A221006)

Chromogenic determination of the Anti-Xa activity of LMWH, UFH, and Orgaran.

	NAME	Reconstitution	Stability	Stabilization in T°
R1	SXa-11 Substrate	6 ml of distilled water (* )	3 months at 2-8°C ** 7 days at room T° Do not freeze	30 mn <i>on board</i> before any use(***)
R2	Factor Xa	7.5 ml of distilled water (* )	3 months at 2-8°C * * 7 days at room T° Do not freeze	30 mn <i>on board</i> before any use (***)
Diluent	Physiological Saline		24 hours on board	

#### 2. Reconstitution of the BIOPHEN Heparin 3 reagent, (Reference A221003)

Chromogenic determination of the Anti-Xa activity of LMWH, UFH, and Orgaran.

	NAME	Reconstitution	Stability	Stabilization in T°
R1	SXa-11 Substrate	3 ml of distilled water (* )	3 months at 2-8°C ** 7 days at room T° Do not freeze	30 mn <i>on board</i> before any use(***)
R2	Factor Xa	3.75 ml of distilled water (* )	3 months at 2-8°C * * 7 days at room T° Do not freeze	30 mn <i>on board</i> before any use (***)
Diluent	Physiological Saline		24 hours on board	

#### Recommendations for reconstitution:

**Reconstitution (\*):** Reconstitute each vial with distilled water. Shake thoroughly (vortex). Let to homogenize for 30 minutes at room temperature (18-25 °C), while shaking the vial from time to time (vortex), until complete dissolution of the content. Check the absence of any solid at the bottom of the vial.

In current practice, in order to allow a good standardization, reconstitute these two reagents the evening before and put them at 2°-8°C following the 30 minutes at room temperature.

**Note:** In all cases, before use, check the absence of solids at the bottom of the vial, which confirms that dissolution is complete. If necessary, incubate for 1 hour at RT or

**better at 37°C, while shaking (vortex) from time to time. If required, then additionally incubate overnight at RT.**

**Storage of reagents: (\*\*)** Stability indicated are validated provided any evaporation or contamination of reconstituted reagents avoided. They must be adjusted according to the actual laboratory working conditions

Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C, and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

If the reagents are kept on the automate board 24H/24H, take care to limit as much as possible any evaporation of the reagents by using a chimney or mini-reducer.

**-Chose of recipients on Olympus instrument: In order to have a minimum breakeven quantity. Pour the reagents into a long tube with a small diameter. Introduce this tube into a plastic recipient specific to the instrument. Cut the vial neck to introduce easily the tube (the tube must be sufficient longer for attending the bottom of the plastic vial)**

**Stabilization of reagents: (\*\*\*)** It is necessary to let the substrate and the Factor Xa temperature to stabilize for at least 30 minutes on the automate before any use. A too low temperature of the reagents can induce an over-estimation. Conversely a too high temperature leads to an under evaluation of heparin

**Homogenize the reagents before each use.**

**Any reagent of biological origin must be handled with all the required cautions, as being potentially infectious.**

**Do not interchange the reagents from different lots.**

**If necessary use microcups when running the assay, to improve management of the reagents volumes by the automate**

### **3. Determination of LMWH. Heparins**

The determination of Low molecular weight Heparin (LMWH) requires its own configuration of the Olympus.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Heparin Calibrator (ref A222001)	1 ml of distilled water (*)	7 days at 2-8°C 48 hours at room T°	30 minutes <i>on board</i> before any use (**)
Quality controls Biophen LMWH Control (ref A223001)	1 ml of distilled water (*)	7 days at 2-8°C 48 hours at room T°	30 minutes <i>on board</i> before any use (**)

### **4. Determination of UFH. Heparins**

The determination of Unfractionated Heparin (UFH) requires its own configuration of the Olympus. This one differs from the configuration used for the *determination* of LMWH only by the denomination of calibrators and controls.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Heparin Calibrator (ref A222001)	1 ml of distilled water (*)	7 days at 2-8°C 48 hours at room T°	30 minutes on board before any use (**)
Quality controls Biophen UFH Control (ref A223101)	1 ml of distilled water (*)	7 days at 2-8°C 48 hours at room T°	30 minutes on board before any use (**)

### 5. Determination of Orgaran:

The determination of Orgaran requires its own configuration of the Olympus. This one differs from the configuration used for the *determination* of UFH/LMWH only by the denomination of calibrators and controls.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Orgaran® Calibrator (ref A222201)	1 ml of distilled water (*)	7 days at room T° 48 hours at room T°	30 minutes on board before any use (**)
Quality controls Biophen Orgaran® Control (ref A223501)	1 ml of distilled water (*)	7 days at room T° 48 hours at room T°	30 minutes on board before any use (**)

#### Reconstitution:

(\*) After reconstitution of calibrators or controls with distilled water, let them to stabilize for 30 minutes at room temperature. It is better to reconstitute calibrators the very day of calibration.

#### Conservation of reagents:

(\*\*) Take care of strictly respecting the 30 minutes temperature stabilization time for *calibrators* and *controls* at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C. Homogenize before each use.

**Nota:** For lyophilized calibrators and controls, following reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature. It is recommended to run the calibration curve with a freshly reconstituted calibrator. It is necessary to let the reagent temperature to stabilize for at least 30 minutes onto the automate before any use. Take care avoiding any contamination or evaporation of the reagents. Stability can be adjusted according to the exact use conditions.

**Homogeneize before each use.**

**Do not freeze calibrators and quality control plasmas.**

**A new calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).**

**Performances may present slight variations according to the instrument used. Validate the expected values in the laboratory working conditions.**

**Check the validity of the series by including quality control plasmas at different levels in each one.**

**Quality controls must be run regularly, and for each new batch of reagents, after an important maintenance of the instrument, or if measured values are not in compliance with the one expected for the method.**

## 6. Results:

- The values obtained for the patients and controls are directly calculated from the calibration curve.
- The results are expressed in IU/ml.
- When Heparin or Orgaran concentrations are out of the working range, assayed plasma must be diluted in normal plasma, appropriately prepared and platelet poor, in order to keep a sufficient concentration of AT III.
- In presence of low AT III concentrations, as it can be the case in young children, an exogenous source of AT III is necessary, in order to correctly measure the heparin concentration.

**Nota:** Biophen Heparin reagents are developed for measuring homogeneously Unfractionated heparin (UFH) and Low Molecular Weight Heparin (LMWH), using the same calibration curve.

The superimposition conditions are susceptible to slightly vary according to the combination of reagents and lots used, and the technical characteristics and specific adjustment of each apparatus.

Therefore, this superimposition has to be verified and validated in the exact laboratory working conditions and for each machine, by establishing the calibration curve using Biophen Heparin calibrators (LMWH, ref A222001), and by assaying the homogeneous system of associated Biophen UFH calibrators (ref A222301) , UFH control plasmas (ref A223101), and LMWH control plasmas at high and low levels (ref A223001, A223701).

The linearity ( $r^2$ ) is expected to be  $\geq 0.98$  for the calibration curve.

Controls must be measured in compliance.

Measured values for calibrators are acceptable when:

Measured value = target value  $\pm 0.10$  IU/ml for levels  $\leq 0.50$  IU/ml (better  $\pm 0.05$  IU/ml)

Measured value = target value  $\pm 0.15$  IU/ml for levels  $> 0.5$  IU/ml (better 90-110% of the target value, or  $\pm 0.10$  IU/ml for levels between 0.5 and 1.0 IU/ml)

Should the superimposition not be obtained, a separate calibration curve must be used for each type of heparin, and validated by assaying the homogeneous quality control system (controls measured in compliance and close to the assigned target value).

**The calibration curve is validated when linearity ( $r^2$ ), as well as measured control values, are in compliance.**

**A new calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).**

**Performances may present slight variations according to the instrument used. Validate the expected values in the laboratory working conditions. Performances, as well as values for each new lot of quality controls used, must then be confirmed (and adjusted if necessary) in the laboratory working conditions.**

**Nota:** The 3 months stability at 2-8°C and the 7 days stability at room temperature (18-25°C) claimed are validated on reconstituted vials kept closed, protected from any evaporation or contamination.

For your information, stability data on STAR board tested in our laboratory working conditions are:

-7 days for reconstituted vials loaded on STAR during working time, using mini reducers or chimneys then closed and stored at 2-8 °C overnight.

-5 days for reconstituted vials left loaded continuously on STAR.

Stabilities of reagents on each automate must be controlled, adjusted and validated under each exact laboratory working conditions.

## **7. Programmation of the OLYMPUS instrument**

- Test name

Selected: (P) Parameter, (Y) Common test parameter, (A) name of the test.

Selected the tab: "Test Name".

Then selected the function ""Set" and choose a free position.

Set the name of the test.

Then selected: (P) Parameter, (Y) Common test parameter, (C) Tray.

Then select the functions "Set" then "Modified".

Click directly on the test to be added for introducing it into the tray.

- Technical sheet programming

Selected: (P) Parameter, (I) Specific test parameters, Tab "General"

Selected the name of the test in the heading "Test name" using the menu.

Then follow the programming:

Specific Test parameters									
General		LIH	ISE	range					
Test Name	HNF	▼	◀	▶	Type	Serum	▼		
Sample :	Volume	25	µl	Dilution	0	µl	Predilution Rate	1	
	R1 Volume	90	µl	Dilution	0	µl	Min OD		Max OD
	R2 Volume	120	µl	Dilution	0	µl		L -2.0000	H 2.5000
Wavelength	Pri.	410	▼	Sec.	Non	▼	Reagent Limit :		
Method	Fixe1		▼				First L	-2.0000	First H 2.5000
Reaction slope		+	▼				Last L	-2.0000	Last H 2.5000
Measuring Point 1	First	11		Last	13		L	0.00	H 1.50
Measuring Point 2	First			Last			Correlation Factor		
Linearity			%				A	1.00000	B 0.0000
No lag Time			▼				On-Board stability		

#### Programming the normal values and the Unit

In the tab "Range" [(P) Parameter, (I) Specific test parameters] program:

- The Normal values                   L: 0.00  
  H: 0.10
- The Unit used                         UI/ml
- The number of decimal             2

#### Programming the rerun

→ In the menu (P) Parameter, (R) repeat parameter, (C) Repeat specific. Do not enter dilution factor (Predilution rate 1)

#### Programming the calibration

→ Programming the number and the position of the calibrators  
Select: (P) Parameters, (B) Calibration, and then (C) Calibrator.

Chose the function "Set" then selected a free position for declaring the calibrator on the yellow support.

BIOPHEN HEPARIN on OLYMPUS

Page 6 sur 8

D.750.07/BI/1006-1003/OLY

Written : 29/07/2009 Version 3

[LLL](http://www.aniara.com)

8580 Gove Court • Mason, OH 45040

Phone: 513.770.1991

Toll Free: 866.783.3797

Fax: 513.573.9241

Email: info@aniara.com

[www.aniara.com](http://www.aniara.com)

→ Programming the type of calibration and the calibrator's values.  
 Selected: (P) Parameter, (B) Calibration, and then (C) Specific calibration.  
 Selected the name of the test in the heading "Test name" using the menu.  
 Selected the function "Set" and follow the programming:

Calibration specific					
General		ISE			
Test Name	HNF	◀	▶	Type :	Serum
Calibration Type	5AB	Formula	Polygonal	Counts	2
		Process	Conc		
	Cal. No.	OD	CONC	Factor/OD-L	Factor/OD-H
Point 1 :	21		0.00*	- 2.0000	2.5000
Point 2 :	22		0.38	- 2.0000	2.5000
Point 3 :	23		0.77	- 2.0000	2.5000
Point 4 :	24		1.14	- 2.0000	2.5000
Point 5 :	25		1.50	- 2.0000	2.5000
Point 6 :					
Point 7 :					
1-Point Cal. Point		<input type="checkbox"/>	With CONC-0	Advanced calibration	Non
MB Type Factor				Calibration Stability Period	

(\*) calibrators position and concentration given as example

Programming reagents

Fix the position of each reagent into the tray:

- Substrat on the tray 1
- Factor Xa on the tray 2

Click on the instrument's icon, Reagent STATUS.

Selected the right position by clicking directly on the reagent support.

Selected the function "Edit"  
 BIOPHEN HEPARIN on OLYMPUS  
 Written : 29/07/2009 Version 3

[LLL](#)

Chose the option "Fixed reagent"  
 Selected the test for each reagents R1 and R2  
 Precise the packaging (Bottle size)

**Calibration curve**

Test name: 71 HP    Type: Serum    Measured Date: 1.07/10/2005 14h50

1	21	0.8684	0.00	R1 :	
2	22	0.5928	0.38	R2 :	
3	23	0.4155	0.77	N° vial	
4	24	0.3046	1.14	R1 :	
5	25	0.2190	1.50	R2 :	
				Pos. R1 :	14
				Pos R2 :	12

Cal Type I: 5AB    Measured Type: Support

Formula: POLYGONAL

Cal Stability: -----    N° seq.:    -    1

Selected in the list printed with the spacebar or clicking the mouse

Help	Exit	Print	Comments	Re-Cal	Factor	Modified axes	N° Test
------	------	-------	----------	--------	--------	---------------	---------

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Dans le menu (P) Paramètre, ® Repassage, © Repassage spécifique ne pas saisir de facteur de prédilution (Zone de Prédilution 1).