Manufactured By: HYPHEN BioMed

# BIOPHEN HEPARIN ON BCS

VARIANT

## Adaptation of BIOPHEN HEPARIN ON BCS

### 1. Reconstitution of the BIOPHEN Heparin 6 reagent, (Reference A221006)

Chromogenic determination of the Anti-Xa activity of LMWH, UFH, and Orgaran.

	NAME	Reconstitution	Stability	Stabilization in T°
R1	SXa-11 Substrate	6 ml of distilled water (* )	3 months at 2-8°C ** 7 days at room T° Do not freeze	30 mn <i>on board</i> before any use(***)
R2	Factor Xa	7.5 ml of distilled water (* )	3 months at 2-8°C ** 7 days at room T° Do not freeze	30 mn <i>on board</i> before any use (***)
Diluent	Physiological Saline		24 hours on board **	

### 2. Reconstitution of the BIOPHEN Heparin 3 reagent, (Reference A221003)

Chromogenic determination of the Anti-Xa activity of LMWH, UFH, and Orgaran.

	NAME	Reconstitution	Stability	Stabilization in T°
R1	SXa-11 Substrate	3 ml of distilled water (* )	3 months at 2-8°C ** 7 days at room T° Do not freeze	30 mn <i>on board</i> before any use(***)
R2	Factor Xa	3.75 ml of distilled water (* )	3 months at 2-8°C ** 7 days at room T° Do not freeze	30 mn <i>on board</i> before any use (***)
Diluent	Physiological Saline		24 hours on board **	

#### Recommendations for reconstitution:

**Reconstitution (\*):** Reconstitute each vial with distilled water. Shake thoroughly (vortex). Let to homogenize for 30 minutes at room temperature (18-25 °C), while shaking the vial from time to time (vortex), until complete dissolution of the content. Check the absence of any solid at the bottom of the vial.

In current practice, in order to allow a good standardization, reconstitute these two reagents the evening before and put them at 2°-8°C following the 30 minutes at room temperature.

Note: In all cases, before use, check the absence of solids at the bottom of the vial, which confirms that dissolution is complete. If necessary, incubate for 1 hour at RT or better at 37°C, while shaking (vortex) from time to time. If required, then additionally incubate overnight at RT.

**Storage of reagents: (\*\*)** Stability indicated are validated provided any evaporation or contamination of reconstituted reagents avoided. They must be adjusted according to the actual laboratory working conditions

Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C, and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

If the reagents are kept on the automate board 24H/24H, take care to limit as much as possible any evaporation of the reagents by using a chimney or mini-reducer.

**Stabilization of reagents:** (\*\*\*) It is necessary to let the substrate and the Factor Xa temperature to stabilize for at least 30 minutes on the automate before any use. A too low temperature of the reagents can induce an over-estimation. Conversely a too high temperature leads to an under evaluation of heparin.

**Homogenize the reagents before each use.**

**Any reagent of biological origin must be handled with all the required cautions, as being potentially infectious.**

**Do not interchange the reagents from different lots.**

**If necessary use microcups when running the assay, to improve management of the reagents volumes by the automate**

### 3. Determination of Heparins UFH and/or LMWH.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Heparin Calibrator (ref A222001)	1 ml of distilled water (*)	7 days at 2-8°C 48 hours at room T°	30 minutes <i>on the automate board</i> before any use (**)
Quality controls Biophen UFH Control (ref A223101) Biophen LMWH Control (ref A223001)	1 ml of distilled water (*)	7 days at 2-8°C 48 hours at room T°	30 minutes <i>on the automate board</i> before any use (**)

### 4. Determination of Organan:

The determination of Organan requires its own configuration of the BCS. This one differs from the configuration used for the *determination* of UFH/LMWH only by the denomination of calibrators and controls.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Organan ® Calibrator (ref A222201)	1 ml of distilled water (*)	7 days at room T° 48 hours at room T°	30 minutes <i>on the automate board</i> before any use (**)
Quality controls Biophen Organan ® Control (ref A223501)	1 ml of distilled water (*)	7 days at room T° 48 hours at room T°	** 30 minutes <i>on the automate board</i> before any use (**)

Reconstitution:

(\*) After reconstitution of calibrators or controls with distilled water, let them to stabilize for 30 minutes at room temperature. It is better to reconstitute calibrators the very day of calibration.

### Conservation of reagents:

(\*\*) Take care of strictly respecting the 30 minutes temperature stabilization time for *calibrators* and *controls* at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C. Homogenize before each use.

**Nota:** For lyophilized calibrators and controls, following reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature. It is recommended to run the calibration curve with a freshly reconstituted calibrator. It is necessary to let the reagent temperature to stabilize for at least 30 minutes onto the automate before any use. Take care avoiding any contamination or evaporation of the reagents. Stability can be adjusted according to the exact use conditions.

**Homogeneize before each use.**

**Do not freeze calibrators and quality control plasmas.**

**A new calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).**

**Performances may present slight variations according to the instrument used. Validate the expected values in the laboratory working conditions.**

**Check the validity of the series by including quality control plasmas at different levels in each one.**

**Quality controls must be run regularly, and for each new batch of reagents, after an important maintenance of the instrument, or if measured values are not in compliance with the one expected for the method.**

### **5. Results:**

- The calibration curve is of the Log (absorbance) - Lin (concentration) type.
- The values obtained for the patients and controls are directly calculated from the calibration curve.
- The results are expressed in IU/ml.
- When Heparin or Orgaran concentrations are out of the working range, assayed plasma must be diluted in normal plasma, appropriately prepared and platelet poor, in order to keep a sufficient concentration of AT III.
- In presence of low AT III concentrations, as it can be the case in young children, an exogenous source of AT III is necessary, in order to correctly measure the heparin concentration.

**Nota:** Biophen Heparin reagents are developed for measuring homogeneously Unfractionated heparin (UFH) and Low Molecular Weight Heparin (LMWH), using the same calibration curve.

The superimposition conditions are susceptible to slightly vary according to the combination of reagents and lots used, and the technical characteristics and specific adjustment of each apparatus.

Therefore, this superimposition has to be verified and validated in the exact laboratory working conditions and for each machine, by establishing the calibration curve using Biophen Heparin calibrators (LMWH, ref A222001), and by assaying the homogeneous system of associated Biophen UFH calibrators (ref A222301), UFH control plasmas (ref A223101), and LMWH control plasmas at high and low levels (ref A223001, A223701).

The linearity (r<sup>2</sup>) is expected to be ≥0.98 for the calibration curve.

Controls must be measured in compliance.

Measured values for calibrators are acceptable when:

Measured value = target value  $\pm 0.10$  IU/ml for levels  $\leq 0.50$  IU/ml (better  $\pm 0.05$  IU/ml)

Measured value = target value  $\pm 0.15$  IU/ml for levels  $> 0.5$  IU/ml (better 90-110% of the target value, or  $\pm 0.10$  IU/ml for levels between 0.5 and 1.0 IU/ml)

Should the superimposition not be obtained, a separate calibration curve must be used for each type of heparin, and validated by assaying the homogeneous quality control system (controls measured in compliance and close to the assigned target value).

The calibration curve is validated when linearity ( $r^2$ ), as well as measured control values, are in compliance.

A new calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).

Performances may present slight variations according to the instrument used. Validate the expected values in the laboratory working conditions. Performances, as well as values for each new lot of quality controls used, must then be confirmed (and adjusted if necessary) in the laboratory working conditions.

**Nota:** The 3 months stability at 2-8°C and the 7 days stability at room temperature (18-25°C) claimed are validated on reconstituted vials kept closed, protected from any evaporation or contamination.

For your information, stability data on STAR board tested in our laboratory working conditions are:

-7 days for reconstituted vials loaded on STAR during working time, using mini reducers or chimneys then closed and stored at 2-8 °C overnight.

-5 days for reconstituted vials left loaded continuously on STAR.

Stabilities of reagents on each automate must be controlled, adjusted and validated under each exact laboratory working conditions.

## 6. Programming the BCS analyser:

### Creation of reagents and diluent files

Return to definition of test. Go in to reagents without data of lot. Create the 2 reagent files, as well as that of Diluent.

<b>Name:</b> Substrate BIOPHEN HEPARIN6	<b>Name:</b> Factor Xa Biophen HEPARIN6
Abbreviation: SUBHEP	Abbreviation: FXAHEP
clean when the reagent is modified high	clean when the reagent is modified special
<input checked="" type="checkbox"/> New calibration with any new lot of reagents	<input checked="" type="checkbox"/> New calibration with any new lot of reagents
<b>Authorized lines</b> from 1 to 14	<b>Authorized lines</b> from 1 to 14
<b>Stability(*)</b> < 15C      > = 15C	<b>Stability (*)</b> < 15C      > = 15C
<b>Mix by jet</b> Intensity 1 Frequency No agitation	<b>Mix by jet</b> Intensity 1 Frequency No agitation
<b>Reference numbers</b> Number add    remove    replace	<b>Reference numbers</b> Number add    remove    replace

<b>Name:</b> PHYSAL
Abbreviation: SP
clean when the reagent is modified Normal/low
<input type="checkbox"/> New calibration with any new lot of reagents
<b>Authorized lines</b> from 1 to 14
<b>Stability (*)</b> < 15C      > = 15C
<b>Mix by jet</b> Intensity 1 Frequency No agitation
<b>Reference numbers</b> Number add    remove    replace

Create in the same way calibration plasmas and control plasmas files.

<b>Name:</b> Biophen HEPARIN Calibrator
Abbreviation: HEPICAL
clean when the reagent is modified Normal/low
<input checked="" type="checkbox"/> New calibration with any new lot of reagents
<b>Authorized lines</b> from 1 to 14
<b>Stability (*)</b> < 15C      > = 15C
<b>Mix by jet</b> Intensity 1 Frequency No agitation
<b>Reference numbers</b> Number add    remove    replace

<b>Name:</b> Biophen plasma UFH Control
Abbreviation: UFHCONT
clean when the reagent is modified Normal/low
<input checked="" type="checkbox"/> New calibration with any new lot of reagents
<b>Authorized lines</b> from 1 to 14
<b>Stability (*)</b> < 15C      > = 15C
<b>Mix by jet</b> Intensity 1 Frequency No agitation
<b>Reference numbers</b> Number add    remove    replace

(\*) To be filled by the user with stabilities data provided into the kit insert.

Repeat these operations as much as necessary to create all calibrator and control plasmas files.

**Assignment of the test to controls**

Assignment of control: UFH	
Test: BIOPHEN HEPARIN	
Type of control	
<input checked="" type="checkbox"/> Exactitude	Precision
Confidence interval	
Allowed variation for the ref.	* <input type="text"/> %
Automatic control measurement of check	
<input type="checkbox"/> No automatic control measurement	
<input type="button" value="cancel"/>	<input type="button" value="record"/>

• **Foot-note**

" **Allowed variation of reference** " must be reported in percentage, in a way adapted for each Quality control used; precision of the method being different according to the measured concentrations, introduce this information with a value in % of the target value claimed for the QC and chosen in order to cover the range of acceptance for each control.

For example, allowed variation value authorized for UFH1 with 0.20IU/ml as target value and 0.12 to 0.28 as acceptance range will be programmed at 40%.

**Creation of the procedure**

Once the various reagents files are created, the procedure file must be created. In order to do it, go to **definition of tests**, procedure and type the following parameters.

<b>General Information</b>			
Former Change	<input type="text"/>		
Test Procedure Number	<input type="text"/>	<input type="text"/>	
Test Procedure Name	BIOPHEN HEPARIN6		
Measurement time stopped after	100 sec	Or	mA (0 = anything time of measure)
	0		
Speed	normal		
Stirring	normal		
Wavelength	405 nm		
Primary absorbance range	131-160 mA		
Duplicate n# for samples and controls (*)	<input type="text"/>	CV Authorized (*)	<input type="text"/> %
Duplicate n# for calibration (*)	<input type="text"/>	CV Authorized (*)	<input type="text"/> %
Duplicate n# of summary of raw data	Arithmetic mean		
Dilution Factor	1 : 2	(1 = no dilution)	
Dilution Buffer	PHYSAL		

(\*) To be filled by the user

Evaluation and Control Method

DeltaA per minute

from  Sec

Sec

Correlation Factor

Absolute

**Pipetting Sequence**

Test Procedure Number

Test Procedure Name

Sampling cycles

Number of Cycle	Arm of transmission	Final Washing	Action of Rotors/transfer	Time Range [sec]	
				Min	Max
1	Arm	normal	No mix. in priority	0	0
2	Reagent Arm	normal	Incubation	160	180
3	Reagent Arm	High	No mix/measure	0	0
4	Reagent Arm	normal	No mix/measure	0	0

Num de transm.	Mean/pipetting goal	Reagent	Speed	Volume (µl)
1	Take air		Slow	20
1	Take the sample		Slow	50
1	In the external cup		Slow	50
2	Take air		Slow	20
2	Take the Reagent	R1 SUPHEP BIO	Slow	100
2	In the central cup		Slow	100
3	Take air		Slow	20
3	Take the Reagent	R2 FXAHEP BIO	Slow	125
3	In the internal cup		Slow	125
4	Take the Reagent	SCS Cleaner	Slow	0

### Test Definition

<b>General Information</b>						
Former Change						
Name		BIOPHEN HEPARIN6		°Test N°		
Result Unit				PC ID t		
Behring	IU/ml	Conversion factor		Abbreviation		
User	IU/ml	1IU/ml	1.0	IU/ml	Place in the list	
PC	IU/ml	1IU/ml	1.0	IU/ml		
					Digit number for results	1000

### Test detailed procedures

Procedure Name test	Method of evaluation
BIOPHEN HEPARIN6	Delta A per minute

### Calculate Formula of raw data

Formula Type	Minimum	
User owns Formula		
Raw data unit	mE/min	

### Detailed test evaluation

Reference range	Determination of result		<b>Standard Concentrations</b>	
	0.0	2.0		IU/ml
	% dN			IU/ml
	at			
Test used for Calibration				
Evaluation	Regression lin/log			
Reference Curve	Calibration curve			
Minimum raw value	0.0			
Upper Extrapolation	1.25	x maxi concentration		
Lower Extrapolation	0.5	x mini concentration		
Calibrator	BIOPHEN HEPARIN calibrator			
Authorized SD respectively to the reference curve	0	%		
Maximal number repeats	0			

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