**1. Reconstitution of the BIOPHEN Heparin 6 reagent, (Reference A221006)**

Chromogenic determination of the Anti-Xa activity of LMWH, UFH, and Orgaran.

	NAME	Reconstitution	Stability	Stabilization in T°
R1	SXa-11 Substrate	7.5 ml of distilled water (*)	3 months at 2-8°C ** 7 days at room T° Do not freeze	30 mn <i>on board</i> before any use(***)
R2	Factor Xa	7.5 ml of distilled water (*)	3 months at 2-8°C ** 7 days at room T° Do not freeze	30 mn <i>on board</i> before any use (***)
Diluent	Physiological Saline		24 hours on board	

**2. Reconstitution of the BIOPHEN Heparin 3 reagent, (Reference A221003)**

Chromogenic determination of the Anti-Xa activity of LMWH, UFH, and Orgaran.

	NAME	Reconstitution	Stability	Stabilization in T°
R1	SXa-11 Substrate	3.75 ml of distilled water (*)	3 months at 2-8°C ** 7 days at room T° Do not freeze	30 mn <i>on board</i> before any use(***)
R2	Factor Xa	3.75 ml of distilled water (*)	3 months at 2-8°C ** 7 days at room T° Do not freeze	30 mn <i>on board</i> before any use (***)
Diluent	Physiological Saline		24 hours on board	

Recommendations for reconstitution:

**Reconstitution (\*):** Reconstitute each vial with distilled water. Shake thoroughly (vortex). Let to homogenize for 30 minutes at room temperature (18-25 °C), while shaking the vial from time to time (vortex), until complete dissolution of the content. Check the absence of any solid at the bottom of the vial.

In current practice, in order to allow a good standardization, reconstitute these two reagents the evening before and put them at 2°-8°C following the 30 minutes at room temperature.

Note: In all cases, before use, check the absence of solids at the bottom of the vial, which confirms that dissolution is complete. If necessary, incubate for 1 hour at RT or better at 37°C, while shaking (vortex) from time to time. If required, then additionally incubate overnight at RT.

**Storage of reagents: (\*\*)** Stability indicated are validated provided any evaporation or contamination of reconstituted reagents avoided. They must be adjusted according to the actual laboratory working conditions

Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C, and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

If the reagents are kept on the automate board 24H/24H, take care to limit as much as possible any evaporation of the reagents by using a chimney or mini-reducer.

**Stabilization of reagents:** (\*\*\*) It is necessary to let the substrate and the Factor Xa temperature to stabilize for at least 30 minutes on the automate before any use. A too low temperature of the reagents can induce an over-estimation. Conversely a too high temperature leads to an under evaluation of heparin.

**Homogenize the reagents before each use.**

Any reagent of biological origin must be handled with all the required cautions, as being potentially infectious.

**Do not interchange the reagents from different lots.**

**If necessary use microcups when running the assay, to improve management of the reagents volumes by the automate**

### 5. Determination of Heparins UFH and/or LMWH.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Heparin Calibrator (ref A222001)	1 ml of distilled water (*)	7 days at 2-8°C 48 hours at room T°	30 minutes <i>on</i> the automate <i>board</i> before any use (**)
Quality controls Biophen UFH Control (ref A223101) Biophen LMWH Control (ref A223001)	1 ml of distilled water (*)	7 days at 2-8°C 48 hours at room T°	30 minutes <i>on</i> the automate <i>board</i> before any use (**)

### 6. Determination of Orgaran:

The determination of Orgaran requires its own configuration of the BCT. This one differs from the configuration used for the *determination* of UFH/LMWH only by the denomination of calibrators and controls.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Orgaran ® Calibrator (ref A222201)	1 ml of distilled water (*)	7 days at room T° 48 hours at room T°	30 minutes <i>on</i> the automate <i>board</i> before any use (**)
Quality controls Biophen Orgaran ® Control (ref A223501)	1 ml of distilled water (*)	7 days at room T° 48 hours at room T°	** 30 minutes <i>on</i> the automate <i>board</i> before any use (**)

Reconstitution:

(\*) After reconstitution of calibrators or controls with distilled water, let them to stabilize for 30 minutes at room temperature. It is better to reconstitute calibrators the very day of calibration.

Conservation of reagents:

(\*\*) Take care of strictly respecting the 30 minutes temperature stabilization time for *calibrators* and *controls* at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C. Homogenize before each use.

**Nota:** For lyophilized calibrators and controls, following reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature. It is recommended to run the calibration curve with a freshly reconstituted calibrator. It is necessary to let the reagent temperature to stabilize for at least 30 minutes onto the automate before any use. Take care avoiding any contamination or evaporation of the reagents. Stability can be adjusted according to the exact use conditions.

**Homogeneize before each use.**

**Do not freeze calibrators and quality control plasmas.**

**A new calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).**

**Performances may present slight variations according to the instrument used. Validate the expected values in the laboratory working conditions.**

**Check the validity of the series by including quality control plasmas at different levels in each one.**

**Quality controls must be run regularly, and for each new batch of reagents, after an important maintenance of the instrument, or if measured values are not in compliance with the one expected for the method.**

## 7. Results:

- The calibration curve is of the Log (absorbance) - Lin (concentration) type.
- The values obtained for the patients and controls are directly calculated from the calibration curve.
- The results are expressed in IU/ml.
- When Heparin or Orgaran concentrations are out of the working range, assayed plasma must be diluted in normal plasma, appropriately prepared and platelet poor, in order to keep a sufficient concentration of AT III.
- In presence of low AT III concentrations, as it can be the case in young children, an exogenous source of AT III is necessary, in order to correctly measure the heparin concentration.

**Nota:** Biophen Heparin reagents are developed for measuring homogeneously Unfractionated heparin (UFH) and Low Molecular Weight Heparin (LMWH), using the same calibration curve.

The superimposition conditions are susceptible to slightly vary according to the combination of reagents and lots used, and the technical characteristics and specific adjustment of each apparatus.

Therefore, this superimposition has to be verified and validated in the exact laboratory working conditions and for each machine, by establishing the calibration curve using Biophen Heparin calibrators (LMWH, ref A222001), and by assaying the homogeneous system of associated Biophen UFH calibrators (ref A222301), UFH control plasmas (ref A223101), and LMWH control plasmas at high and low levels (ref A223001, A223701).

The linearity ( $r^2$ ) is expected to be  $\geq 0.98$  for the calibration curve.

Controls must be measured in compliance.

Measured values for calibrators are acceptable when:

Measured value = target value  $\pm 0.10$  IU/ml for levels  $\leq 0.50$  IU/ml (better  $\pm 0.05$  IU/ml)

Measured value = target value  $\pm 0.15$  IU/ml for levels > 0.5 IU/ml (better 90-110% of the target value, or  $\pm 0.10$  IU/ml for levels between 0.5 and 1.0 IU/ml)

Should the superimposition not be obtained, a separate calibration curve must be used for each type of heparin, and validated by assaying the homogeneous quality control system (controls measured in compliance and close to the assigned target value).

The calibration curve is validated when linearity ( $r^2$ ), as well as measured control values, are in compliance.

A new calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).

Performances may present slight variations according to the instrument used. Validate the expected values in the laboratory working conditions. Performances, as well as values for each new lot of quality controls used, must then be confirmed (and adjusted if necessary) in the laboratory working conditions.

**Nota:** The 3 months stability at 2-8°C and the 7 days stability at room temperature (18-25°C) claimed are validated on reconstituted vials kept closed, protected from any evaporation or contamination.

For your information, stability data on STAR board tested in our laboratory working conditions are:

-7 days for reconstituted vials loaded on STAR during working time, using mini reducers or chimneys then closed and stored at 2-8 °C overnight.

-5 days for reconstituted vials left loaded continuously on STAR.

Stabilities of reagents on each automate must be controlled, adjusted and validated under each exact laboratory working conditions.

## 8. Programming the BCT analyzer:

### Creation of reagents and diluent files

Return to definition of test. Go in to reagents without data of lot. Create the 2 reagent files, as well as that of Diluent.

#### A. Create Reagents

##### Substrat

Reagent Parameter	
Name	HNF/ SUB BIOPHEN
Abbreviation	UFH/LMWH
Type	Normal
Version :	
<b>Reagents Data -</b>	
Stability :	99 h
<input checked="" type="checkbox"/> New calibration with any new lot of reagents <input type="checkbox"/> Use as enzyme blank <input type="checkbox"/> Rerun enzyme blank with any new lot	
<b>Washing</b>	
<input checked="" type="checkbox"/> Normal      Wash Time	10 s
<input type="checkbox"/> Intensive <input type="checkbox"/> Supplementary washing step	
<b>MIX</b>	
<input checked="" type="checkbox"/> Disabled <input type="checkbox"/> Enabled <input type="checkbox"/> Interval	

##### FXa

Reagent Parameter	
Name	HNF/HBPM FXa BIOPHEN
Abbreviation	UFH/LMWH FXa
Type	Normal
Version :	
<b>Reagents Data -</b>	
Stability :	99 h
<input checked="" type="checkbox"/> New calibration with any new lot of reagents <input type="checkbox"/> Use as enzyme blank <input type="checkbox"/> Rerun enzyme blank with any new lot	
<b>Washing</b>	
<input type="checkbox"/> Normal      Wash Cycle	5
<input checked="" type="checkbox"/> Intensive <input type="checkbox"/> Supplementary washing step	
<b>MIX</b>	
<input checked="" type="checkbox"/> Disabled <input type="checkbox"/> Enabled <input type="checkbox"/> Interval	

##### Buffer

Reagent Parameter	
Name	PHYSAL
Abbreviation	SP
Type	Buffer
Version :	
<b>Reagents Data -</b>	
Stability :	24 h
<input type="checkbox"/> New calibration with any newlot of reagents <input type="checkbox"/> Use as enzyme blank <input type="checkbox"/> Rerun enzyme blank with any new lot	
<b>Washing</b>	
<input checked="" type="checkbox"/> Normal      Wash Cycle	1.0 s
<input type="checkbox"/> Intensive <input type="checkbox"/> Supplementary washing step	
<b>MIX</b>	
<input checked="" type="checkbox"/> Disabled <input type="checkbox"/> Enabled <input type="checkbox"/> Interval	

## Create the program for UFH/LMWH

Test Configuration	
Name :	UFH/LMWH BIOPHEN
Abbreviation :	UFH/LMWH
PC ID	
Evaluation Mode:	
OD KINETIC without enzyme blank	

Click on *Pipetting sequence and write:*

Run Parameter - UFH/LMWH			
Sample Vol :		50	µl
Buffer :	SP	1+	1
<input type="checkbox"/> Sample starts reading			
Deficient Plasma	without	0	µl
Reagent 1	without	0	µl
<input type="checkbox"/> Pre-incubation		0	s
Reagent 2	UFH/LMWH SUB	125	µl
Reagent 3		0	µl
Incubation time :		180	s
Start Reagent :	UFH/LMWH F.Xa	125	µl
Reading Time		50	s
Determination :	<input checked="" type="checkbox"/> Single	<input type="checkbox"/> Double	
Max Deviation :	10	%	

Click on *Rerun mode*

Rerun Mode UFH/LMWH		
Evaluation error	Rerun Mode	Data
1. Too high Deviation	no rerun	
2. No intersection	no rerun	
3. Extrapolation up	no rerun	
4. Extrapolation down	no rerun	
5. No base line	no rerun	
6. Invalid Parameters	no rerun	
7. Too high or low OD	no rerun	
8. Incorrect OD	no rerun	

Rerun Modes enabled

Click on **Evaluation** and write.

OD KINETIC without enzyme blank	
Latence phase:	<input type="text" value="30"/> s
Late phase :	<input type="text" value="5"/> s
<input type="checkbox"/> To negative values	
<input type="checkbox"/> To screen inverse graph in research mode	

Click on **Units** and write the parameters and valid.

Units and values UFH/LMWH	
<b>Units</b>	
Raw Data :	mE/min
result Unit :	<input type="text" value="IU/ml"/>
Decimal Number:	<input type="text" value="2"/>
<b>Calibration</b>	
<input type="checkbox"/> No Calibration	
<input checked="" type="checkbox"/> Normal calibration	
<input type="checkbox"/> Calibration Ratio	
<input type="checkbox"/> Ratio/INR Calibration	
<input type="text"/>	
<b>Normal Values</b>	
<input checked="" type="checkbox"/> Normal range for result	
Minimum :	<input type="text" value="0,01"/>
Maximum :	<input type="text" value="2,00"/>

Click on Calibration and enter following parameters. Refer to the certificate of analysis for all information needed.

UFH/LMWH Calibration					
Type :	Plasma Calibration				
Calibration data					
Durability	Not limited			name	Conc IU /ml
Alarm hour			1	CAL 1	
Points number	5		2	CAL2	
Determination	2		3	CAL 3	
Time	50	S	4	CAL 4	
Max deviation	10	%	5	CAL 5	
Approximation					
Mode	Regression				
	LOG-LIN				
Extrapolation	Low		1		
	High		1		

Click on **Technic**, write and valid the process by OK.

Technical Parameters		
Photometer		
Wavelength :	<input type="text" value="405 nm"/>	
<input type="checkbox"/> High sensitivity		
<input type="checkbox"/> High resolution		
Mix		
Initial speed:	<input type="text" value="400"/>	1/mm
Change after :	<input type="text" value="4"/>	s
Stop speed	<input type="text" value="0"/>	1/mm
Pipeting		
Reflux	<input type="text" value="3.0"/>	µl
Separation Bubble :	<input type="text" value="5.0"/>	µl
Diluted speed distribution	<input type="text" value="4000"/>	Phases/s

Click on the icon set up software and create the program for quality control.

Quality control test creation

Click on quality control then on the table; manage the different levels of quality control with name, lot, ect. (Refer to certificate of analysis)

click on data and give a test number for different level of controls

Select UFH/LMWH test with following parameters:

-Quality mode

-Rerun mode

-Work method: upon request

For all information concerning target values and limits of controls refer to the documentation supplied with the kit (insert, certificate of analysis)

Do the same thing for each controls of the kit.

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