Manufactured By: HYPHEN BioMed

Biophen Heparin (AT+) A221007 on ACL7000 (research software)

NOTE : reconstitution, preparation, calibration and controls, expression of results: usual cautions and realization according to the specific technical insert D.750.02/BI/1007.

Reagents are reconstituted by 7.5ml of distilled water for R1/R2 and by 2.5ml of distilled water for R3, well equilibrated, and preheated at 37°C

• **For UFH or LMWH assay:**

PHASE	CONDITION	OPERATION	PARAMETERS
1 st loading step	YES	Reagent loading (pos 18)	0µl
		Washing after reagent loading (pos 18)	0
		Diluent (pos DIL) Reagent R3 AT(h)	16µL
		Specimen	16µl
		Reagent (nacelle 2) Reagent R1 Substrate	80µL
		Washing between specimen	3
		Optical ref. (cycle abs.)	
		Diluent (pos DIL)	80µL
		Reagent (nacelle 2)	80µL
		Washing (at the end of loading)	3
2 nd loading step	YES	Incubation time (T1)	120 s
		Diluent (pos DIL)	0µL
		Specimen	0µL
		Washing between specimen	0
		Reagent (nacelle 3) Reagent R2 FXa	80µL
		Washing after reagent loading	3
		Mix time (T2)	120 s
3 rd loading step	NO	Activation time (T3)	0 s
		Reagent (pos 17)	0µL
		Washing after reagent loading	0
		Mix time (T4)	0 s
Washing	NO	Activation time (T5)	0 s
		Reagent (pos 1)	0µL
		Reagent (pos 16)	0µL
Levels detection	NO	Cycles number	1
		Washing at the end of loading	1
Acquisition parameters	YES	Ramp	yes
		Inter-ramp interval (T6)	1s
		Delay before acquisition (T7)	5s
		Acquisition time (T8)	95s (*)
		Canal used	405nm
		Rotation speed	1200 rpm

(*) the acquisition time is limited to 90 sec by programming (T8) at 95 sec.

Note: Biophen Heparin AT(+) reagents are developed for measuring homogeneously Unfractionated heparin (UFH) and Low Molecular Weight Heparin (LMWH), using the same calibration curve.

The superimposition conditions are susceptible to slightly vary according to the combination of reagents and lots used, and the technical characteristics and specific adjustment of each apparatus.

Therefore, this superimposition has to be verified and validated in the exact laboratory working conditions and for each machine, by establishing the calibration curve using Biophen Heparin calibrators (LMWH, ref A222001), and by assaying the homogeneous system of associated Biophen UFH calibrators (ref A222301), UFH control plasmas (ref A223101), and LMWH control plasmas at high and low levels (ref A223001, A223701).

The linearity (r^2) is expected to be ≥ 0.98 for the calibration curve.

Controls must be measured in compliance.

Measured values for calibrators are acceptable when:

Measured value = target value ± 0.10 IU/ml for levels ≤ 0.50 IU/ml (better ± 0.05 IU/ml)

Measured value = target value ± 0.15 IU/ml for levels > 0.5 IU/ml (better 90-110% of the target value, or ± 0.10 IU/ml for levels between 0.5 and 1.0 IU/ml)

Should the superimposition not be obtained, a separate calibration curve must be used for each type of heparin, and validated by assaying the homogeneous quality control system (controls measured in compliance and close to the assigned target value).

Example of calibration curve obtained:

	Acquisition 90 sec
IU/ml LMWH	A405
0	1.01
0.46	0.63
0.93	0.41
1.46	0.28
1.90	0.19
r^2 (logA405-lin conc)	>0.99

Note: The stability at 2-8°C and at room temperature (18-25°C) claimed are validated on reconstituted vials kept closed, protected from any evaporation or contamination.

Stabilities of reagents on each automate must be controlled, adjusted and validated under each exact laboratory working conditions.

- **For Arixtra assay:**

For assaying Arixtra in plasma samples in the range 0 to about 1.5µg/ml Arixtra, it is recommended to predilute the sample 1:2 in physiological saline before loading it on the ACL7000 for the test.

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