#### Determination of anti-Xa chromogenic activity of heparin and analogues

**NOTE :** reconstitution, preparation, calibration and controls, expression of results: usual cautions and realization according to the specific technical insert (D.750.02/BI/1011) and to the specific calibrators and controls used for the assayed molecule. (Ref. A221011)

#### Reagents are preheated at 37°C

PHASE	CONDITION	OPERATION	PARAMETERS
1 <sup>st</sup> loading step	YES	Reagent loading (pos 18)	0µl
		Washing after reagent loading (pos 18)	0
		<b>Diluant (pos DIL) Physiological saline</b>	<b>20µL</b>
		<b>Tested Specimen (loaded undiluted)</b>	<b>20µl</b>
		<b>Reagent (nacelle 2) Substrate R1</b>	<b>100µL</b>
		Washing between specimen	3
		Optical ref. (cycle abs.)	
		Diluent (pos DIL)	80µL
		Reagent (nacelle 2)	80µL
		Washing (at the end of loading)	3
		Incubation time (T1)	120 s
2 <sup>nd</sup> loading step	YES	Diluent (pos DIL)	0µL
		Specimen	0µL
		Washing between specimen	0
		<b>Reagent (nacelle 3) FXa R2</b>	<b>100µL</b>
		Washing after reagent loading	3
		Mix time (T2)	120 s
		Activation time (T3)	0 s
3 <sup>rd</sup> loading step	NO	Reagent (pos 17)	0µL
		Washing after reagent loading	0
		Mix time (T4)	0 s
		Activation time (T5)	0 s
Washing	NO	Reagent (pos 1)	0µL
		Reagent (pos 16)	0µL
		Cycles number	1
		Washing at the end of loading	1
Levels detection	NO		
Acquisition parameters	YES	Ramp	yes
		Inter-ramp interval (T6)	1s
		Delay before acquisition (T7)	<b>5s</b>
		Acquisition time (T8)	<b>65s</b>
		Canal used	405nm
		Rotation speed	1200 rpm

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