Manufactured By: HYPHEN BioMed

## BIOPHEN HEPARIN LRT

### ON BCS XP

### Adaptation of BIOPHEN HEPARIN LRT ON BCS XP

#### 1. Preparation of BIOPHEN Heparin (LRT) reagent, (Reference A221011)

Both reagents, chromogenic substrate SXa-11 and Factor Xa, included in Biophen Heparin LRT kit, are ready to use, for the assay of UFH/LMWH (see chapter II) or Orgaran® (see chapter III). A calibration curve and a specific configuration will be use for each type of antithrombotic substance.

<b>R1 : SXa-11 Substrate</b> (liquid, ready to use).
<b>R2: Factor Xa</b> (liquid, ready to use).
Diluent (not provided): physiological saline.

- **Recommendation for stabilization of the reagents:**

Before use, let the reagents to stabilize for 30 minutes at room temperature (18-25°C), while shaking from time to time.

**It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate board before any use.**

- **Stability data indicated on the insert are validated on vials kept closed, and provided any evaporation or contamination of reagents is avoided. They must be verified and adjusted if required according to the actual laboratory working conditions**

**Storage of reagents:** Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C, and of well homogenising and strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

If the reagents are kept on the automate board, take care and use reducers to limit as much as possible any evaporation of the reagents.

- **Homogenize the reagents before each use.**
- **Any reagent of biological origin must be handled with all the required cautions, as being potentially infectious.**
- **Do not interchange the reagents from different lots.**

## 2. Determination of Heparins UFH and/or LMWH.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Heparin Calibrator (ref A222001)	1 ml of distilled water (*)	7 days at 2-8°C  48 hours at room T°	30 minutes <i>on board</i> before any use (**)
Quality controls Biophen UFH Control (ref A223101) Biophen LMWH Control (ref A223001)	1 ml of distilled water (*)	7 days at 2-8°C  48 hours at room T°	30 minutes <i>on BCS board</i> before any use (**)

## 3. Determination of Orgaran:

The determination of Orgaran requires its own configuration of the BCS. This one differs from the configuration used for the *determination* of UFH/LMWH only by the denomination of calibrators and controls.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Orgaran ® Calibrator (ref A222201)	1 ml of distilled water (*)	7 days at room T°  48 hours at room T°	30 minutes <i>on board</i> before any use (**)
Quality controls Biophen Orgaran ® Control (ref A223501)	1 ml of distilled water (*)	7 days at room T°  48 hours at room T	30 minutes <i>on board</i> before any use (**)

### Reconstitution:

(\*) After reconstitution of calibrators or controls with distilled water, let them to stabilize for 30 minutes at room temperature. It is better to reconstitute calibrators the very day of calibration.

### Stabilization

(\*\*) Take care of strictly respecting the 30 minutes temperature stabilization time for *calibrators* and *controls* at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C. Homogenize before each use.

**Nota:** For lyophilized calibrators and controls, following reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature. It is recommended to run the calibration curve with a freshly reconstituted calibrator. It is necessary to let the reagent temperature to stabilize for at least 30 minutes onto the automate before any use. Take care avoiding any contamination or evaporation of the reagents. Stability can be adjusted according to the exact use conditions.

**Homogenize before each use.**

**Do not freeze calibrators and quality control plasmas.**

**A new calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).**

**Performances may present slight variations according to the instrument used. Validate the expected values in the laboratory working conditions.**

**Check the validity of the series by including quality control plasmas at different levels in each one.**

**Quality controls must be run regularly, and for each new batch of reagents, after an important maintenance of the instrument, or if measured values are not in compliance with the one expected for the method.**

#### 4. **Results:**

- The calibration curve is of the Log (absorbance) - Lin (concentration) type.
- The values obtained for the patients and controls are directly calculated from the calibration curve.
- The results are expressed in IU/ml.
- When Heparin or Orgaran concentrations are out of the working range, assayed plasma must be diluted in normal plasma, appropriately prepared and platelet poor, in order to keep a sufficient concentration of AT III.
- In presence of low AT III concentrations, as it can be the case in young children, an exogenous source of AT III is necessary, in order to correctly measure the heparin concentration.

**Nota:** Biophen Heparin LRT reagents are developed for measuring homogeneously Unfractionated heparin (UFH) and Low Molecular Weight Heparin (LMWH), using the same calibration curve.

The superimposition conditions are susceptible to slightly vary according to the combination of reagents and lots used, and the technical characteristics and specific adjustment of each apparatus.

Therefore, this superimposition has to be verified and validated in the exact laboratory working conditions and for each machine, by establishing the calibration curve using Biophen Heparin calibrators (LMWH, ref A222001), and by assaying the homogeneous system of associated Biophen UFH calibrators (ref A222301), UFH control plasmas (ref A223101), and LMWH control plasmas at high and low levels (ref A223001, A223701).

The linearity ( $r^2$ ) is expected to be  $\geq 0.98$  for the calibration curve.

Controls must be measured in compliance.

Measured values for calibrators are acceptable when:

Measured value = target value  $\pm 0.10$  IU/ml for levels  $\leq 0.50$  IU/ml (better  $\pm 0.05$  IU/ml)

Measured value = target value  $\pm 0.15$  IU/ml for levels  $> 0.5$  IU/ml (better 90-110% of the target value, or  $\pm 0.10$  IU/ml for levels between 0.5 and 1.0 IU/ml)

**Should the superimposition not be obtained, a separate calibration curve must be used for each type of heparin, and validated by assaying the homogeneous quality control system (controls measured in compliance and close to the assigned target value).**

**The calibration curve is validated when linearity ( $r^2$ ), as well as measured control values, are in compliance.**

**A new calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).**

**Performances may present slight variations according to the instrument used. Validate the expected values in the laboratory working conditions. Performances, as well as values for each new lot of quality controls used, must then be confirmed (and adjusted if necessary) in the laboratory working conditions.**

**Nota:** The stability at 2-8°C and at room temperature (18-25°C) claimed are validated on reconstituted vials kept closed, protected from any evaporation or contamination. Stabilities of reagents on each automate must be controlled, adjusted and validated under each exact laboratory working conditions.

**Programming the BCS XP analyser:**  
**Creation of reagents and diluent files**

Return to definition of test. Go into reagents without data of lot. Create the 2 reagent files, as well as that of Diluent.

<b>Name:</b>	Substrate BIOPHEN HEPARIN LRT		
Abbreviation:	SUBHEPLRT		
clean when the reagent is modified	high		
<input checked="" type="checkbox"/>	New calibration with any new lot of reagents		
<b>Authorized lines</b>	from 1 to 14		
<b>Stability (*)</b>	< 15C > = 15C		
<b>Mix by jet</b>	Intensity 1 Frequency No agitation		
<b>Reference numbers</b>	Number add remove replace		

<b>Name:</b>	Factor Xa Biophen HEPARIN LRT		
Abbreviation:	FXAHEPLRT		
clean when the reagent is modified	special		
<input checked="" type="checkbox"/>	New calibration with any new lot of reagents		
<b>Authorized lines</b>	from 1 to 14		
<b>Stability (*)</b>	< 15C > = 15C		
<b>Mix by jet</b>	Intensity 1 Frequency No agitation		
<b>Reference numbers</b>	Number add remove replace		

<b>Name:</b>	Physal		
Abbreviation:	SP		
clean when the reagent is modified	Normal/low		
<input type="checkbox"/>	New calibration with any new lot of reagents		
<b>Authorized lines</b>	from 1 to 14		
<b>Stability (*)</b>	< 15C > = 15C		
<b>Mix by jet</b>	Intensity 1 Frequency No agitation		
<b>Reference numbers</b>	Number add remove replace		

Create in the same way calibration plasmas and control plasmas files.

<b>Name:</b>	Biophen HEPARIN Calibrator		
Abbreviation:	HEPCAL		
clean when the reagent is modified	Normal/low		
<input checked="" type="checkbox"/>	New calibration with any new lot of reagents		
<b>Authorized lines</b>	from 1 to 14		
<b>Stability (*)</b>	< 15C > = 15C		
<b>Mix by jet</b>	Intensity 1 Frequency No agitation		
<b>Reference numbers</b>	Number add remove replace		

<b>Name:</b>	Biophen plasma UFH Control		
Abbreviation:	UFHCONT		
clean when the reagent is modified	Normal/low		
<input checked="" type="checkbox"/>	New calibration with any new lot of reagents		
<b>Authorized lines</b>	from 1 to 14		
<b>Stability (*)</b>	< 15C > = 15C		
<b>Mix by jet</b>	Intensity 1 Frequency No agitation		
<b>Reference numbers</b>	Number add remove replace		

(\*) To be filled by the user with stabilities data provided into the kit insert.

Repeat these operations as much as necessary to create all calibrator and control plasmas files.

**Assignment of the test to controls**

Assignment of control: UFH	
Test: BIOPHEN HEPARIN	
Type of control	
<input checked="" type="checkbox"/> Exactitude	Precision
Confidence interval	
Allowed variation for the ref.	* %
Automatic control measurement of check	
<input type="checkbox"/> No automatic control measurement	
cancel	record

• **Foot-note**

" **Allowed variation of reference** " must be reported in percentage, in a way adapted for each Quality control used; precision of the method being different according to the measured concentrations, introduce this information with a value in % of the target value claimed for the QC and chosen in order to cover the range of acceptance for each control.

For example, allowed variation value authorized for UFH1 with 0.20IU/ml as target value and 0.12 to 0.28 as acceptance range will be programmed at 40%.

**Creation of the procedure**

Once the various reagents files are created, the procedure file must be created. In order to do it, go to **definition of tests**, procedure and type the following parameters.

General Information			
Former Change			
Test Procedure Number			
Test Procedure Name	BIOPHEN HEPARIN LRT		
Measurement time stopped after	100 sec	Or	mA (0 = anything time of measure)
	0		
Speed	normal		
Stirring	normal		
Wavelength	405 nm		
Primary absorbance range	131-160 mA		
Duplicate n# for samples and controls (*)		CV Authorized (*)	0.0 %
Duplicate n# for calibration (*)	2	CV Authorized (*)	10 %
Duplicate n# of summary of raw data	Arithmetic mean		
Dilution Factor	1 : 2	(1 = no dilution)	
Dilution Buffer	PHYSAL		

\*testing in duplicate is recommended at least for the calibration curve (with CV expected <10%, to be adjusted by the user).

Evaluation and Control Method			
Delta A per minute			
from	20	Sec	50
Correlation Factor			0.95
			Sec
			Absolute

AmaxMin: 10 sec; Minimal absorption limit 5000mA.

Pipetting Sequence					
Test Procedure Number					
Test Procedure Name			BIOPHEN HEPARIN LRT		
Sampling cycles					
Number of Cycle	Arm of transmission	Final Washing	Action of Rotors/transfer	Time Range [sec]	
				Min	Max
1	Sample Arm	normal	No mix. in priority	0	0
2	Reagent Arm	normal	Incubation	160	180
3	Reagent Arm	High	No mix/measure	0	0
Transmission number	Mean/pipetting goal	Reagent	Speed	Volume (µl)	
1	Take air		Slow	20	
1	Take the sample		Slow	50	
1	In the external cup		Slow	50	
2	Take air		Slow	20	
2	Take the Reagent	R1 SUBHEPLRT	Slow	125	
2	In the central cup		Slow	125	
3	Take air		Slow	20	
3	Take the Reagent	R2 FXAHEPLRT	Slow	125	
3	In the internal cup		Slow	125	

### Test Definition

<b>General Information</b>				
Former Change		<input type="text"/>		
Name		BIOPHEN HEPARIN LRT		°Test N°
Result Unit				PC ID t
Behring	IU/ml	Conversion Factor		Abbreviation
User	IU/ml	1IU/ml	1.0	IU/ml
PC	IU/ml	1IU/ml	1.0	IU/ml
				Place in the list
				Digit number for results

### Test detailed procedures

Procedure Name test	Method of evaluation
BIOPHEN HEPARIN LRT	Delta A per minute

### Calculate Formula of raw data

Formula Type	Minimum	<input type="text"/>
User owns Formula	<input type="text"/>	
Raw data unit	mE/min	<input type="text"/>

### Detailed test evaluation

Reference range	Determination of result		<b>Standard Concentrations</b>  IU/ml
	0.0* % dN at	2.0* IU/ml	
Test used for Calibration	<input type="text"/>		
Evaluation	Regression lin/log		
Reference Curve	Calibration curve		
Minimum raw value	0.0		
Upper Extrapolation	1.0	x maxi concentration	
Lower Extrapolation	1.0	x mini concentration	
Calibrator	BIOPHEN HEPARIN calibrator		
Authorized SD respectively to the reference curve	0	%	
Maximal number repeats	0		

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