Manufactured by Hyphen BioMed.

BIOPHEN HEPARIN LRT ON SYSMEX CA 1500

1. Preparation of BIOPHEN Heparin (LRT) reagent, (Reference A221011)

Both reagents, chromogenic substrate SXa-11 and Factor Xa, included in Biophen Heparin LRT kit, are ready to use, for the assay of UFH/LMWH (see chapter II) or Orgaran® (see chapter III). A calibration curve and a specific configuration will be use for each type of antithrombotic substance.

R1 : SXa-11 Substrate (liquid, ready to use).
R2: Factor Xa (liquid, ready to use).
Diluent (not provided): physiological saline.

- **Recommendation for stabilization of the reagents:**

Before use, let the reagents to stabilize for 30 minutes at room temperature (18-25°C), while shaking from time to time.

It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate board before any use.

- **Stability data indicated on the insert are validated on vials kept closed, and provided any evaporation or contamination of reagents is avoided. They must be verified and adjusted if required according to the actual laboratory working conditions**

Storage of reagents: Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C, and of well homogenising and strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

If the reagents are kept on the automate board, take care and use reducers to limit as much as possible any evaporation of the reagents.

- **Homogenize the reagents before each use.**
- **Any reagent of biological origin must be handled with all the required cautions, as being potentially infectious.**
- **Do not interchange the reagents from different lots.**

2. Determination of Heparins UFH and/or LMWH.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Heparin Calibrator (ref A222001)	1 ml of distilled water (*)	7 days at 2-8°C 48 hours at room T°	30 minutes on board SYSMEX CA 1500 before any use (**)
Quality controls Biophen UFH Control (ref A223101) Biophen LMWH Control (ref A223001)	1 ml of distilled water (*)	7 days at 2-8°C 48 hours at room T°	30 minutes on board SYSMEX CA 1500 before any use (**)

3. Determination of Orgaran:

The determination of Orgaran requires its own configuration of the CA 1500. This one differs from the configuration used for the *determination* of UFH/LMWH only by the denomination of calibrators and controls.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Orgaran® Calibrator (ref A222201)	1 ml of distilled water (*)	7 days at room T° 48 hours at room T°	30 minutes on board SYSMEX CA 1500 before any use (**)
Quality controls Biophen Orgaran® Control (ref A223501)	1 ml of distilled water (*)	7 days at room T° 48 hours at room T°	30 minutes on board SYSMEX CA 1500 before any use (**)

Reconstitution:

(*) After reconstitution of calibrators or controls with distilled water, let them to stabilize for 30 minutes at room temperature. It is better to reconstitute calibrators the very day of calibration.

Conservation of reagents:

(**) Take care of strictly respecting the 30 minutes temperature stabilization time for *calibrators* and *controls* at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C. Homogenize before each use.

Foot-note: Do not freeze calibrators or controls.

Foot-note: A calibration curve must be carried out for each new batch of reagents.

4. Results:

- The calibration curve is of the Log (absorbance) - Lin (concentration) type.
- The values obtained for the patients and controls are directly calculated from the calibration curve.
- The results are expressed in IU/ml.
- When Heparin or Orgaran concentrations are out of the working range, assayed plasma must be diluted in normal plasma, appropriately prepared and platelet poor, in order to keep a sufficient concentration of AT III.
- In presence of low AT III concentrations, as it can be the case in young children, an exogenous source of AT III is necessary, in order to correctly measure the heparin concentration.

5. Instrument setting for sysmex CA 1500:

Click on the window Set up software for the “manager program” and create the program according to: Chromogen for AT3

Change the data with the data indicated on the table below

Use the program and create the name for the reagents, Calibration plasmas and control plasmas

System		Ready		
Parameter : HEP LRT Chr Para Code				
Sample.Vol			18	µl
Diluent Vol	Dist Water		36	µl
Wash	Without			
Second Dilution	No	0µl		
Diluent Vol	No	0µl		
Wash	Without			
Deficient Plasma	No			
Wash	Without	x0/	Without	x0
First reagent	SUB	90µl	20 sec	
Diluent Vol	No	0µl		
Wash	Without	x0/	Clean I x 1	
Second reagent	Xa	90µl	250 sec	Norm
Diluent Vol	No	0µl		
Wash	Without	x0/ Off / CleanI x 1		
Third reagent	Without			
Diluent Vol	No Util	0µl		
Wash	Without	x0/	Without x0	
Detector	Chromogenic for AT III			
Sens/wavelength	Low sens /405			
Read time	70 sec			
Select Tests	Name Tests	Special	↑	↓

Detection settings : Analysis range 20 - 50 sec

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