Manufactured by Hyphen BioMed.

## BIOPHEN HEPARIN LRT ON SYSMEX CA 7000

### 1. Preparation of BIOPHEN Heparin (LRT) reagent, (Reference A221011)

Both reagents, chromogenic substrate SXa-11 and Factor Xa, included in Biophen Heparin LRT kit, are ready to use, for the assay of UFH/LMWH (see chapter II) or Orgaran® (see chapter III). A calibration curve and a specific configuration will be use for each type of antithrombotic substance.

<b>R1 : SXa-11 Substrate</b> (liquid, ready to use).
<b>R2: Factor Xa</b> (liquid, ready to use).
Diluent (not provided): physiological saline.

- **Recommendation for stabilization of the reagents:**

Before use, let the reagents to stabilize for 30 minutes at room temperature (18-25°C), while shaking from time to time.

**It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate board before any use.**

- Stability data indicated on the insert are validated on vials kept closed, and provided any evaporation or contamination of reagents is avoided. They must be verified and adjusted if required according to the actual laboratory working conditions

**Storage of reagents:** Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C, and of well homogenising and strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

If the reagents are kept on the automate board, take care and use reducers to limit as much as possible any evaporation of the reagents.

- **Homogeneize the reagents before each use.**
- **Any reagent of biological origin must be handled with all the required cautions, as being potentially infectious.**
- **Do not interchange the reagents from different lots.**

## 2. Determination of Heparins UFH and/or LMWH.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Heparin Calibrator (réf A222001)	1 ml of distilled water (*)	7 days at 2-8°C 48 hours at room T° Do not freeze	30 minutes <i>on board</i> before any use (**)
Quality controls Biophen UFH Control (réf A223101) Biophen LMWH Control (réf A223001)	1 ml of distilled water (*)	7 days at 2-8°C 48 hours at room T° Do not freeze	30 minutes <i>on board</i> before any use (**)

## 3. Determination of Orgaran:

The determination of Orgaran requires its own configuration of the CA 7000. This one differs from the configuration used for the *determination* of UFH/LMWH only by the denomination of calibrators and controls.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Orgaran® Calibrator (réf A222201)	1 ml of distilled water (*)	7 days at room T° 48 hours at room T°	30 minutes <i>on board</i> before any use (**)
Quality controls Biophen Orgaran® Control (réf A223501)	1 ml of distilled water (*)	7 days at room T° 48 hours at room T°	30 minutes <i>on board</i> before any use (**)

### Reconstitution:

(\*) After reconstitution of calibrators or controls with distilled water, let them to stabilize for 30 minutes at room temperature. It is better to reconstitute calibrators the very day of calibration.

### Conservation of reagents:

(\*\*) Take care of strictly respecting the 30 minutes temperature stabilization time for *calibrators* and *controls* at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C. Homogenize before each use.

**Nota:** For lyophilized calibrators and controls, following reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature. It is recommended to run the calibration curve with a freshly reconstituted calibrator. It is necessary to let the reagent temperature to stabilize for at least 30 minutes onto the automate before any use. Take care avoiding any contamination or evaporation of the reagents. Stability can be adjusted according to the exact use conditions.

**Homogenize before each use.**

**Do not freeze calibrators and quality control plasmas.**

**A new calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).**

**Performances may present slight variations according to the instrument used. Validate the expected values in the laboratory working conditions.**

**Check the validity of the series by including quality control plasmas at different levels in each one.**

**Quality controls must be run regularly, and for each new batch of reagents, after an important maintenance of the instrument, or if measured values are not in compliance with the one expected for the method.**

#### 4. Results:

- The calibration curve is of the Log (absorbance) - Lin (concentration) type.
- The values obtained for the patients and controls are directly calculated from the calibration curve.
- The results are expressed in IU/ml.
- When Heparin or Orgaran concentrations are out of the working range, assayed plasma must be diluted in normal plasma, appropriately prepared and platelet poor, in order to keep a sufficient concentration of AT III.
- In presence of low AT III concentrations, as it can be the case in young children, an exogenous source of AT III is necessary, in order to correctly measure the heparin concentration.

**Nota:** Biophen Heparin LRT reagents are developed for measuring homogeneously Unfractionated heparin (UFH) and Low Molecular Weight Heparin (LMWH), using the same calibration curve.

The superimposition conditions are susceptible to slightly vary according to the combination of reagents and lots used, and the technical characteristics and specific adjustment of each apparatus.

Therefore, this superimposition has to be verified and validated in the exact laboratory working conditions and for each machine, by establishing the calibration curve using Biophen Heparin calibrators (LMWH, ref A222001), and by assaying the homogeneous system of associated Biophen UFH calibrators (ref A222301), UFH control plasmas (ref A223101), and LMWH control plasmas at high and low levels (ref A223001, A223701).

The linearity ( $r^2$ ) is expected to be  $\geq 0.98$  for the calibration curve.

Controls must be measured in compliance.

Measured values for calibrators are acceptable when:

Measured value = target value  $\pm 0.10$  IU/ml for levels  $\leq 0.50$  IU/ml (better  $\pm 0.05$  IU/ml)

Measured value = target value  $\pm 0.15$  IU/ml for levels  $> 0.5$  IU/ml (better 90-110% of the target value, or  $\pm 0.10$  IU/ml for levels between 0.5 and 1.0 IU/ml)

Should the superimposition not be obtained, a separate calibration curve must be used for each type of heparin, and validated by assaying the homogeneous quality control system (controls measured in compliance and close to the assigned target value).

**The calibration curve is validated when linearity ( $r^2$ ), as well as measured control values, are in compliance.**

**A new calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).**

**Performances may present slight variations according to the instrument used. Validate the expected values in the laboratory working conditions. Performances, as well as values for each new lot of quality controls used, must then be confirmed (and adjusted if necessary) in the laboratory working conditions.**

**Nota:** The stability at 2-8°C and at room temperature (18-25°C) claimed are validated on reconstituted vials kept closed, protected from any evaporation or contamination.

Stabilities of reagents on each automate must be controlled, adjusted and validated under each exact laboratory working conditions.

## Instrument setting for Sysmex CA 7000

### Create the reagents

Go to Analysis Configuration

Then Reagent information

Create the reagents

**Example of creation for Biophen Heparin LRT reference A221011:**

ID	NAME	LOT	Code Bar	Expir	Limit Init Vol.	Limit Alert	Intern Limit	Vial
1011	Sub HepLRT	00000-1*		00/00*	7.5 ml	1.02		GW5
1012	FXa Hep LRT	00000-2*		00/00*	7.5 ml	1.02		GW5
1008	Eau Phy	0		0	5 ml	1.02		GW5
1009	Cal 1....	4000*		00/00*	1 ml	0.2		Cup2

\*user defined

**ID** can be created with 4 numbers.

**Name:** Sub Hep LRT or FXa Hep LRT or eau Phy (NaCl must be in vial of 5 ml as the other reagents).

**Limit Alert:** 1.02 ml for the vial containing 5 ml and more.

**Vial** = GW5 for the vial of 5 ml or more and Cup2 for Standards and controls (vial of 1 ml).

**Go to « Protocol test » and write the following informations:**

Click on the window Set up software for the “manager program” and create the program according to: Chromogen for HEPARIN

Change the data with the data indicated on the table below

Use the program and create the name for the reagents, Calibration plasmas and control plasmas

**PROTOCOL FOR 333 TESTS:**

System		Ready		
				Emergency
<b>Analysis protocol</b>				
Parameter : UFHLMWHLRT      Para Code ...				
Sample Vol.			18**	µl
Vol. diluent	NaCl		36**	µl
Washing	Without			Detection method
Second dilution	No	0µl		1
Diluent Vol.	Without	0µl		7
Washing	Without			8
<b>Deficient Plasma</b>	Without	0 µl	0 sec	9
Vol Diluent	No Util.	No	0 µl	4
Washing	Without	x0/Off/Without x0		5
<b>First Reagent</b>	SUBhepLRT 90µl** 12 sec. Norm			6
Diluent Vol	Needle R A2	No	0 µl	1
Washing	Without	x0/	Clean I x 1	2
<b>Second Reagent</b>	XaHepLRT 90µl** 252 sec. Norm			3
Diluent Vol	Needle B1	No	0µl	0
Washing	Without	x0/	Clean I x3	Enter
<b>Third Reagent</b>	Without	0µl	0 sec	
Diluent Vol	No util	No	0µl	
Washing	Without	x0/Off/Without x0		
Detector (*)	Chromogen For AT III			
Wavelength (*)	Bas sens /405nm		Inc	
Detection setting (*)	20 sec		60 sec	
Select Tests	Test Name	Special	↑	↓
				Return

(\*) For CA7000 ⇒ This information is mentioned at the beginning of the table

**PROTOCOL FOR 375 TESTS: (\*\*) 16µl sample Vol, 32 µl Diluent Vol, 80µl SubLRT, 80µl XaLRT**

**Go to « Analysis Configuration »:**

Click on “Calcul Parameter” and complete the following table:

System		Ready analysis	HC GP		
			Emergency		

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