Manufactured By: HYPHEN BioMed

Biophen AT (#A221102/A221105) on ACL7000 (research software)

NOTE : reconstitution, preparation, calibration and controls, expression of results: usual cautions and realization according to the specific technical insert D.750.02/BI/1102 or /1105.

Reagents are reconstituted as per the insert, well equilibrated, and preheated at 37°C

PHASE	CONDITION	OPERATION	PARAMETERS
1 st loading step	YES	Reagent loading (pos 18)	0µl
		Washing after reagent loading (pos 18)	0
		Diluent (pos DIL)	0µL
		Specimen (working dilution 1 : 20 in physiological saline)	40µl
		Reagent (nacelle 2) Reagent R1 FXa	100µL
		Washing between specimen	2
		Optical ref. (cycle abs.) Diluent (pos DIL) Reagent (nacelle 2)	80µL 80µL
		Washing (at the end of loading)	2
		Incubation time (T1)	60 s
		2 nd loading step	YES
Specimen	0µL		
Washing between specimen	0		
Reagent (nacelle 3) Reagent R2 substrate	100µL		
Washing after reagent loading	2		
Mix time (T2)	60 s		
Activation time (T3)	0 s		
3 rd loading step	NO	Reagent (pos 17)	0µL
		Washing after reagent loading	0
		Mix time (T4)	0 s
Washing	NO	Reagent (pos 1)	0µL
		Reagent (pos 16)	0µL
		Cycles number	1
		Washing at the end of loading	1
Levels detection	NO		
Acquisition parameters	YES	Ramp	yes
		Inter-ramp interval (T6)	1s
		Delay before acquisition (T7)	5s
		Acquisition time (T8)	65s
		Canal used	405nm
		Rotation speed	1200 rpm

Note: The stability at 2-8°C and at room temperature (18-25°C) claimed are validated on reconstituted vials kept closed, protected from any evaporation or contamination. Stabilities of reagents on each automate must be controlled, adjusted and validated under each exact laboratory working conditions.

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