Manufactured By: HYPHEN BioMed

# Biophen AT on AMAX DESTINY

## 1. Reconstitution of Biophen AT Reagent (Ref A221102-A221105) :

Name	Reconstitution Biophen AT 5	Reconstitution Biophen AT 2.5	Stability***	Stabilization in T°
R1 AT (FXa)	5.00 ml of R3*	2.50ml of R3*	Refer insert***	30 minutes on board before use (**)
R2 AT (Substrate)	5.00 ml of Distilled water*	2.50ml of Distilled water*	Refer insert***	30 minutes on board before use (**)
R3 AT (Buffer)	Ready to use (5.0 ml)	Ready to use (2.50 ml)	Refer insert***	30 minutes on board before use (**)

\*The stability data claimed on the insert were obtained on reconstituted vials, kept closed, protected from and provided any contamination or evaporation is avoided. Stability must be controlled, and can be adjusted and validated if required, according to the exact use conditions for each laboratory.

**Reconstitution:** (\*) After reconstitution, let the reagent to stabilize for 30 minutes at room temperature (18-25°C).

**Stabilization of reagents:** (\*\*) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate board before any use.

**Storage of reagents:** Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C, and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

If the reagents are kept on the automate board, take care and use reducers to limit as much as possible any evaporation of the reagents.

**Homogenize the reagents before each use.**

Any reagent of biological origin must be handled with all the required cautions, as being potentially infectious.

**Do not interchange the reagents from different lots.**

### Reagents required but not provided:

- System Saline (9g/l NaCl; physiological saline)
- Destiny Pro Wash (DPW)
- Plasma calibrator and quality control titrated for AT (eg A222101,A223201,A223301)

## 2.Preparation of the calibrator, quality controls and samples to test :

### I. Reconstitution of the Calibrator for the calibration curve:

Name	Reconstitution	Stability	Stabilization in T°
<b>Biophen plasma calibrator</b> Ref # A222101	1 ml distilled water	Refer insert***	<b>30 minutes on board</b> before use (**)

- The exact AT concentration ( %) of the plasma calibrator is already indicated on the specific flyer for the lot used. The Amax Destiny machine realizes the dilutions in physiological saline to perform the standard calibration curve (following the instructions programmed by the user)  
The results' curve is given in AT rate (%) = f (mEquivalent/min).

**Note: For more information on the unit, see the Amax Destiny manual on chapter 10.**

### II. Reconstitution of controls

Name	Reconstitution	Stability	Stabilization in T°
<b>Biophen Normal Control</b> Ref # A223201	1 ml distilled water	Refer insert***	<b>30 minutes on board</b> before use (**)
<b>Biophen Abnormal Control</b> Ref # A223301	1 ml distilled water	Refer insert***	<b>30 minutes on board</b> before use (**)

Nota : For lyophilized calibrators and controls, following reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature. It is recommended to run the calibration curve with a freshly reconstituted calibrator. It is necessary to let the reagent temperature to stabilize for at least 30 minutes onto the automate before any use. Take care avoiding any contamination or evaporation of the reagents. Stability can be adjusted according to the exact use conditions.

**Homogeneize before each use.**

**Do not freeze calibrators and quality control plasmas.**

**Quality controls must be run regularly, and for each new batch of reagents, after an important maintenance of the instrument, or if measured values are not in compliance with the one expected for the method.**

### III. Samples

Samples, as controls, are loaded undiluted on the machine, and the dilution is directly managed by the automate according to the programming instructions.

### **3) Results :**

- The calibration curve is of the Lin-Lin type.
- The values obtained for patients and controls are directly calculated from the calibration curve (when the standard dilution is used for the test).
- The value of the standard calibration curve, the samples and controls is rendered in mEquivalent and the results are expressed in **AT concentration (%)**.

The calibration curve is validated when linearity (r2), as well as measured control values, are in compliance.

A new calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).

Performances may present variations according to the instrument used, and to detection sensitivity adjustment. Validate the expected values in the laboratory working conditions. Performances, as well as values for each new lot of quality controls used, must then be confirmed (and adjusted if necessary) in the laboratory working conditions.

### **4) Operations step by step on Amax Destiny (software version 1.1.5) Test Set up:**

**Consult the Amax Destiny user manual and the reagent insert.**

Before use, do the maintenance operation (cleaning, empty waste bin, control the liquid jet).

1) In the general menu enter in the **programming icon** « % R Sec µl g/l »

Enter the instructions in the various menus **Title – Reagents Procedure - Measure - Calibration - Curves** (See items into the programming part next page). Bind parameters into the tree

2) In the general menu enter in the **reagent icon** « symbolized by vials »

In the scroll text menu enter the reagents and all the necessary auxiliary reagents to perform the test. Load the reagents on the Amax

3) In the general menu enter in the **sample icon** « symbolized by a board »

Load calibrators, controls and samples on the machine and run the tests.

**Note:** The values measured for controls must be in compliance with expected results for the laboratory. As for samples, a complementary study must be conducted for each unexpected result.

---

## PROGRAMMING MENU

---

-----  
**First enter the codes chosen by the user (example AT) and the title (example Biophen AT)**

« Save » allows to save parameters in the little screen at the bottom of the left side.  
-----

**Enter the reagent area to define the parameters of the reagents bound**

-----//R1 AT//-----

<b>Name:</b> R1 AT	<b>Lot #:</b> User Defined		
<b>Producer:</b> HBM	<b>Expires:</b> User Defined		
<b>Barcode:</b>	<b>Stability:</b> 7 days ***		
<b>Type:</b> REAGENT	<b>Wash:</b> 4	<b>Requires Stirring:</b> No	
<b>Min. Volume:</b> User defined	<b>Max. Volume:</b> User defined	<b>Diameter:</b> User defined	<b>Bottom:</b> User Defined

« Save » allows to save parameters in the little screen at the bottom of the left side.

-----//R2 AT//-----

<b>Name:</b> R2 AT	<b>Lot #:</b> User Defined		
<b>Producer:</b> HBM	<b>Expires:</b> User Defined		
<b>Barcode:</b>	<b>Stability:</b> 7 days ***		
<b>Type:</b> starter	<b>Wash:</b> 4	<b>Requires Stirring:</b> No	
<b>Min. Volume:</b> User defined	<b>Max. Volume:</b> User defined	<b>Diameter:</b> User defined	<b>Bottom:</b> User Defined

« Save » allows to save parameters in the little screen at the bottom of the left side.

-----// System Saline//-----

<b>Name:</b> System Saline	<b>Lot #:</b> User Defined		
<b>Producer:</b> User defined	<b>Expires :</b> User Defined		
<b>Barcode:</b>	<b>Stability :</b> User Defined		
<b>Type:</b> System Buffer	<b>Wash:</b> 2	<b>Requires Stirring:</b> No	
<b>Min. Volume:</b> User defined	<b>Max. Volume:</b> User defined	<b>Diameter:</b> User defined	<b>Bottom:</b> User Defined

« Save » allows to save parameters in the little screen at the bottom of the left side.

-----// DPW//-----

<b>Name:</b> DPW	<b>Lot #:</b> User Defined		
<b>Producer:</b> Trinity	<b>Expires:</b> User Defined		
<b>Barcode:</b>	<b>Stability :</b> User Defined		
<b>Type:</b> Cleaner	<b>Wash:</b> 8	<b>Requires Stirring:</b> No	
<b>Min. Volume:</b> User defined	<b>Max. Volume:</b> User defined	<b>Diameter:</b> User defined	<b>Bottom:</b> User Defined

« Save » allows to save parameters in the little screen at the bottom of the left side.

-----  
**Enter the procedure area to define the Biophen AT test**  
 -----

-----// AT Procedure//-----

**Name:** AT procedure

Name : AT

Chromogenical

Dispense 50 µl of sample prediluted 1/40 with System Saline Add 100 µl of R1 AT Decontaminate with DPW Incubate 240 seconds Start Measure adding 100 µl R2 AT				
Step	Action	Val. I	Reagent/Buffer	Val. II
1	Predilute	50	System Saline	40
2	Add reagent	100	R1 AT	
3	Decontaminate		DPW	
4	Incubate	240		
5	Start Measure	100	R2 AT	

**SAVE .The AT procedure appears in the list of the available procedures.**

-----  
**Enter the measure area to define the measuring parameters.**  
 -----

-----// AT Measure//-----

<b>Name:</b> AT	<b>Duplicate Enable:</b> User defined	10 %dev
<b>Measure Mode:</b> Chromogenical	<b>Baseline Minimum:</b> NA	mEq
<b>Lag Time:</b> 10	<b>Baseline Maximum:</b> NA	mEq
<b>Timeout:</b> 60	<b>Linearity Check Enable:</b> no	
<b>Extended time out:</b>	<b>Negative Slope Enable:</b> NA	
<b>Activated:</b> no NA sec		

Save AT measure definition. It will appear in the list of the available measures.

-----  
**Enter the calibrator definition area to define calibrators (vial by vial)**  
 -----

-----// AT calib //-----

**Name:** AT

**code:** automatic

**Producer :** user defined

**lot#** user defined

Assigned and measured values:

**Exp:** user defined

Assay	dil	Date/time	rate	Measure	edited
AT	0	calibrator lot /calibrator number	***		

\*\*\*user defined.

Save the calibrator definition. It will appear in the list of the available calibrators.

### Enter in the curve definition area

-----// AT curve //-----

**Name:** AT

**curve:** lin/lin

**Assay:** AT

**Producer** user defined

**Unit %**

Id/lot	dil	conc	read
Calibrator lot	1:1		
Calibrator lot	1:2		
Calibrator lot	1:4.5		
Calibrator lot	1:7.5		
Calibrator lot	1:98		

**Enter five times the** above defined calibrator using the scroll text menu, and the five required dilutions to perform the standard calibration curve.

The software indicates the corresponding rates **when the test is run**.

The machine then performs the dilutions

The measured values in mExtinction appear in the column "read".

Save the curve definition. It will appear in the list of available calibration curves.

### Programming menu

In the higher left screen click on « system » to suppress the arborescence.


Click on the flash icon.

### Push on the button title

In the lower left screen the various titles appear.

### Click twice on the title AT

The various elements are bound in the following order: AT title ; AT curve, AT calibration , AT measure, AT procedure, AT Reagent .

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-----//R1 FXa//-----

<b>Name:</b> R1 AT	<b>Lot #:</b> User Defined		
<b>Producer:</b> HBM	<b>Expires:</b> User Defined		
<b>Barcode:</b>	<b>Stability:</b> 7 Days ***		
<b>Type:</b> REAGENT	<b>Wash:</b> 4	<b>Requires Stirring:</b> No	
<b>Min. Volume:</b> User defined	<b>Max. Volume:</b> User defined	<b>Diameter:</b> User defined	<b>Bottom:</b> UserDefined

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-----//R2 AT//-----

<b>Name:</b> R2 AT	<b>Lot #:</b> User Defined		
<b>Producer:</b> HBM	<b>Expires:</b> User Defined		
<b>Barcode:</b>	<b>Stability:</b> 7days ***		
<b>Type</b> starter	<b>Wash:</b> 4	<b>Requires Stirring:</b> No	
<b>Min. Volume:</b> User defined	<b>Max. Volume:</b> User defined	<b>Diameter:</b> User defined	<b>Bottom:</b> UserDefined

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-----// System Saline//-----

<b>Name:</b> System Saline	<b>Lot #:</b> User Defined		
<b>Producer:</b> User defined	<b>Expires:</b> User Defined		
<b>Barcode:</b>	<b>Stability:</b> User Defined		
<b>Type:</b> System Buffer	<b>Wash:</b> 2	<b>Requires Stirring:</b> No	
<b>Min. Volume:</b> User defined	<b>Max. Volume:</b> User defined	<b>Diameter:</b> User defined	<b>Bottom:</b> UserDefined

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-----// DPW//-----

<b>Name:</b> DPW	<b>Lot #:</b> User Defined		
<b>Producer:</b> Trinity	<b>Expires:</b> User Defined		
<b>Barcode:</b>	<b>Stability:</b> User Defined		
<b>Type:</b> Cleaner	<b>Wash:</b> 8	<b>Requires Stirring:</b> No	
<b>Min. Volume:</b> User defined	<b>Max. Volume:</b> User defined	<b>Diameter:</b> User defined	<b>Bottom:</b> UserDefined

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**Name:** AT procedure

Nom : AT

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Dispense 50 µl of sample prediluted 1/40 with System Saline Add 100 µl of R1 AT Decontaminate with DPW Incubate 240 seconds Start Measure adding 100 µl R2 AT				
Step	Action	Val. I	Reagent/Buffer	Val. II
1	Sample	50	Systeme Saline	40
2	Add reagent	100	R1 AT	
3	Decontaminate		DPW	
4	Incubate	240		
5	Start Measure	100	R2 AT	

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