Manufactured By: HYPHEN BioMed

BIOPHEN ANTITHROMBIN ON BCS

Adaptation of BIOPHEN ANTITHROMBIN ON BCS

1. Reconstitution of BIOPHEN ANTITHROMBIN 2.5 (Ref A221102) reagent:

Chromogenic determination of the Antithrombin.

	NAME	Reconstitution	Stability	Stabilization in T°
R1	Factor Xa	2.5 ml TRIS Buffer (*)	3 months at 2-8°C (**) 7 days at room T° (**) Do not freeze	30 mn before any use (***)
R2	SXa-11 Substrate	2.5 ml of distilled water (*)	3 months at 2-8°C (**) 7 days at room T° (**) Do not freeze	30 mn before any use (***)
Diluent	Physiological Saline		24 hours <i>on BCS board</i> (**)	

2. Reconstitution of BIOPHEN ANTITHROMBIN 5 (Ref A221105) reagent:

	NAME	Reconstitution	Stability	Stabilization in T°
R1	Factor Xa	5 ml TRIS Buffer (*)	3 months at 2-8°C (**) 7 days at room T° (**) Do not freeze	30 mn before any use (***)
R2	SXa-11 Substrate	5 ml of distilled water (*)	3 months at 2-8°C (**) 7 days at room T° (**) Do not freeze	30 mn before any use (***)
Diluent	Physiological Saline		24 hours <i>on BCS board</i> (**)	

Reconstitution:

(*) Following reconstitution with distilled water let the R1 and R2 reagents to stabilize for 30 minutes at room temperature.

Conservation of reagents Stability:

()** Provided any contamination or evaporation is avoided. Stability can be adjusted according to the exact use conditions.

Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

If the substrate and Factor Xa are kept 24h/24h on automate board, take care and use the chimneys to avoid or limit as much as possible any evaporation.

Stabilization of reagents:

(*)** It is necessary to let the substrate (R2) and the Xa Factor (R1) temperature to stabilize for at least 30 minutes on the automate before any use

Foot-note: Do not interchange the reagents from different lots.

3. Determination of ANTITHROMBIN:

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Plasma Calibrator (ref A222101)	1 ml of distilled water (*)	24 hours at 2-8°C (**) 8 hours at room T° (**)	30 minutes on BCS board before any use (***)
Quality controls Biophen Normal Control (ref A223201) Biophen Abnormal Control (ref A223301)	1 ml of distilled water (*)	24 hours at 2-8°C (**) 8 hours at room T° (**)	30 minutes on BCS board before any use (***)

Reconstitution:

(*) Following reconstitution of calibrators or controls with distilled water let them to stabilize for 30 minutes at room temperature.

Foot-note: A calibration curve must be carried out for each new batch of reagents

Conservation of reagents :

()** Provided any contamination or evaporation is avoided. Stability can be adjusted according to the exact use conditions.

Stabilization of reagents

(*)** Take care of strictly respecting the 30 minutes temperature stabilization time for *calibrators* and *controls* at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C. Homogenize before each use.

Foot-note : Do not freeze calibrators or controls.

Quality controls must be run regularly, and for each new batch of reagents, after an important maintenance of the instrument, or if measured values are not in compliance with the one expected.

4. Results:

- The calibration curve is of the Lin (absorbance) - Lin (concentration) type.
- The values obtained for the patients and controls are directly calculated from the calibration curve.
- The results are expressed in % activity.

5. Programming the BCS analyzer:

Creation of reagents and diluent files

Return to definition of test. Go in reagents without data of lot. Create both reagents files, and also the Diluent file.

Name: Substrate BIOPHEN AT			
Abbreviation: SXA11			
clean when the reagent is modified			high
<input checked="" type="checkbox"/> New calibration with any new lot of reagents			
Authorized lines			
from		1	to 14
Stability (*)			
< 15C			> = 15C
Mix by jet			
Intensity		1	
Frequency		No agitation	
Reference numbers			
	Number		
	add	remove	replace

Name: Factor Xa Biophen AT			
Abbreviation: FXAAT			
clean when the reagent is modified			special
<input checked="" type="checkbox"/> New calibration with any new lot of reagents			
Authorized lines			
from		1	to 14
Stability (*)			
< 15C			> = 15C
Mix by jet			
Intensity		1	
Frequency		No agitation	
Reference numbers			
	Number		
	add	remove	replace

Name: Physiological saline			
Abbreviation: PHYSAL			
clean when the reagent is modified			Normal/low
New calibration with any new lot of reagents			
Authorized lines			
from		1	to 14
Stability (*)			
< 15C			> = 15C
Mix by jet			
Intensity		1	
Frequency		No agitation	
Reference numbers			
	Number		
	add	remove	replace

Create in the same way calibration plasmas and control plasmas files.

Name: Biophen plasma Calibrator			
Abbreviation: ATCAL			
clean when the reagent is modified			Normal/low
<input checked="" type="checkbox"/> New calibration with any new lot of reagents			
Authorized lines			
from		1	to 14
Stability (*)			
< 15C			> = 15C
Mix by jet			
Intensity		1	
Frequency		No agitation	
Reference numbers			
	Number		
	add	remove	replace

Name: Biophen plasma Control			
Abbreviation: AT Cont			
clean when the reagent is modified			Normal/low
<input checked="" type="checkbox"/> New calibration with any new lot of reagents			
Authorized lines			
from		1	to 14
Stability (*)			
< 15C			> = 15C
Mix by jet			
Intensity		1	
Frequency		No agitation	
Reference numbers			
	Number		
	add	remove	replace

(*) To be filled by the user with stabilities data provided into the kit insert

Creation of the procedure

Once the various reagents files are created, the procedure file must be created. In order to do it, go to **definition of tests**, procedure and type the following parameters.

General Information

Former Change			
Test Procedure Number			
Test Procedure Name	BIOPHEN AT		
Measurement time stopped after	100 sec	Or	mA
	0		
Speed	normal		
Stirring	normal		
Wavelength	405 nm		
Primary absorbance range	131-160 mA		
Duplicate n# for samples and controls (*)		CV Authorized (*)	0.0 %
Duplicate n# for calibration (*)		CV Authorized (*)	0.0 %
Duplicate n# of summary of raw data	Arithmetic mean		
Dilution Factor	1 : 4	(1 = no dilution)	
Dilution Buffer	Physiological saline		

(*) To be filled by the user

Evaluate and Control Method			
Delta A per minute			
from	10	Sec	60
			Sec
Correlation Factor		0.9	Absolute

Pipetting Sequence

Test Procedure Number

Test Procedure Name

BIOPHEN AT

Sampling cycles

Number of Cycle	Arm of transmission	Final Washing	Action of Rotors/transfer	Time Range [sec]	
				Min	Max
1	Arm	normal	No mix. in priority	0	0
2	Reagent Arm	normal	Incubation	50	60
3	Reagent Arm	normal	No mix/measure	0	0
4	Reagent Arm	normal	No mix/measure	0	0

Num de transm.	Mean/pipetting goal	Reagent	Speed	Volume (µl)
1	Take air		Slow	20
1	Take the Reagent	PHYSAL	Slow	45
1	Take the sample		Slow	5
1	In the external cup		Slow	50
2	Take air		Slow	20
2	Take the Reagent	R1AT FactorXa	Slow	100
2	In the central cup		Slow	100
3	Take air		Slow	20
3	Take the Reagent	R2 AT Substrate	Slow	100
3	In the central cup		Slow	100
4	Take the Reagent	SCS Cleaner	Slow	0

Test Definition

General Information				
Former Change		<input type="text"/>		
Name		<input type="text" value="BIOPHEN AT"/>		°Test N°
Result Unit		PC ID t		
Behring	% dN	Conversion factor		Abbreviation
User	% dN	1%dN	1.0	1%dN
PC	% dN	1 %dN	1.0	1%dN
Place in the list				<input type="text" value="AT"/>
Digit number for results				<input type="text" value="1000"/>

Test detailed procedures

Procedure Name test	Evaluation Method
<input type="text" value="BIOPHEN AT"/>	<input type="text" value="Delta A per minute"/>

Calculate Formula of raw data

Formula Type	Minimum	<input type="text"/>
User own Formula	<input type="text"/>	
Unit	mE/min	<input type="text"/>

Detailed test evaluation

Usual Values		Determination of result		Standard Concentrations
		<input type="text" value="70"/> % dN	at <input type="text" value="140"/> %dN	
Test used for Calibration		<input type="text"/>		
Evaluation	<input type="text" value="Regression lin/lin"/>			
Reference Curve	<input type="text" value="Calibration curve"/>			
Minimum raw value	<input type="text" value="0.0"/>			
Upper Extrapolation	<input type="text" value="1"/>	x maxi concentration		
Lower Extrapolation	<input type="text" value="1"/>	x mini concentration		
Calibrator	Standard plasma calibrator			
Authorized SD respectively to the reference curve	<input type="text" value="0"/>	%		
Maximal number repeats	<input type="text" value="0"/>			

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