Manufactured By: HYPHEN BioMed

BIOPHEN ANTITHROMBIN 5 ON STA-R

Adaptation of BIOPHEN ANTITHROMBIN 5 ON STA-R

1. Reconstitution of BIOPHEN ANTITHROMBIN 5 (Ref A221105) reagents:

Chromogenic determination of the Antithrombin.

	NAME	Reconstitution	Stability	Stabilization in T°
R1	Factor Xa	5 ml TRIS Buffer (*)	3 months at 2-8°C 7 days at room T° Do not freeze	30 mn before any use (**)
R2	SXa-11 Substrate	5 ml of distilled water (*)	3 months at 2-8°C 7 days at room T° Do not freeze	30 mn before any use (**)
Diluent	Physiological Saline		24 hours <i>on STA-R board</i>	

Reconstitution:

(*) Following reconstitution with distilled water, let the R1 and R2 reagents to stabilize for 30 minutes at room temperature.

Conservation of reagents:

Take care of putting up the specific caps back on the bottles before storing them at 2°-8°C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Stabilization of reagents:

(**) It is necessary to let the substrate (R2) and the Xa Factor (R1) temperature to stabilize for at least 30 minutes on the automate before any use.

Foot-note: Do not interchange the reagents from different lots.

2. Determination of ANTITHROMBIN.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Plasma Calibrator (réf A222101)	1 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T°	30 minutes <i>on STA-R board</i> before any use (**)
Quality controls Biophen Normal Control (réf A223201) Biophen Abnormal Control (réf A223301)	1 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T°	30 minutes <i>on STA-R board</i> before any use (**)

Reconstitution:

(*) After reconstitution of calibrators or controls with distilled water, let them to stabilize for 30 minutes at room temperature.

Foot-note: A calibration curve must be carried out for each new batch of reagents.

Conservation of reagents:

(**) Take care of strictly respecting the 30 minutes temperature stabilization time for *calibrators* and *controls* at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C. Homogénéize before each use.

Foot-note: Do not freeze calibrators or controls.

3. Results:

- The calibration curve is of the Lin (absorbance) - Lin (concentration) type.
- The values obtained for the patients and controls are directly calculated from the calibration curve.
- The results are expressed in % activity.
The 100% activity is that of reference normal pooled citrated plasma.

4. Programming the STA-R analyzer:

Click on the icon **set up software** for the **manager program** and create the program according to:

TESTS						Configuration		
Method								
IDENTIFICATION								
Abbreviation		AT 5		Last up date				
Name		BIOPHEN AT5		Method		Chro-Kinetics OD /min		
SAMPLE			DILUENT					
Vol. μ l	Incu.	Dilution	Id.	Name		Stab. h	Continued	
100 μ l	0 sec	1/50	PS	PHYSAL		24		
REAGENTS								
	Id.	Name		Stab. h	Vol. μ l	Incu. sec	Prec	Continued
Ra	AT XA	AT XA		168	100	240	<input type="checkbox"/>	
Rb							<input type="checkbox"/>	
Rc							<input type="checkbox"/>	
Rd	AT SUB	AT SUB		168	100			

Enter data concerning washing

	Id.	Name	Stab. h	Vol. μ l	Incu. sec	Prec.	Vial. ml	Min Vol ml
Diluent	PS	PHYSAL	24				5	1
Ra	AT XA	AT XA	168	100	240	<input type="checkbox"/>	5	0.5
Rb						<input type="checkbox"/>		
Rc						<input type="checkbox"/>		
Rd	AT SUB	AT SUB	168	100			5	0.5
Washing								
	Ra		Rb		Rc		Rd	
Before	No						No	
After	Special						Intensive	
Wash	12227							
Name	STA-DESORB U							
Stab. H	120							
Vial ml	15							
Vol in. ml	0.9							

Click on the icon **Result software** for the **manager program** and create the program according to:

TESTS		Configuration							
		Result							
METHOD Chro KINETICS OD/mn		RESULT							
Linearity Min.	0.95	Primary Units	%						
First Point	10	Corrector	1.00						
Second Point	60	Determination	<input checked="" type="radio"/> Single <input type="radio"/> Duplicate						
		Precision (%)	5.00						
<input checked="" type="checkbox"/> Rd Heating <input type="checkbox"/> Stirring		Redilution	Condition						
		1/	<						
		1/	>						
		<table border="1"> <thead> <tr> <th colspan="2">Validation (%)</th> </tr> <tr> <th>Min.</th> <th>Max.</th> </tr> </thead> <tbody> <tr> <td>0</td> <td>140</td> </tr> </tbody> </table>		Validation (%)		Min.	Max.	0	140
Validation (%)									
Min.	Max.								
0	140								

Write calibration data.

TESTS		Configuration	
		Calibration	
MODE :	Linear	Determination	
		<input checked="" type="checkbox"/> Duplicate <input type="checkbox"/> Single	
ASSAYS		Concentrations	
SCALE : <input checked="" type="checkbox"/> Linear <input type="checkbox"/> Log		<input checked="" type="checkbox"/> Linear <input type="checkbox"/> Inverse <input type="checkbox"/> Log	
	Id.	Key	Name
	Stab. H	Dilution	
Calibrator 1	ZERO		ZEROCAL
Calibrator 2	BPC		BPCAL
Calibrator 3	BPC		BPCAL
Calibrator 4	BPC		BPCAL
Calibrator 5			
Ctrl, Niv. 1			
Ctrl, Niv. 2			

Click on the icon **Printout/Transmission software** for the **manager program** and create the program according to:

TESTS		Configuration				
		Printout /Trans.				
PARAMETERS						
	Unit	Conver Factors.	Print	Transmission Test number	Usual Values Min. Max.	
Main	%		<input checked="" type="checkbox"/>	0	70	140
Aux1	OD/min		<input checked="" type="checkbox"/>			
Aux2			<input type="checkbox"/>			
Aux3			<input type="checkbox"/>			
Printout Limits						
Min.	<input type="text" value="0"/>	Max.	<input type="text" value="140"/>			

Show the last page and enter data concerning the quality control.

TESTS		Configuration				
		Q.C				
LEVEL1						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol. ml
Control	BNC		BN CONT	8	1	0.2
LEVEL 2						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control	BAC		BA CONT	8	1	0.2
LEVEL3						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control						
Period : <input type="text"/> hours <input type="text"/> <input type="checkbox"/> Tests <input type="checkbox"/> Vial						

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