Adaptation of BIOPHEN AT (LRT) ON STA-R

1. BIOPHEN AT (LRT)(Ref A221111) reagent:

Chromogenic determination of the Antithrombin.

Reagent of BIOPHEN AT (LRT) kit are liquid and ready to use.
BIOPHEN AT (LRT) kit must be stored at 2 - 8 °C, in their original packaging box.
They are then stable until the expiration date printed on the box.

	NAME	Packaging	Stability after opening	Stabilization in T°
R1	Factor Xa	4 x 11ml ready for use	3 months at 2-8°C (*) 7 days at room T° (*) Do not freeze	30 mn before any use (**)
R2	SXa-11 Substrate	4 x 4 ml ready for use	3 months at 2-8°C (*) 7 days at room T° (*) Do not freeze	30 mn before any use (**)
Diluent	Physiological Saline		24 hours <i>on board</i> STA	

(*) . If any contamination is avoided

Conservation of reagents:

Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Stabilization of reagents:

()** It is necessary to let the substrate R2 and the Xa Factor R1 temperature to stabilize for at least 30 minutes on the automate before any use.

Foot-note: Do not interchange the reagents from different lots.

2. Determination of ANTITHROMBIN.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Plasma Calibrator (ref A222101)	1 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T°	30 minutes <i>on STA board</i> before any use (**)
Quality controls Biophen Normal Control (ref A223201) Biophen Abnormal Control (ref A223301)	1 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T°	30 minutes <i>on STA board</i> before any use (**)

Reconstitution:

(*) Following reconstitution of calibrators or controls with distilled water, let them to stabilize for 30 minutes at room temperature.

Foot-note: A calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range.

Storage of reagents:

(**) Take care of strictly respecting the 30 minutes temperature stabilization time for *calibrators* and *controls* at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C. Homogenize before each use.

Foot-note: Do not freeze calibrators or controls.

3. Results:

- The calibration curve is of the Lin (absorbance) - Lin (concentration) type.
- The values obtained for the patients and controls are directly calculated from the calibration curve.
- The results are expressed in % activity.

4. Programming the STA-R analyzer:

Click on the icon **set up software** for the **manager program** and create the program according to:

TESTS		Configuration						
Method								
IDENTIFICATION								
Abbreviation		AT LRT		Last up date				
Name		BIOPHEN AT LRT		Method		Colo-Kinetics OD /min		
SAMPLE				DILUENT				
Vol. μ l	Incu.	Dilution	Id.	Name			Stab. h	Continued
25 μ l	0 sec	1/10	PS	PHYSAL			24	
REAGENTS								
	Id.	Name		Stab. h	Vol. μ l	Incu. sec	Prec	Continued
Ra	LR1	R1 LRT XA		168	200	90	<input type="checkbox"/>	
Rb							<input type="checkbox"/>	
Rc							<input type="checkbox"/>	
Rd	LR2	R2 LRT SUB		168	75			

Enter data concerning washing

	Id.	Name	Stab. h	Vol. μ l	Incu. sec	Prec.	Vial. ml	Min Vol ml
Diluent	PS	PHYSAL	24				5	1
Ra	LR1	R1 LRT XA	168	200	90	<input type="checkbox"/>	11	0.5
Rb						<input type="checkbox"/>		
Rc						<input type="checkbox"/>		
Rd	LR2	R2 LRT SUB	168	75			4	0.5
Washing								
	Ra		Rb		Rc		Rd	
Before	No						No	
After	Special						Intensive	
Wash	12227							
Name	STA-DESORB U							
Stab. H	120							
Vial ml	15							
Vol in. ml	0.9							

Click on the icon **Result software** for the **manager program** and create the program according to:

TESTS		Configuration			
		Result			
METHOD Chro KINETICS OD/mn		RESULT			
Linearity Min.	0.95	Primary Units	%		
First Point	10	Corrector	1.00		
Second Point	60	Determination	<input checked="" type="radio"/> Single		
			<input type="radio"/> Duplicate		
		Precision (%)	5.00	Validation (%)	
		Redilution	1/	Min.	Max.
<input checked="" type="checkbox"/> Rd Heating			<	0	120
<input type="checkbox"/> Stirring			>		
		Condition			

Write calibration data.

TESTS		Configuration			
		Calibration			
MODE :		Linear	Determination		
			<input checked="" type="checkbox"/> Duplicate		<input type="checkbox"/> Single
ASSAYS		Concentrations			
SCALE : <input checked="" type="checkbox"/> Linear <input type="checkbox"/> Log		<input checked="" type="checkbox"/> Linear <input type="checkbox"/> Inverse <input type="checkbox"/> Log			
	Id.	Key	Name	Stab. H	Dilution
Calibrator 1	LRTCAL	1/10	LRTCAL	8	1/10
Calibrator 2	LRTCAL	1/20	LRTCAL	8	1/20
Calibrator 3	LRTCAL	1/40	LRTCAL	8	1/40
Calibrator 4	LRTCAL	1/80	LRTCAL	8	1/80
Calibrator 5	ZEROCAL	1/1	ZEROCAL	8	1/1
Ctrl, Niv. 1					
Ctrl, Niv. 2					

Click on the icon **Printout/Transmission software** for the **manager program** and create the program according to:

TESTS		Configuration				
		Printout /Trans.				
PARAMETERS						
	Unit	Conver Factors.	Print	Transmission Test number	Usual Values Min. Max.	
Main	%		<input checked="" type="checkbox"/>	0	70	120
Aux1	OD/min		<input checked="" type="checkbox"/>			
Aux2			<input type="checkbox"/>			
Aux3			<input type="checkbox"/>			
Printout Limits						
Min.	<input type="text" value="0"/>	Max.	<input type="text" value="140"/>			

Show the last page and enter data concerning the quality control.

TESTS		Configuration				
		Q.C				
LEVEL1						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol. ml
Control	BNC		BN CONT	8*	1*	0.2*
LEVEL 2						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control	BAC		BA CONT	8*	1*	0.2*
LEVEL3						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control						
Period : <input type="text"/> hours <input type="text"/> <input type="checkbox"/> Tests <input type="checkbox"/> Vial						

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