Manufactured By: HYPHEN BioMed

Biophen PC on AMAX DESTINY

1. Reconstitution of Biophen PC Reagents (Ref A221202-A221205) :

Name	Reconstitution Bi. PC 5.0	Reconstitution Bi. PC 2.5	Stability***	Stabilization at RT
R1 PC (PROTAC)	5.0 ml of Distilled water*	2.5ml of Distilled water*	Refer insert***	30 minutes on board before use (**)
R2 PC (Substrate)	5.0 ml of Distilled water*	2.5ml of Distilled water*	Refer insert***	30 minutes on board before use (**)

*The stability data claimed on the insert were obtained on reconstituted vials, kept closed, protected from and provided any contamination or evaporation is avoided. Stability must be controlled, and can be adjusted and validated if required, according to the exact use conditions for each laboratory.

Reconstitution: (*) After reconstitution, let the reagent to stabilize for 30 minutes at room temperature (18-25°C).

Stabilization of reagents: (**) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate board before any use.

Storage of reagents: Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C, and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

If the reagents are kept on the automate board, take care and use reducers to limit as much as possible any evaporation of the reagents.

Homogenize the reagents before each use.

Any reagent of biological origin must be handled with all the required cautions, as being potentially infectious.

Do not interchange the reagents from different lots.

Reagents required but not provided:

- System Saline (9g/l NaCl; physiological saline)
- Destiny Pro Wash (DPW)
- Plasma calibrator and quality control titrated for PC (eg A222101,A223201,A223301)

2. Preparation of the calibrator, quality controls and samples to test :

I. Reconstitution of the Calibrator for the calibration curve:

Name	Reconstitution	Stability	Stabilization in T°
Biophen plasma calibrator Ref # A222101	1 ml distilled water	Refer insert***	30 minutes on board before use (**)

The exact PC concentration (%) of the plasma calibrator is already indicated on the specific flyer for the lot used. The Amax Destiny machine realizes the dilutions in physiological saline to perform the standard calibration curve (following the instructions programmed by the user)

The results' curve is given in PC rate (%) = f (mEquivalent/min).

Note: For more information on the unit, see the manual Amax Destiny on the chapter 10

I. Reconstitution of controls

Name	Reconstitution	Stability	Stabilization in T°
Biophen Normal Control Ref # A223201	1 ml distilled water	Refer insert***	30 minutes on board before use (**)
Biophen Abnormal Control Ref # A223301	1 ml distilled water	Refer insert***	30 minutes on board before use (**)

Nota : For lyophilized calibrators and controls, following reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature. It is recommended to run the calibration curve with a freshly reconstituted calibrator. It is necessary to let the reagent temperature to stabilize for at least 30 minutes onto the automate before any use. Take care avoiding any contamination or evaporation of the reagents. Stability can be adjusted according to the exact use conditions.

Homogenize before each use.

Do not freeze calibrators and quality control plasmas.

Quality controls must be run regularly, and for each new batch of reagents, after an important maintenance of the instrument, or if measured values are not in compliance with the one expected for the method.

III. Samples

Samples, as controls, are loaded undiluted on the machine, and the dilution, if required, is directly managed by the automate according to the programming instructions.

3) Results :

- The calibration curve is of the Lin-Lin type.
- The values obtained for patients and controls are directly calculated from the calibration curve (when the standard dilution is used for the test).
- The value of the standard calibration curve, the samples and controls is rendered in mEquivalent and the results are expressed in **PC concentration (%)**.

The calibration curve is validated when linearity (r^2), as well as measured control values, are in compliance.

A new calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).

Performances may present variations according to the instrument used, and to detection sensitivity adjustment. Validate the expected values in the laboratory working conditions. Performances, as well as values for each new lot of quality controls used, must then be confirmed (and adjusted if necessary) in the laboratory working conditions.

4) Operations step by step on Amax Destiny (software version 1.1.5) Test Set up:

Consult the Amax Destiny user manual and the reagent insert.

Before use, do the maintenance operation (cleaning, empty waste bin, control the liquid jet) if required.

1) In the general menu enter in the **programming icon** « % R Sec μ l g/l »

Enter the instructions in the various menus **Title – Reagents Procedure -Measure - Calibration - Curves** (See items into the programming part next page). Bind parameters into the tree

2) In the general menu enter in the **reagent icon** « symbolized by vials »

In the scroll text menu enter the reagents and all the necessary auxiliary reagents to perform the test. Load the reagents on the Amax.

3) In the general menu enter in the **sample icon** « symbolized by a board ».

Load calibrators, controls and samples on the machine and run the tests.

Note: The values measured for controls must be in compliance with expected results for the laboratory. As for samples, a complementary study must be conducted for each unexpected result.

PROGRAMMING MENU

First enter the codes chosen by the user (example PC) and the title (example Biophen PC)

« Save » allows to save parameters in the little screen at the bottom of the left side.

Enter the reagent area to define the parameters of the reagents bound

-----//R1 PC//-----

Name: R1 PC	Lot #: User Defined		
Producer: HBM	Expires: User Defined		
Barcode:	Stability: 3 days ***		
Type: REAGENT	Wash: 8	Requires Stirring: No	
Min. Volume: User defined	Max. Volume: User defined	Diameter: User defined	Bottom: User Defined

« Save » allows to save parameters in the little screen at the bottom of the left side.

-----//R2 PC//-----

Name: R2 PC	Lot #: User Defined		
Producer: HBM	Expires: User Defined		
Barcode:	Stability: 3 days ***		
Type starter	Wash: 4	Requires Stirring: No	
Min. Volume: User defined	Max. Volume: User defined	Diameter: User defined	Bottom: User Defined

« Save » allows to save parameters in the little screen at the bottom of the left side.

-----// System Saline//-----

Name: System Saline	Lot #: User Defined		
Producer: User defined	Expires : User Defined		
Barcode:	Stability : User Defined		
Type: System Buffer	Wash: 2	Requires Stirring: No	
Min. Volume: User defined	Max. Volume: User defined	Diameter: User defined	Bottom: User Defined

« Save » allows to save parameters in the little screen at the bottom of the left side.

-----// DPW//-----

Name: DPW	Lot #: User Defined		
Producer: Trinity	Expires: User Defined		
Barcode:	Stability: : User Defined		
Type: Cleaner	Wash: 8	Requires Stirring: No	
Min. Volume: User defined	Max. Volume: User defined	Diameter: User defined	Bottom: User Defined

« Save » allows to save parameters in the little screen at the bottom of the left side.

Enter the procedure area to define the Biophen PC test

-----// PC Procedure//-----

Name: PC procedure

Name : PC

Chromogenical

Dispense 20 µl of sample Add 80 µl of R1 PC Decontaminate with DPW Incubate 300 seconds Start Measure adding 80 µl R2 PC Decontaminate with DPW				
Step	Action	Val. I	Reagent/Buffer	Val. II
1	Add sample	20		
2	Add reagent	80	R1 PC	
3	Decontaminate		DPW	
4	Incubate	300		
5	Start Measure	80	R2 PC	
6	Decontaminate		DPW	

SAVE .The PC procedure appears in the list of the available procedures.

Enter the measure area to define the measuring parameters.

-----// PC Measure//-----

Name: PC	Duplicate Enable: User defined (10 %dev)
Measure Mode: Chromogenical	Baseline Minimum: NA mEq
Lag Time: 20	Baseline Maximum: NA mEq
Timeout: 140	Linearity Check Enable: Yes (15 %dev)
Extended time out:	Negative Slope Enable: NA
Activated: no NA sec	

Save PC measure definition. It will appear in the list of available measures

Enter the calibrator definition area to define calibrators (vial by vial)

-----// PC calib //-----

Name: PC

code: automatic

Producer : user defined

lot# user defined

Assigned and measured values:

Exp: user defined

Assay	dil	Date/time	rate	Measure	edited
PC	0	calibrator lot /calibrator number	***		

***user defined.

Save calibrator definition. It will appear in the list of available calibrators.

Enter the curve definition area

-----// **PC curve** //-----

Name: PC

curve: lin/lin

Assay: PC

Producer user defined

Unit %

Id/lot	dil	conc	read
Calibrator lot	1:1		
Calibrator lot	1:2		
Calibrator lot	1:4		
Calibrator lot	1:8		

Enter four times the above defined calibrator by the scroll text menu and the four required dilutions to perform the standard calibration curve.

The software indicates the corresponding rates **when the test is run.**

The machine then performs the dilutions

The measured values in mExtinction appear in the column read .

Save the curve definition. It will appear in the list of available calibration curves.

Programming menu

In the higher left screen click on « system » to suppress the arborescence


Click on the flash icon.

Push on the button title

In the lower left screen the various titles appear.

Click twice on the title PC

The various elements are bound in the following order: PC title ; PC curve, PC calibration, PC measure, PC procedure, PC Reagents .

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-----//R1 PC //-----

Name: R1 PC	Lot #: User Defined		
Producer: HBM	Expires: User Defined		
Barcode:	Stability: 3 Days ***		
Type: REAGENT	Wash: 8	Requires Stirring: No	
Min. Volume: User defined	Max. Volume: User defined	Diameter: User defined	Bottom: UserDefined

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-----//R2 PC //-----

Name: R2 PC	Lot #: User Defined		
Producer: HBM	Expires: User Defined		
Barcode:	Stability: 3days ***		
Type starter	Wash: 4	Requires Stirring: No	
Min. Volume: User defined	Max. Volume: User defined	Diameter: User defined	Bottom: UserDefined

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-----// System Saline//-----

Name: System Saline	Lot #: User Defined		
Producer: User defined	Expires: User Defined		
Barcode:	Stability: User Defined		
Type: System Buffer	Wash: 2	Requires Stirring: No	
Min. Volume: User defined	Max. Volume: User defined	Diameter: User defined	Bottom: UserDefined

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-----// DPW//-----

Name: DPW	Lot #: User Defined		
Producer: Trinity	Expires: User Defined		
Barcode:	Stability: User Defined		
Type: Cleaner	Wash: 8	Requires Stirring: No	
Min. Volume: User defined	Max. Volume: User defined	Diameter: User defined	Bottom: UserDefined

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Name: PC procedure

Nom : PC

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Dispense 20 µl of sample Add 80 µl of R1PC Decontaminate with DPW Incubate 300 seconds Start Measure adding 80 µl R2 PC Decontaminate with DPW				
Step	Action	Val. I	Reagent/Buffer	Val. II
1	Sample	20 µl		
2	Add reagent	80 µl	R1 PC	
3	Decontaminate		DPW	
4	Incubate	300 sec		
5	Start Measure	80 µl	R2 PC	
6	Decontaminate		DPW	

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