### **Adaptation of BIOPHEN Protein C 2.5 ON BCS**

#### **1. Reconstitution of BIOPHEN Protein C 2.5 (Ref A221202) reagents.**

Chromogenic determination of Protein C.

	<b>NAME</b>	<b>Reconstitution</b>	<b>Stability</b>	<b>Stabilization in T°</b>
R1	Protac	2.5 ml of distilled water (*)	3 months at 2-8°C 3 days at room T° Do not freeze	30 mn before any use (**)
R2	SAPC-21 Substrate	2.5 ml of distilled water (*)	3 months at 2-8°C 3 days at room T° Do not freeze	30 mn before any use (**)

#### Reconstitution:

(\*) Following reconstitution with distilled water, let the R1 and R2 reagents to stabilize for 30 minutes at room temperature.

#### Conservation of reagents:

Take care of putting up the specific caps back on the bottles before storing them at 2°-8°C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

#### Stabilization of reagents:

(\*\*) It is necessary to let the substrate (R2) and the Protac (R1) temperature to stabilize for at least 30 minutes on the automate before any use.

**Foot-note: Do not interchange the reagents from different lots.**

## 2. Determination of Protein C.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Plasma Calibrator (ref A222101)	1 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T°	30 minutes <i>on</i> STA <i>board</i> before any use (**)
Quality controls Biophen Normal Control (ref A223201) Biophen Abnormal Control (ref A223301)	1 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T°	30 minutes <i>on</i> STA <i>board</i> before any use (**)

### Reconstitution:

(\*) Following reconstitution of calibrators or controls with distilled water, let them to stabilize for 30 minutes at room temperature.

**Foot-note: A calibration curve must be carried out for each new batch of reagents.**

### Conservation of reagents:

(\*\*) Take care of strictly respecting the 30 minutes temperature stabilization time for *calibrators* and *controls* at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C. Homogeneize before each use.

**Foot-note: Do not freeze calibrators or controls.**

## 3. Results.

- The calibration curve is of the Lin (absorbance) - Lin (concentration) type.
- The values obtained for the patients and controls are directly calculated from the calibration curve.
- The results are expressed in % activity.  
The 100% activity is that of reference normal pooled citrated plasma.

#### 4. Programming the BCS analyzer.

##### Creation of reagents and diluent files

Return to definition of test. Go to “reagents without data of lot”. Create both reagents files, and also the diluent file.

<b>Name:</b> PROTAC			
Abbreviation: PROTAC			
clean when the reagent is modified			high
<input checked="" type="checkbox"/> New calibration with any new lot of reagents			
<b>Authorized lines</b>			
from 1		to 14	
<b>Stability</b>			
< 15C unlimited		> = 15C	
<b>Mix by jet</b>			
Intensity		1	
Frequency		No agitation	
<b>Reference numbers</b>			
		Number	
		add	remove
		replace	

<b>Name:</b> SAPC21			
Abbreviation: SAPC21			
clean when the reagent is modified			special
<input checked="" type="checkbox"/> New calibration with any new lot of reagents			
<b>Authorized lines</b>			
from 1		to 14	
<b>Stability</b>			
< 15C unlimited		> = 15C	
<b>Mix by jet</b>			
Intensity		1	
Frequency		No agitation	
<b>Reference numbers</b>			
		Number	
		add	remove
		replace	

<b>Name:</b> Physiological saline			
Abbreviation: PHYSAL			
clean when the reagent is modified			Normal/low
New calibration with any new lot of reagents			
<b>Authorized lines</b>			
from 1		to 14	
<b>Stability (*)</b>			
< 15C		> = 15C	
<b>Mix by jet</b>			
Intensity		1	
Frequency		No agitation	
<b>Reference numbers</b>			
		Number	
		add	remove
		replace	

Create in the same way calibration plasmas and control plasmas files.

<b>Name:</b> Biophen plasma Calibrator			
Abbreviation: PCCAL			
clean when the reagent is modified			Normal/low
<input checked="" type="checkbox"/> New calibration with any new lot of reagents			
<b>Authorized lines</b>			
from 1		to 14	
<b>Stability (*)</b>			
< 15C		> = 15C	
<b>Mix by jet</b>			
Intensity		1	
Frequency		No agitation	
<b>Reference numbers</b>			
		Number	
		add	remove
		replace	

<b>Name:</b> Biophen plasma Control			
Abbreviation: PC Cont			
clean when the reagent is modified			Normal/low
<input checked="" type="checkbox"/> New calibration with any new lot of reagents			
<b>Authorized lines</b>			
from 1		to 14	
<b>Stability (*)</b>			
< 15C		> = 15C	
<b>Mix by jet</b>			
Intensity		1	
Frequency		No agitation	
<b>Reference numbers</b>			
		Number	
		add	remove
		replace	

(\*) To be filled by the user

**Creation of the procedure**

Once the various reagent files are created, the procedure file must be created. In order to do it, go to definition of tests, procedure and type the following parameters.

**General Information**

Former Change			
Test Procedure Number			
Test Procedure Name	BIOPHEN PC 2.5		
Measurement time stopped after	120 sec	Or	mA (0 = anything time of measure)
	0		
Speed	normal		
Stirring	normal		
Wavelength	405 nm		
Primary absorbance range	131-160 mA		
Duplicate n# for samples and controls (*)		CV Authorized (*)	0.0 %
Duplicate n# for calibration (*)		CV Authorized (*)	0.0 %
Duplicate n# of summary of raw data	Arithmetic mean		
Dilution Factor	1 : 4 (1 = no dilution)		
Dilution Buffer	Physiological saline		

**(\*) To be filled by the user**

Evaluate and Control Method			
Delta A per minute			
from	20	Sec	120
			Sec
Correlation Factor		0.9	Absolute
Curve Type	Decrease		

## Pipetting Sequence

Test Procedure Number

Test Procedure Name

BIOPHEN PC 2.5

### Sampling cycles

Number of Cycle	Arm of transmission	Final Washing	Action of Rotors/transfer	Printout limits [sec]	
				Min	Max
1	Arm	normal	No mix. in priority	0	0
2	Reagent Arm	normal	Incubation	280	310
3	Reagent Arm	normal	No mix/measure	0	0

Num de transm.	Mean/pipetting goal	Reagent	Speed	Volume (µl)
1	Take air		Slow	20
1	Take the sample		Slow	50
1	In the external cup		Slow	50
2	Take air		Slow	20
2	Take the Reagent	R1 PROTAC	Slow	100
2	In the central cup		Slow	100
3	Take air		Slow	20
3	Take the Reagent	R2 SAPC21	Slow	100
3	In the central cup		Slow	100
4	Take the Reagent	SCS Cleaner	Slow	0

### Test Definition

<b>General Information</b>				
Former Change		[ ]		
Name		BIOPHEN PC 2.5		Test N°
Result Unit		PC ID t		
Behring	% dN	Conversion factor		Abbreviation
User	% dN	1%dN	1.0	1%dN
PC	% dN	1 %dN	1.0	1%dN
Place in the list				PC2.5
Digit number for results				1000

### Test detailed procedures

Procedure Name test	Evaluation Method
BIOPHEN PC 2.5	Delta A fix

### Calculate Formula of raw data

Formula Type	Minimum	
User own Formula	[ ]	
Unit	mE/min	

### Detailed evaluation test

Usual Values		Determination of result		<b>Standard Concentrations</b>	
		70	% dN		at
		140	%dN		% dN
Test used for Calibration		[ ]			100
Evaluation		Regression lin/lin			50
Reference Curve		Calibration curve			25
					12.5
					6.25
Minimum raw value		0.0			
Upper Extrapolation		1.5			
Lower Extrapolation		0.01			
Calibrator		Standard plasma calibrator			
Authorized SD respectively to the reference curve		0		%	
Maximal number repeats		0			

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