Manufactured By: HYPHEN BioMed

# BIOPHEN Protein C 2.5 ON STA

## Adaptation of BIOPHEN PROTEIN C 2.5 ON STA

### 1. Reconstitution of BIOPHEN PROTEIN C 2.5 (Ref A221202) reagents.

Chromogenic determination of Protein C.

	NAME	Reconstitution	Stability	Stabilization in T°
R1	Protac	2.5 ml of distilled water (*)	3 months at 2-8°C 3 days at room T° Do not freeze	30 mn before any use (**)
R2	SAPC-21 Substrate	2.5 ml of distilled water (*)	3 months at 2-8°C 3 days at room T° Do not freeze	30 mn before any use (**)

#### Reconstitution:

(\*) Following reconstitution with distilled water, let the R1 and R2 reagents to stabilize for 30 minutes at room temperature.

#### Conservation of reagents:

Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

#### Stabilization of reagents:

(\*\*) It is necessary to let the substrate (R2) and the Protac (R1) temperature to stabilize for at least 30 minutes on the automate before any use.

**Foot-note: Do not interchange the reagents from different lots.**

## 2. Determination of Protein C.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Plasma Calibrator (ref A222101)	1 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T°	30 minutes <i>on</i> STA <i>board</i> before any use (**)
Quality controls Biophen Normal Control (ref A223201) Biophen Abnormal Control (ref A223301)	1 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T°	30 minutes <i>on</i> STA <i>board</i> before any use (**)

### Reconstitution:

(\*) Following reconstitution of calibrators or controls with distilled water, let them to stabilize for 30 minutes at room temperature.

**Foot-note: A calibration curve must be carried out for each new batch of reagents.**

### Conservation of reagents:

(\*\*) Take care of strictly respecting the 30 minutes temperature stabilization time for *calibrators* and *controls* at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C. Homog eneize before each use.

**Foot-note: Do not freeze calibrators or controls.**

## 3. Results.

- The calibration curve is of the Lin (absorbance) - Lin (concentration) type.
- The values obtained for the patients and controls are directly calculated from the calibration curve.
- The results are expressed in % activity.  
The 100% activity is that of reference normal pooled citrated plasma.

#### 4. Programming STA Analyzer

Click on the icon **set up software** for the **manager program** and create the program according to:

Identification									
Abbreviation		Name			Method			Date	
PC 2.5		BIOPHEN PC2.5			Chro-kinetics OD/min				
Sample				Diluent					
Volume	Incu	Dil.	Id.	Name			vial.	stab	Min. Vol.
30µl	0 sec	1/1	PH	PHYSAL			5 ml	8 h	0.2 ml
reagents									
Id.	Name	Incu sec	Vol. µl	Vial ml	Stab h	Min.Vol ml	Washing		
							Before	after	
Ra	PROTAC	300	125	2.5	168	0.5	non	Spécial	
Rb									
Rc									
Rd	SAPC21		125	2.5	168	0.5	non	Intensive	
Analysis		Result				Validation (%)			
Linearity min	0.80	Primary Units		%		Standard	min	max	
First point	15 sec	Corrector		1.00			0	120	
Second point	80 sec	Single/Duplicate		single					
Rd Heating	Yes	Precision		5%					
Stirring	Yes	Redil. Condition							
		<							
		>							

Enter data concerning washing.

Washing				
	Ra	Rb	Rc	Rd
Before	No			No
After	Special			Intensive
Wash	12227			
Name	STA-DESORB U			
Stab. H	120			
Vial ml	15			
Vol in. ml	0.9			

Show the second page and write calibration data.

**CONFIGURATION – Page 2 / 3**

Calibration										
Mode graph – linear regression						Scale				
Calibrators				Vial ml	Stab. H	Min Vol ml	Raw data Concentration Single/Duplicate	Lin Lin Single		
Id.	Key	Name	Dil.							
BPC		BPCAL	1/1	1	8	0.2				
BPC		BPCAL	1/2	1	8	0.2				
BPC		BPCAL	1/4	1	8	0.2				
BPC		BPCAL	1/8	1	8	0.2				
ZERO		ZEROCAL	1/1	1	8	0.2				
Control										
	Id.	Key	Name	Vial. ml	Stab. H	Min Vol ml				
<b>Level 1</b>										
<b>Level 2</b>										

Show the last page and enter data concerning the quality control.

The STA can manage 3 levels of quality control. For example: the following configuration is possible.

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