Adaptation of BIOPHEN Protein C 5 ON BCS

1. Reconstitution of BIOPHEN Protein C 5 (Ref A221205) reagents.

Chromogenic determination of Protein C.

	NAME	Reconstitution	Stability	Stabilization in T°
R1	Protac	5 ml of distilled water (*)	3 months at 2-8°C 3 days at room T° Do not freeze	30 mn before any use (**)
R2	SAPC-21 Substrate	5 ml of distilled water (*)	3 months at 2-8°C 3 days at room T° Do not freeze	30 mn before any use (**)

Reconstitution:

(*) Following reconstitution with distilled water, let the R1 and R2 reagents to stabilize for 30 minutes at room temperature.

Conservation of reagents:

Take care of putting up the specific caps back on the bottles before storing them at 2°-8°C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Stabilization of reagents:

(**) It is necessary to let the substrate (R2) and the Protac (R1) temperature to stabilize for at least 30 minutes on the automate before any use.

Foot-note: Do not interchange the reagents from different lots.

2. Determination of Protein C.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Plasma Calibrator (ref A222101)	1 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T°	30 minutes <i>on</i> STA <i>board</i> before any use(**)
Quality controls Biophen Normal Control (ref A223201) Biophen Abnormal Control (ref A223301)	1 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T°	30 minutes <i>on</i> STA <i>board</i> before any use (**)

Reconstitution:

(*) Following reconstitution of calibrators or controls with distilled water, let them to stabilize for 30 minutes at room temperature.

Foot-note: A calibration curve must be carried out for each new batch of reagents.

Conservation of reagents:

(**) Take care of strictly respecting the 30 minutes temperature stabilization time for *calibrators* and *controls* at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C. Homogenize before each use.

Foot-note: Do not freeze calibrators or controls.

3. Results.

- The calibration curve is of the Lin (absorbance) - Lin (concentration) type.
- The values obtained for the patients and controls are directly calculated from the calibration curve.
- The results are expressed in % activity.
The 100% activity is that of reference normal pooled citrated plasma.

4. Programming the BCS analyzer.

Creation of reagents and diluent files

Return to definition of test. Go to reagents without data of lot. Create both reagent files, and also the diluent file.

Name:	PROTAC		
Abbreviation:	PROTAC		
clean when the reagent is modified	high		
<input checked="" type="checkbox"/> New calibration with any new lot of reagents			
Authorized lines			
from	1	to	14
Stability			
< 15C	unlimited	> = 15C	
Mix by jet			
Intensity	1		
Frequency	No agitation		
Reference numbers			
	Number		
	add	remove	replace

Name:	SAPC21		
Abbreviation:	SAPC21		
clean when the reagent is modified	special		
<input checked="" type="checkbox"/> New calibration with any new lot of reagents			
Authorized lines			
from	1	to	14
Stability			
< 15C	unlimited	> = 15C	
Mix by jet			
Intensity	1		
Frequency	No agitation		
Reference numbers			
	Number		
	add	remove	replace

Name:	Physiological saline		
Abbreviation:	PHYSAL		
clean when the reagent is modified	Normal/low		
New calibration with any new lot of reagents			
Authorized lines			
from	1	to	14
Stability			
< 15C		> = 15C	
Mix by jet			
Intensity	1		
Frequency	No agitation		
Reference numbers			
	Number		
	add	remove	replace

Create in the same way calibration plasmas and control plasmas files.

Name:	Biophen plasma Calibrator		
Abbreviation:	PC CAL		
clean when the reagent is modified	Normal/low		
<input checked="" type="checkbox"/> New calibration with any new lot of reagents			
Authorized lines			
from	1	to	14
Stability			
< 15C		> = 15C	
Mix by jet			
Intensity	1		
Frequency	No agitation		
Reference numbers			
	Number		
	add	remove	replace

Name:	Biophen plasma Control		
Abbreviation:	PC Cont		
clean when the reagent is modified	Normal/low		
<input checked="" type="checkbox"/> New calibration with any new lot of reagents			
Authorized lines			
from	1	to	14
Stability			
< 15C		> = 15C	
Mix by jet			
Intensity	1		
Frequency	No agitation		
Reference numbers			
	Number		
	add	remove	replace

(*) To be filled by the user

Creation of the procedure

Once the various reagents files are created, the procedure file must be created. In order to do it, go to **definition of tests**, procedure and type the following parameters.

General Information

Former Change			
Test Procedure Number			
Test Procedure Name	BIOPHEN PC 5		
Measurement time hopped after	120 sec	Or	mA (0 = anything time of measure)
	0		
Speed	normal		
Mix	normal		
Wavelength	405 nm		
Primary absorbance range	131-160 mA		
Duplicate n# for samples and controls (*)		CV Authorized (*)	0.0 %
Duplicate n# for calibration (*)		CV Authorized (*)	0.0 %
Duplicate n# of summary of raw data	Arithmetic mean		
Dilution Factor	1 : 4	(1 = no dilution)	
Dilution Buffer	Physiological saline		

(*) To be filled by the user

Evaluate and Control Method			
Delta A per minute			
from	20	Sec	120
			Sec
Correlation Factor	0.9		Absolute
Curve Type	Decrease		

Pipetting Sequence

Test Procedure Number

Test Procedure Name

BIOPHEN PC 5

Sampling cycles

Number of Cycle	Arm of transmission	Final Washing	Action of Rotors/transfer	Time Range [sec]	
				Min	Max
1	Arm	normal	No mix. in priority	0	0
2	Reagent Arm	normal	Incubation	280	310
3	Reagent Arm	normal	No mix/measure	0	0

Num de transm.	Mean/pipetting goal	Reagent	Speed	Volume (µl)
1	Take air		Slow	20
1	Take the sample		Slow	50
1	In the external cup		Slow	50
2	Take air		Slow	20
2	Take the Reagent	R1 PROTAC	Slow	100
2	In the central cup		Slow	100
3	Take air		Slow	20
3	Take the Reagent	R2 SAPC21	Slow	100
3	In the central cup		Slow	100
4	Take the Reagent	SCS Cleaner	Slow	0

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