Manufactured By: HYPHEN BioMed

BIOPHEN Protein C 5 ON STA

Adaptation of BIOPHEN PROTEIN C 5 ON STA

1. Reconstitution of BIOPHEN PROTEIN C 5 (Ref A221205) reagent.

Chromogenic determination of Protein C.

| | NAME | Reconstitution | Stability | Stabilization in T° |
|----|----------------------|--------------------------------|---------------------------------------------------------|------------------------------|
| R1 | Protac | 5 ml of distilled water (*) | 3 months at 2-8°C 3 days at room T° Do not freeze | 30 mn before any use (**) |
| R2 | SAPC-21 Substrate | 5 ml of distilled water (*) | 3 months at 2-8°C 3 days at room T° Do not freeze | 30 mn before any use (**) |

Reconstitution:

(*) Following reconstitution with distilled water, let the R1 and R2 reagents to stabilize for 30 minutes at room temperature.

Conservation of reagents:

Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Stabilization of reagents:

(**) It is necessary to let the substrate (R2) and the Protac (R1) temperature to stabilize for at least 30 minutes on the automate before any use.

Foot-note: Do not interchange the reagents from different lots.

2. Determination of Protein C.

| NAME | Reconstitution | Stability | Stabilization in T° |
|----------------------------------------------------------------------------------------------------------|-----------------------------|-----------------------------------------|----------------------------------------------------|
| Calibration Biophen Plasma Calibrator (ref A222101) | 1 ml of distilled water (*) | 24 hours at 2-8°C 8 hours at room T° | 30 minutes <i>on STA board</i> before any use(**) |
| Quality controls Biophen Normal Control (ref A223201) Biophen Abnormal Control (ref A223301) | 1 ml of distilled water (*) | 24 hours at 2-8°C 8 hours at room T° | 30 minutes <i>on STA board</i> before any use (**) |

Reconstitution:

(*) Following reconstitution of calibrators or controls with distilled water, let them to stabilize for 30 minutes at room temperature.

Foot-note: A calibration curve must be carried out for each new batch of reagents.

Conservation of reagents:

(**) Take care of strictly respecting the 30 minutes temperature stabilization time for *calibrators* and *controls* at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C. Homog enize before each use.

Foot-note: Do not freeze calibrators or controls.

3. Results.

- The calibration curve is of the Lin (absorbance) - Lin (concentration) type.
- The values obtained for the patients and controls are directly calculated from the calibration curve.
- The results are expressed in % activity.
The 100% activity is that of reference normal pooled citrated plasma.

4. Programming STA Analyzer

Click on the icon **set up software** for the **manager program** and create the program according to:

| Identification | | | | | | | | | |
|----------------|--------|------------------|---------|---------|----------------------|----------------|---------|-----------|---------|
| Abbreviation | | Name | | | Method | | | Date | |
| PC 5 | | BIOPHEN PC5 | | | Chro-kinetics OD/min | | | | |
| Sample | | | | Diluent | | | | | |
| Volume | Incu | Dil. | Id. | Name | | | Vial. | stab | Min Vol |
| 30µl | 0 sec | 1/1 | PH | PHYSAL | | | 5 ml | 8 h | 0.2 ml |
| reagents | | | | | | | | | |
| Id. | Name | Incu sec | Vol. µl | Vial ml | Stab h | Min Vol ml | Washing | | |
| | | | | | | | before | after | |
| Ra | PROTAC | 300 | 125 | 5 | 168 | 0.5 | non | Special | |
| Rb | | | | | | | | | |
| Rc | | | | | | | | | |
| Rd | SAPC21 | | 125 | 5 | 168 | 0.5 | non | Intensive | |
| Analysis | | Result | | | | Validation (%) | | | |
| Linearity min | 0.80 | Primary Units | | % | | Standard | min | max | |
| First point | 15 sec | Corrector | | 1.00 | | | 0 | 120 | |
| Second point | 80 sec | Single/Duplicate | | single | | | | | |
| Rd Heating | Yes | Precision | | 5% | | | | | |
| Stirring | Yes | Redil. Condition | | | | | | | |
| | | < | | | | | | | |
| | | > | | | | | | | |

Enter data concerning washing.

| Washing | | | | |
|------------|--------------|----|----|-----------|
| | Ra | Rb | Rc | Rd |
| Before | No | | | No |
| After | Special | | | Intensive |
| Wash | 12227 | | | |
| Name | STA-DESORB U | | | |
| Stab. H | 120 | | | |
| Vial ml | 15 | | | |
| Vol in. ml | 0.9 | | | |

Show the second page and write calibration data.

CONFIGURATION – Page 2 / 3

| Calibration | | | | | | | | |
|--------------------------------|-----|---------|------|----------|---------|-------------|-----------------------------------------------|----------------------|
| Mode graph – linear regression | | | | | | Scale | | |
| Calibrators | | | | Vial ml | Stab. H | Min Vol ml | Raw data Concentration Single/Duplicate | Lin Lin Single |
| Id. | Key | Name | Dil. | | | | | |
| BPC | | BPCAL | 1/1 | 1 | 8 | 0.2 | | |
| BPC | | BPCAL | 1/2 | 1 | 8 | 0.2 | | |
| BPC | | BPCAL | 1/4 | 1 | 8 | 0.2 | | |
| BPC | | BPCAL | 1/8 | 1 | 8 | 0.2 | | |
| ZERO | | ZEROCAL | 1/1 | 1 | 8 | 0.2 | | |
| Control | | | | | | | | |
| | Id. | Key | Name | Vial. ml | Stab. H | Min Vol. ml | | |
| Level 1 | | | | | | | | |
| Level 2 | | | | | | | | |

Show the last page and enter data concerning the quality control.

The STA can manage 3 levels of quality control. For example: the following configuration is possible.

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