Manufactured By: HYPHEN BioMed

BIOPHEN PROTEIN C 5 ON STA-R

Adaptation of BIOPHEN Protein C 5 ON STA-R

1. Reconstitution of BIOPHEN Protein C 5 (Ref A221205) reagent.

	NAME	Reconstitution	Stability	Stabilization in T°
R1	Protac	5 ml of distilled water (*)	3 months at 2-8°C 3 days at room T° Do not freeze	30 mn before any use (**)
R2	SAPC-21 Substrate	5 ml of distilled water (*)	3 months at 2-8°C 3 days at room T° Do not freeze	30 mn before any use (**)

Reconstitution:

(*) Following reconstitution with distilled water, let the R1 and R2 reagent to stabilize for 30 minutes at room temperature.

Conservation of reagents:

Take care of putting up the specific caps back on the bottles before storing them at 2°-8°C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Stabilization of reagents:

(**) It is necessary to let the substrate (R2) and the Protac (R1) temperature to stabilize for at least 30 minutes on the automate before any use.

Foot-note: Do not interchange the reagents from different lots.

2. Determination of Protein C :

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Plasma Calibrator (ref A222101)	1 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T°	30 minutes <i>on</i> STA -Rboard before any use(**)
Quality controls Biophen Normal Control (ref A223201) Biophen Abnormal Control (ref A223301)	1 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T°	30 minutes <i>on</i> STA-R board before any use (**)

Reconstitution:

(*) Following reconstitution of calibrators or controls with distilled water, let them to stabilize for 30 minutes at room temperature.

Foot-note: A calibration curve must be carried out for each new batch of reagents.

Conservation of reagents:

(**) Take care of strictly respecting the 30 minutes temperature stabilization time for *calibrators* and *controls* at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C. Homog eneize before each use.

Foot-note: Do not freeze calibrators or controls.

3. Results.

- The calibration curve is of the Lin (absorbance) - Lin (concentration) type.
- The values obtained for the patients and controls are directly calculated from the calibration curve.
- The results are expressed in % activity.

The 100% activity is that of reference normal pooled citrated plasma.

4. Programming STA-R Analyzer

Click on the icon **set up software** for the **manager program** and create the program according to:

TESTS						Configuration		
Method								
IDENTIFICATION								
Abbreviation		PC 5		Last up date				
Name		BIOPHEN PC 5		Method		Chro-Kinetics OD /min		
SAMPLE			DILUENT					
Vol. μ l	Incu.	Dilution	Id.	Name			Stab. h	Continued
30 μ l	0 sec	1/1	PH	PHYSAL			24	
REAGENTS								
	Id.	Name		Stab. h	Vol. μ l	Incu. sec	Prec	Continued
Ra	PROTAC	PROTAC		168	125	300	<input type="checkbox"/>	
Rb							<input type="checkbox"/>	
Rc							<input type="checkbox"/>	
Rd	SAPC 21	SAPC 21		168	125			

Enter data concerning washing

	Id.	Name	Stab. h	Vol. μ l	Incu. sec	Prec.	Vial. ml	Min Vol. ml
Diluent	PH	PHYSAL	24				5.0	1
Ra	PROTAC	PROTAC	168	125	300	<input type="checkbox"/>	5.0	0.5
Rb						<input type="checkbox"/>		
Rc						<input type="checkbox"/>		
Rd	SAPC 21	SAPC 21	168	125			5.0	0.5
Washing								
	Ra		Rb		Rc		Rd	
Before	No						No	
After	Special						Intensive	
Wash	12227							
Name	STA-DESORB U							
Stab. H	120							
Vial ml	15							
Vol in. ml	0.9							

Click on the icon **Result software** for the **manager program** and create the program according to:

TESTS		Configuration	
		Result	
METHOD CHRO-KINETICS OD/mn		RESULT	
Linearity Min.	0.80	Primary Units.	%
First Point	15	Corrector	1.00
Second Point	80	Single/Duplicate	<input checked="" type="radio"/> Single <input type="radio"/> Duplicate
		Precision (%)	5.00
<input checked="" type="checkbox"/> Rd Heating <input checked="" type="checkbox"/> Stirring		Redilution	Condition
		1/	<
		1/	>
		Validation (%)	
		Min.	Max.
		0	120

Write calibration data.

TESTS		Configuration	
		Calibration	
MODE :	Linear	Determination	
		<input checked="" type="checkbox"/> Duplicate <input type="checkbox"/> Single	
ASSAYS		Concentrations	
SCALE :	<input checked="" type="checkbox"/> Linearity <input type="checkbox"/> Log	<input checked="" type="checkbox"/> Linear <input type="checkbox"/> Inverse <input type="checkbox"/> Log	
	Id.	Key	Name
	Stab. H	Dilution	
Calibrator 1	CAL1		BPCAL
Calibrator 2	CAL2		BPCAL
Calibrator 3	CAL3		BPCAL
Calibrator 4	CAL4		BPCAL
Calibrator 5	CAL5		ZEROCAL
Ctrl, Niv. 1			
Ctrl, Niv. 2			

Click on the icon **Printout/Transmission software** for the **manager program** and create the program according to:

TESTS		Configuration				
		Printout/Trans.				
PARAMETERS						
	Units	Conver Factors	Print	Transmission Test number	Usual Values Min. Max.	
Main	%		<input checked="" type="checkbox"/>	0	70	120
Aux1	OD/min		<input checked="" type="checkbox"/>			
Aux2			<input type="checkbox"/>			
Aux3			<input type="checkbox"/>			
Printout Limits						
Min.	<input type="text" value="0"/>	Max.	<input type="text" value="120"/>			

Show the last page and enter data concerning the quality control.

TESTS		Configuration				
		Q.C				
LEVEL1						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control	BNC		BN CONT	8	1	0.2
LEVEL 2						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control	BAC		BA CONT	8	1	0.2
LEVEL3						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control						
Period : <input type="text"/> hours <input type="text"/> <input type="checkbox"/> Tests <input type="checkbox"/> Vial						

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