Biophen FVIII (#A221402/A221406)

on ACL7000 (research software)

Manufactured By: HYPHEN BioMed

**NOTE :** reconstitution, preparation, calibration and controls, expression of results: usual cautions and realization according to the specific technical insert D.750.02/BI/1402 or /1406.

**Reagents are reconstituted as per the insert, well equilibrated, and preheated at 37°C**

PHASE	CONDITION	OPERATION	PARAMETERS
1 <sup>st</sup> loading step	YES	Reagent loading (pos 18) <b>Reagent 1</b>	<b>50µl</b>
		Washing after reagent loading (pos 18)	2
		Diluent (pos DIL)	0µL
		<b>Specimen (working dilution : 1:40 in R4 for High Range, OR 1:10 in R4 for Low Range protocol).***</b>	<b>50µl</b>
		<b>Reagent (nacelle 2) Reagent R2</b>	<b>50µL</b>
		Washing between specimen	2
		Optical ref. (cycle abs.)	
		Diluent (pos DIL)	80µL
		Reagent (nacelle 2)	80µL
		Washing (at the end of loading)	1
		Incubation time (T1)	<b>300 s</b>
2 <sup>nd</sup> loading step	YES	Diluent (pos DIL)	0µL
		Specimen	0µL
		Washing between specimen	0
		<b>Reagent (nacelle 3) Reagent R3 substrate</b>	<b>50µL</b>
		Washing after reagent loading	1
		Mix time (T2)	<b>300 s*</b>
		Activation time (T3)	0 s
3 <sup>rd</sup> loading step	NO	Reagent (pos 17)	0µL
		Washing after reagent loading	0
		Mix time (T4)	0 s
		Activation time (T5)	0 s
Washing	NO	Reagent (pos 1)	0µL
		Reagent (pos 16)	0µL
		Cycles number	1
		Washing at the end of loading	1
Levels detection	NO		
Acquisition parameters	YES	Ramp	yes
		Inter-ramp interval (T6)	1s
		Delay before acquisition (T7)	<b>5s</b>
		Acquisition time (T8)	<b>300s*</b>
		Canal used	405nm
		Rotation speed	1200

\*could be reduced if required, but lower OD values expected.

Biophen FVIII:C on ACL7000 (research software)

1/2

D.750.07/BI/1402/1406/ACL7000

HYPHEN BioMed Version 1 04/02/2010

LLL

6560 Gove Court - Mason, OH 45040

Phone: 513.770.1991

Toll Free: 866.783.3797

Fax: 513.573.9241

Email: info@aniara.com

www.aniara.com

**Note \*\*\* :**

- **High range protocol :**

Calibrator at C% of FVIII:c or normal plasma pool considered at 100% FVIII :C (refer insert).  
The working dilution is 1:40 in R4 buffer, and corresponds by definition to C% or 100% FVIII :C.

Example of calibration curve obtained for the high range:

% FVIII :c	A405
156%	1.08
78% (or « C% », tested at the 1 :40 dilution)	0.62
39%	0.28
9.8%	0.036
$r^2$ ( logA405-log conc)	>0.98

- **Low range protocol :**

Calibrator at C% of FVIII:c or normal plasma pool considered at 100% FVIII :C, diluted 1:4 (or adequate dilution) in FVII:C deficient plasma to get a calibrator at 25% of FVIII:C (refer insert).  
The working dilution is 1:10 in R4 buffer, and corresponds by definition to 25% FVIII :C.

Note: The stability at 2-8°C and at room temperature (18-25°C) claimed are validated on reconstituted vials kept closed, protected from any evaporation or contamination.  
Stabilities of reagents on each automate must be controlled, adjusted and validated under each exact laboratory working conditions.

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