BIOPHEN FACTOR VIII ON BCS

Provisory Not VALIDATED

Adaptation of BIOPHEN Factor VIII on BCS

1. Reconstitution of the BIOPHEN Factor VIII reagents.

Chromogenic determination of Factor VIII:C with Biophen Factor VIII:C.

BIOPHEN Factor VIII (Ref A221402) reagents

	NAME	Reconstitution	Stability*	Stabilization in T°
R1	F.X	2.5 ml of distilled water *	Refer insert	** 30 mn before any use
R2	Activator	2.5 ml of distilled water *		** 30 mn before any use
R3	SXa-11	2.5 ml of distilled water		** 30 mn before any use
R4	TRIS-BSA	Ready to use		** 30 mn before any use

BIOPHEN Factor VIII (6) (Ref A221406) reagents

	NAME	Reconstitution	Stability*	Stabilization in T°
R1	F.X	6.0 ml of distilled water *	Refer insert	** 30 mn before any use
R2	Activator	6.0 ml of distilled water *		** 30 mn before any use
R3	SXa-11	6.0 ml of distilled water		** 30 mn before any use
R4	TRIS-BSA	Ready to use		** 30 mn before any use

*The stability data claimed on the insert were obtained on reconstituted vials, kept closed, protected from and provided any contamination or evaporation is avoided. Stability must be controlled, and can be adjusted and validated if required, according to the exact use conditions for each laboratory.

Reconstitution:

(*) Following reconstitution with distilled water, leave the R1, R2 and R3 reagents to stabilize for 30 minutes at room temperature.

Provisory Not VALIDATED

Conservation of reagents:

Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Stabilization of reagents:

(**) It is necessary to leave the substrate (R3), (R2) (R1) and (R4) temperature to stabilize for at least 30 minutes on the automate before any use

Foot-note: Do not interchange the reagents from different lots.

2.Preparation of the high range concentration calibrators

- Calibration is performed with normal pooled citrated plasma with the assigned value of 100% Factor VIII: C. The assay includes a standard plasma dilution of 1:20. By definition, this later dilution of the pool represents the 100% Factor VIII: C activity. The dynamic range is from 0 to 100% Factor VIII: C.

Or

-Calibration is performed with plasma calibrator commercially available, with a known factor VIII: C concentration. The 1:20 dilution corresponds to the indicated Factor VIII: C concentration.

Tested plasma:

In these conditions, the tested plasmas and controls are assayed at the 1:20 dilution.

Preparation of the low range concentration calibrators

- Calibration is performed with normal pooled citrated plasma with the assigned value of 100% Factor VIII: C. The normal pooled citrated plasma is diluted 1:4 in a Factor VIII:C deficient plasma in order to obtain a concentration of 25% Factor VIII: C (1 volume of normal pooled citrated + 3 volume of deficient factor VIII:C plasma).

This Normal citrated pooled plasma, diluted 1:4 with VIII: C deficient plasma is then diluted at 1:10. By definition this later dilution of the pool represents the 25% Factor VIII: C activity. The dynamic range is from 0 to 25 % Factor VIII: C activity.

Or

- Calibration is performed with a commercially available plasma calibrator, with a known factor VIII: C concentration. Following reconstitution, the calibrator must be appropriately diluted in Factor VIII deficient plasma in order to obtain a concentration of 25% Factor VIII: C.

This calibrator is then diluted 1:10. By definition this later dilution of the calibrator represents the 25% Factor VIII: C activity. The dynamic range is from 0 to 25 % Factor VIII: C activity.

Tested plasma:

In these conditions, the tested plasmas and controls are assayed at the 1:10 dilution.

A calibration curve must be carried out for each new batch of reagents.

3.Controls

The control is performed with commercially available control plasmas, titrated for Factor VIII: C.

Nota : Following reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature. It is necessary to let the reagent temperature to stabilize for at least 30 minutes onto the automate before any use.

Do not freeze quality control plasmas

Quality controls must be run regularly, and for each new batch of reagents, after an important maintenance of the instrument, or if measured values are not in compliance with the one expected for the method.

4. Results:

- **HIGH CALIBRATION RANGE (0 to 100 %):** The calibration curve (working range) is of the Log (absorbance) – Log (concentration) type.
- **LOW CALIBRATION RANGE (0 to 25 %):** The calibration curve (low concentration working range) is of the Lin (absorbance) - Lin (concentration) type.
- The values obtained for the patients and controls are directly calculated from the calibration curve.
- The results are expressed as % factor VIII:C activity (100 % being equivalent to 1 Unit/ml)

The calibration curve is validated when linearity (r^2), as well as measured control values, are in compliance.

A new calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).

Performances may present variations according to the instrument used. Validate the expected values in the laboratory working conditions. Performances, as well as values for each new lot of quality controls used, must then be confirmed (and adjusted if necessary) in the laboratory working conditions.

5. Programming the BCS analyzer:

Click on the icon **set up software** for the **manager program** and create the program according to:

High Factor VIII: C concentration working range:

This calibration allows measuring concentrations in the range 0 -200% of Factor VIII: C

Creation of reagents and diluent files

Return to definition of test. Go into reagents without data of lot. Create the 3 reagent files, as well as that of Diluent.

<p>Name: F.X</p> <p>Abbreviation: R1</p> <p>clean when the reagent is modified Normal/low</p> <p><input checked="" type="checkbox"/> New calibration with any new lot of reagents</p> <p>Authorized lines from 1 to 14</p> <p>Stability < 15C * > = 15C</p> <p>Mix by jet Intensity 1 Frequency No agitation</p> <p>Reference numbers Number add remove replace</p>	<p>Name: Activator</p> <p>Abbreviation: R2</p> <p>clean when the reagent is modified special</p> <p><input checked="" type="checkbox"/> New calibration with any new lot of reagents</p> <p>Authorized lines from 1 to 14</p> <p>Stability < 15C * > = 15C</p> <p>Mix by jet Intensity 1 Frequency No agitation</p> <p>Reference numbers Number add remove replace</p>
<p>Name: SXa-11</p> <p>Abbreviation: R3</p> <p>clean when the reagent is modified High</p> <p><input checked="" type="checkbox"/> New calibration with any new lot of reagents</p> <p>Authorized lines from 1 to 14</p> <p>Stability < 15C * > = 15C</p> <p>Mix by jet Intensity 1 Frequency No agitation</p> <p>Reference numbers Number add remove replace</p>	<p>Name: Buffer Tris-BSA</p> <p>Abbreviation: R4</p> <p>clean when the reagent is modified Normal/low</p> <p><input type="checkbox"/> New calibration with any new lot of reagents</p> <p>Authorized lines from 1 to 14</p> <p>Stability < 15C * > = 15C</p> <p>Mix by jet Intensity 1 Frequency No agitation</p> <p>Reference numbers Number add remove replace</p>

*user defined; **Provisory Not VALIDATED**

Creation of the procedure

Once the various reagents files are created, the procedure file must be created. In order to do it,

General Information			
Former Change			
Test Procedure Number			
Test Procedure Name	BIOPHEN F.VIII C		
Measurement time hopped after	250 sec	Or	mA (0 = anything time of measure)
	0		
Speed	normal		
Mix	normal		
Wavelength	405 nm		
Primary absorbance range	131-160 mA		
Duplicate n# for samples and controls	1	CV Authorized	5.0 %
Duplicate n# for calibration	2	CV Authorized	5.0 %
Duplicate n# of summary of raw data	Arithmetic mean		
Factor dilution	1 : 2	(1 = no dilution)	
Buffer dilution	Tris R4		

Evaluate and Control Method			
Delta A per minute			
from	10	Sec	240
			Sec
Correlation Factor			0.9
			Absolute

Provisory Not VALIDATED

Pipetting Sequence

Test Procedure Number

Test Procedure Name

BIOPHEN F.VIII C

Sampling cycles

Number of Cycle	Arm of transmission	Final Washing	Action of Rotors/transfer	Time Range [sec]	
				Min	Max
1	Arm	normal	No mix. in priority	0	0
2	Reagent Arm	special	No mix. in priority	0	0
3	Reagent Arm	High	Incubation	290	310
4	Reagent Arm	normal	No mix/measure	0	0

Num de transm.	Mean/pipetting goal	Reagent	Speed	Volume (µl)
1	Take up air		Slow	20
1	Take up reagent	R4	Slow	45
1	Take up sample		Slow	5
1	Dispense outer cuvette		Slow	50
2	Take up air		Slow	20
2	Take up Reagent	R1	Slow	50
2	Dispense center cuvette		Slow	50
3	Take up air		Slow	20
3	Take the Reagent	R2	Slow	50
3	Dispense center cuvette		Slow	50
4	Take up air		Slow	20
4	Take the Reagent	R3	Slow	50
4	Dispense inner cuvette		Slow	50

Provisory Not VALIDATED

Test Definition					
General Information					
Former Change					
Name		BIOPHEN F.VIII C		°Test N°	
Result Unit				PC ID t	
Behring	%	Conversion factor		Abbreviation	
User	%	1.0	%	Place in the list	
PC	%	1.0	%	Digit number for results	
				F8C	
				100,0	

Test detailed procedures

Procedure Name test	Method of evaluation
BIOPHEN F.VIII C	Delta A per minute

Calculate Formula of raw data

Formula Type	Minimum	
User own's Formula		
Raw data unit	mE/min	

Detailed test evaluation

Determination of result		Standard Concentrations
Reference range	0.0 % dN at 100 %	
Test used for Calibration		
Evaluation	Regression log/log	
Reference Curve	Calibrate curve	
Minimum raw value	0.0	
Upper Extrapolation	1.1 x maxi concentration	
Lower Extrapolation	1.0 x mini concentration	
Calibrator	FVIII calibrator	
Authorized SD respectively to the reference curve	0 %	
Maximal number repeats	0	

Provisory Not VALIDATED

Click on the icon **Edition/Transmission software** for the **manager program** and create the program according to:

Low concentration Factor VIII: C working range:

This calibration allows measuring concentration between 0 and 25 % of Factor VIII: C activity

General Information			
Former Change			
Test Procedure Number			
Test Procedure Name	BIOPHEN F.VIII C Low		
Measurement time hopped after	250 sec	Or	mA (0 = anything time of measure)
	0		
Speed	normal		
Mix	normal		
Wavelength	405 nm		
Primary absorbance range	131-160 mA		
Duplicate n# for samples and controls	1	CV Authorized	5.0 %
Duplicate n# for calibration	2	CV Authorized	5.0 %
Duplicate n# of summary of raw data	Arithmetic mean		
Factor dilution	1 : 1	(1 = no dilution)	
Buffer dilution	Tris R4		

Evaluate and Control Method			
Delta A per minute			
from	10	Sec	240
			Sec
Correlation Factor			0.9
			Absolute

Provisory Not VALIDATED

Pipetting Sequence					
Test Procedure Number					
Test Procedure Name			BIOPHEN F.VIII C Low		
Sampling cycles					
Number of Cycle	Arm of transmission	Final Washing	Action of Rotors/transfer	Time Range [sec]	
				Min	Max
1	Arm	normal	No mix. in priority	0	0
2	Reagent Arm	special	No mix. in priority	0	0
3	Reagent Arm	High	Incubation	290	310
4	Reagent Arm	normal	No mix/measure	0	0
Num de transm.	Mean/pipetting goal	Reagent	Speed	Volume (µl)	
1	Take up air		Slow	20	
1	Take up reagent	R4	Slow	45	
1	Take up sample		Slow	5	
1	Dispense outer cuvette		Slow	50	
2	Take up air		Slow	20	
2	Take up Reagent	R1	Slow	50	
2	Dispense center cuvette		Slow	50	
3	Take up air		Slow	20	
3	Take the Reagent	R2	Slow	50	
3	Dispense center cuvette		Slow	50	
4	Take up air		Slow	20	
4	Take the Reagent	R3	Slow	50	
4	Dispense inner cuvette		Slow	50	

Provisory Not VALIDATED

Test Definition					
General Information					
Former Change					
Name		BIOPHEN F.VIII C Low		°Test N°	
Result Unit				PC ID t	
Behring		%		Conversion factor	
User		%		1.0 %	
PC		%		1.0 %	
				Abbreviation	
				F8CL	
				Place in the list	
				Digit number for results	
				100,0	

Test detailed procedures

Procedure Name test	Method of evaluation
BIOPHEN F.VIII c	Delta A per minute

Calculate Formula of raw data

Formula Type	Minimum	
User own's Formula		
Raw data unit	mE/min	

Detailed test evaluation

		Determination of result		
Reference range	0.0	% dN	at	25 %
Test used for Calibration				Standard Concentrations
Evaluation	Regression lin/lin			25%
Reference Curve	Calibrate curve			12.5
				6.25
				3.12
				1.56
Minimum raw value	0.0			
Upper Extrapolation	1.0			x maxi concentration
Lower Extrapolation	0.5			x mini concentration
Calibrator	FVIII calibrator (low)			
Authorized SD respectively to the reference curve	0			%
Maximal number repeats	0			

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