Adaptation of BIOPHEN Factor VIII on STA and STA-C

1. Reconstitution of the BIOPHEN Factor VIII reagents.

Chromogenic determination of Factor VIII:C with Biophen Factor VIII:C.

BIOPHEN Factor VIII (Ref A221402) reagents

	NAME	Reconstitution	Stability*	Stabilization in T°
R1	F.X	2.5 ml of distilled water *	72 hours at 2-8°C * 24 hours at room T°	** 30 mn before any use
R2	Activator	2.5 ml of distilled water *	72 hours at 2-8°C * 24 hours at room T°	** 30 mn before any use
R3	SXa-11	2.5 ml of distilled water	3 months at 2-8°C * 7 days at room T°	** 30 mn before any use
R4	TRIS-BSA	Ready to use	Stable until the expiration date at 2-8°C	** 30 mn before any use

BIOPHEN Factor VIII (6) (Ref A221406) reagents

	NAME	Reconstitution	Stability*	Stabilization in T°
R1	F.X	6.0 ml of distilled water *	72 hours at 2-8°C * 24 hours at room T°	** 30 mn before any use
R2	Activator	6.0 ml of distilled water *	72 hours at 2-8°C * 24 hours at room T°	** 30 mn before any use
R3	SXa-11	6.0 ml of distilled water	3 months at 2-8°C * 7 days at room T°	** 30 mn before any use
R4	TRIS-BSA	Ready to use	Stable until the expiration date at 2-8°C	** 30 mn before any use

*The stability data claimed on the insert were obtained on reconstituted vials, kept closed, protected from and provided any contamination or evaporation is avoided. Stability must be controlled, and can be adjusted and validated if required, according to the exact use conditions for each laboratory.

Reconstitution:

(*) After reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature.

Stabilization of reagents:

(**) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate before any use.

Conservation of reagents:

Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Foot-note: Do not interchange the reagents from different lots.

NOTE: Use microcups when running the assay, to improve management of the reagents volumes by the automate.

2.Preparation of the high range concentration calibrators

- Calibration is performed with normal pooled citrated plasma with the assigned value of 100% Factor VIII: C. The assay includes a standard plasma dilution of 1:40. By definition, this later dilution of the pool represents the 100% Factor VIII: C activity. The dynamic range is from 0 to 200% Factor VIII: C. The 200% Factor VIII activity is the 1:20 dilution of the plasma pool.

Or

-Calibration is performed with plasma calibrator commercially available, with a known factor VIII: C concentration. The 1:40 dilution corresponds to the indicated Factor VIII:C concentration, and the 1:20 to twice this concentration.

Tested plasma:

In these conditions, the tested plasmas and controls are assayed at the 1:40 dilution.

Preparation of the low range concentration calibrators

- Calibration is performed with normal pooled citrated plasma with the assigned value of 100% Factor VIII: C. The normal pooled citrated plasma is diluted 1:4 in a Factor VIII:C deficient plasma in order to obtain a concentration of 25% Factor VIII: C (1 volume of normal pooled citrated + 3 volume of deficient factor VIII:C plasma).

This Normal citrated pooled plasma is then diluted at 1:10. By definition this later dilution of the pool represents the 25% Factor VIII:C activity. The dynamic range is from 0 to 25 % Factor VIII: C activity.

Or

- Calibration is performed with a commercially available plasma calibrator, with a known factor VIII: C concentration. Following reconstitution, the calibrator must be appropriately diluted in Factor VIII deficient plasma in order to obtain a Factor VIII: C concentration of 25%.

This calibrator is then diluted 1:10. By definition this later dilution of the calibrator represents the 25% Factor VIII: C activity. The dynamic range is from 0 to 25 % Factor VIII:C activity.

Tested plasma:

In these conditions, the tested plasmas and controls are assayed at the 1:10 dilution.

A calibration curve must be carried out for each new batch of reagents.

3.Controls

The control is performed with commercially available control plasmas, titrated for Factor VIII: C.

Nota : Following reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature. It is necessary to let the reagent temperature to stabilize for at least 30 minutes onto the automate before any use.

Do not freeze quality control plasmas

Quality controls must be run regularly, and for each new batch of reagents, after an important maintenance of the instrument, or if measured values are not in compliance with the one expected for the method.

4. Results:

- **HIGH CALIBRATION RANGE (0 to 200 %):** The calibration curve (working range) is of the Log (absorbance) – Log (concentration) type.
- **LOW CALIBRATION RANGE (0 to 25 %):** The calibration curve (low concentration working range) is of the Lin (absorbance) - Lin (concentration) type.
- The values obtained for the patients and controls are directly calculated from the calibration curve.
- The results are expressed as % factor VIII:C activity (100 % being equivalent to 1 Unit/ml)

The calibration curve is validated when linearity (r²), as well as measured control values, are in compliance.

A new calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).

Performances may present variations according to the instrument used. Validate the expected values in the laboratory working conditions. Performances, as well as values for each new lot of quality controls used, must then be confirmed (and adjusted if necessary) in the laboratory working conditions.

5. Programming the STA analyser:

Click on the icon **set up software** for the **manager program** and create the program according to:

1. High Factor VIII:C concentration working range :

This calibration allows measuring concentrations in the range 0 -200% of Factor VIII: C

CONFIGURATION – Page 1 / 3

Identification										
Abbreviation		Name			Method			Date		
B F.VIII		BIOPHEN F.VIII			Two points Chromo					
Sample				Diluent						
Volume	Incu	Dil.	Id.	Name			Vial.	stab	Min.Vol.	
75µl	0 sec	1/40	R4	R4 VIII			5 ml	24 h	0.2 ml	
reagents										
Id.		Name		Incu sec	Vol. µl	Vial ml	Stab h	Min.Vol. ml	Washing	
									before	after
Ra	R1VIII	R1 VIII		0	75	3*	24	0.2*	no	Normal
Rb	R2VIII	R2 VIII		300	75	3*	24	0.2*	no	Special
Rc										
Rd	R3VIII	R3 VIII			75	3*	24	0.2*	no	Intensive
Analysis			Result				Validation (%)			
Rd Heating	yes	Primary Units			%		Standard	min	max	
Stirring	yes	Corrector			1.00			5	200	
First point	20	Single/Duplicate			single					
Last point	140	Precision			5%					
		Redil. Condition								
<										
>										

*vial 3ml and min vol 0.2ml for A221402 ; Vial 6ml and min vol 0.2ml for A221406.

Enter data concerning washing.

Washing	Ra	Rb	Rc	Rd
Front	No	No		No
After	Normal	Special		Intensive
Wash		12227		
Name		STA-DESORB U		
Stab. H		120		
Vial ml		15		
Vol in. ml		0.9		

Show the second page and write calibration data.

CONFIGURATION – Page 2 / 3

Calibration								
Mode graph – linear regression					Scale			
Calibrators				Vial* ml	Stab.* H	Min.Vol.* ml	Raw Data Concentration Single/Duplicate	Log Log Single
Id.	Key	Name	Dil.					
VIIICAL		VIIICAL	1/20	1	8	0.2		
VIIICAL		VIIICAL	1/40	1	8	0.2		
VIIICAL		VIIICAL	1/80	1	8	0.2		
VIIICAL		VIIICAL	1/200	1	8	0.2		
VIIICAL		VIIICAL	1/400	1	8	0.2		
Control								
	Id.	Key	Name	Vial. ml	Stab. H	Min.Vol. ml		
Level 1								
Level 2								

(*) To be filled by the user . Homogenize before each use

Show the last page and enter data concerning the quality control.

The STA can manage 3 levels of quality control. For example: the following configuration is possible.

CONFIGURATION – Page 3 / 3

Printout / Transmission							
Units	Factor convers.	Print	Transmission Test number	Usual Values	%		
Main	%	Yes	0	STANDARD	Min.	Max.	
Aux.1	OD/min	Yes	0		5	200	
Aux.2							
Aux.3							
Printout limits		min	0				
		max	200				
Controls							
	Id.	Key	Name	Period H	Vial. ml	Stab. H	Min.Vol. ml
Level 1	(*)	(*)	(*)	(*)	(*)	(*)	(*)
Level 2	(*)	(*)	(*)	(*)	(*)	(*)	(*)
Level 3	(*)	(*)	(*)	(*)	(*)	(*)	(*)

(*) To be filled by the user . Homogenize before each use

2. Low Factor VIII:C concentration working range:

This calibration allows measuring concentration between 0 and 25 % of Factor VIII:C activity

CONFIGURATION – Page 1 / 3

Identification									
Abbreviation		Name			Method			Date	
B F.VIII		BIOPHEN F.VIII			Two points Chrom				
Sample				Diluent					
Volume	Incu	Dil.	Id.	Name			vial.	stab	Min.Vol
75µl	0 sec	1/10	R4	R4 VIII			5 ml	8 h	0.2 ml
reagents									
Id.	Name	Incu sec	Vol. µl	Vial ml	Stab h	Min.Vol ml	Washing		
							before	after	
Ra R1VIII	R1 VIII	0	75	3*	24	0.2*	no	Normal	
Rb R2VIII	R2 VIII	300	75	3*	24	0.2*	no	Special	
Rc									
Rd R3VIII	R3 VIII		75	3*	24	0.2*	no	Intensive	
Analysis		Result				Validation (%)			
Rd Heating Stirring First point Last point	yes	Primary Units		%		Standard	min	max	
	yes	Corrector		1.00			0	25	
	20 sec	Single/Duplicate		single					
	300 sec	Precision		5%					
		Redil. Condition							
		<							
		>							

*vial 3ml and min vol 0.2ml for A221402 ; Vial 6ml and min vol 0.2ml for A221406.

Enter data concerning washing.

Washing	Ra	Rb	Rc	Rd
Before	No	No		No
After	Normal	Special		Intensive
Wash		12227		
Name		STA-DESORB U		
Stab. H		120		
Vial ml		15		
Vol in. ml		0.9		

Show the second page and write calibration data.

CONFIGURATION – Page 2 / 3

Calibration								
Mode graph – linear regression						Scale		
Calibrators				Vial* ml	Stab.* H	Min.Vol* ml	Raw Data Concentration Single/Duplicate	Lin Lin Single
Id.	Key	Name	Dil.					
VIIICALB		VIIICALB	1/10	1	8	0.2		
VIIICALB		VIIICALB	1/20	1	8	0.2		
VIIICALB		VIIICALB	1/40	1	8	0.2		
VIIICALB		VIIICALB	1/80	1	8	0.2		
VIIICALB		VIIICALB	1/160	1	8	0.2		
Control								
	Id.	Key	Name	Vial. ml	Stab. H	Min Vol		
Level 1								
Level 2								

(*) To be filled by the user . Homogenize before each use

Show the last page and enter data concerning the quality control.

The STA can manage 3 levels of quality control. For exemple : the following configuration is possible.

CONFIGURATION – Page 3 / 3

Printout / Transmission							
Units	Factor Conver	Print	Transmission Test number	Usual Values	%		
Main	%	Yes	0	STANDARD	Min	Max.	
Aux.1	OD/min	Yes	0		0	25	
Aux.2							
Aux.3							
Range data		min	0				
		max	25				
Controls							
	Id.	Key	Name	Period H	Vial. ml	Stab. H	Min Vol ml
Level 1	(*)	(*)	(*)	(*)	(*)	(*)	(*)
Level 2	(*)	(*)	(*)	(*)	(*)	(*)	(*)
Level 3	(*)	(*)	(*)	(*)	(*)	(*)	(*)

(*) To be filled by the user . Homogenize before each use

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