Manufactured By: HYPHEN BioMed

BIOPHEN Plasminogen ON BCS

Adaptation of BIOPHEN Plasminogen ON BCS

1. Reconstitution of BIOPHEN Plasminogen (Ref A221502) reagents:

Chromogenic determination of plasminogen activity.

	NAME	Reconstitution	Stability	Stabilization in T°
R1	Streptokinase (activator)	2.5ml of distilled water (*)	1 month at 2-8°C 7 days at room T° Do not freeze	30 mn before any use (**)
R2	Spm41 Substrate	2.5 ml of distilled water (*)	1 month at 2-8°C 7 days at room T° Do not freeze	30 mn before any use (**)
Diluent	Physiological Saline /sodium chloride		24 hours <i>on BCS board</i>	

Reconstitution:

(*) Following reconstitution with distilled water, let the reagents to stabilize for 30 minutes at room temperature.

Storage of reagents:

Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Stabilization of reagents:

(**) It is necessary to let the reagents temperature to stabilize for at least 30 minutes on the automate before any use.

Foot-note: Do not interchange the reagents from different lots.

NOTE: If necessary; use microcups when running the assay, to improve management of the reagents volumes by the automate.

2. Determination of Plasminogen:

A normal human plasma pool, titrating by definition 100% plasminogen activity, can be used to establish the calibration curve.

Alternatively, the calibration and control is performed with commercially available control plasmas, titrated for plasminogen activity("C"%).

(example: Biophen Plasma calibrator #A222101, Biophen Normal Control Plasma #A223201, Biophen Abnormal control Plasma #A223301):

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Plasma Calibrator (ref A222101)	1 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T°	30 minutes on BCS board before any use (**)
Quality controls Biophen Normal Control (ref A223201) Biophen Abnormal Control (ref A223301)	1 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T°	30 minutes on BCS board before any use (**)

Note : samples are loaded undiluted on the instrument. The 1:20 working dilution is directly managed by the instrument, and corresponds to the indicated concentration ("C").

Note : if lyophilized calibrator and controls are used:

- Following reconstitution of calibrators or controls with distilled water, let them to stabilize for 30 minutes at room temperature.
- Take care of strictly respecting the 30 minutes temperature stabilization time for *calibrators* and *controls* at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C. Homogenize before each use.

Foot-note: A calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, and if QC values are out of the acceptance range for the method.

Foot-note: Do not freeze calibrators or controls.

3. Results:

- The calibration curve is of the Lin (absorbance) - Lin (concentration) type.
- The values obtained for the patients and controls are directly calculated from the calibration curve, and the results are expressed in % activity. .

4. Programming the BCS analyzer:

Creation of reagents and diluent files

Return to definition of test. Go in reagents without data of lot. Create both reagents files, and also the Diluent file.

Name: Streptokinase Plg BIOPHEN			
Abbreviation: Plg Act		clean when the reagent is modified special	
<input checked="" type="checkbox"/> New calibration with any new lot of reagents			
Authorized lines			
from	1	to	14
Stability			
< 15C	unlimited	> = 15C	
Mix by jet			
Intensity	1		
Frequency	No agitation		
Reference numbers			
	Number		
	add	remove	replace

Name: Substrate Plg BIOPHEN			
Abbreviation: Plg Sub		clean when the reagent is modified high	
<input checked="" type="checkbox"/> New calibration with any new lot of reagents			
Authorized lines			
from	1	to	14
Stability			
< 15C	unlimited	> = 15C	
Mix by jet			
Intensity	1		
Frequency	No agitation		
Reference numbers			
	Number		
	add	remove	replace

Name: Physiological saline / sodium chloride			
Abbreviation: Sodium chloride		clean when the reagent is modified Normal/low	
<input type="checkbox"/> New calibration with any new lot of reagents			
Authorized lines			
from	1	to	14
Stability (*)			
< 15C		> = 15C	
Mix by jet			
Intensity	1		
Frequency	No agitation		
Reference numbers			
	Number		
	add	remove	replace

Create in the same way calibration plasmas and control plasmas files.

Name: Calibrator			
Abbreviation: CAL		clean when the reagent is modified Normal/low	
<input checked="" type="checkbox"/> New calibration with any new lot of reagents			
Authorized lines			
from	1	to	14
Stability (*)			
< 15C		> = 15C	
Mix by jet			
Intensity	1		
Frequency	No agitation		
Reference numbers			
	Number		
	add	remove	replace

Name: Control			
Abbreviation: CT		clean when the reagent is modified Normal/low	
<input checked="" type="checkbox"/> New calibration with any new lot of reagents			
Authorized lines			
from	1	to	14
Stability (*)			
< 15C		> = 15C	
Mix by jet			
Intensity	1		
Frequency	No agitation		
Reference numbers			
	Number		
	add	remove	replace

(*) To be filled by the user

Creation of the procedure

Once the various reagents files are created, the procedure file must be created. In order to do it, go to definition of tests, procedure and type the following parameters.

General Information

Former Change			
Test Procedure Number			
Test Procedure Name	BIOPHEN Plasminogen		
Measurement time hopped after	90 sec	Or	mA
	0		
Speed	normal		
Stirring	normal		
Wavelength	405 nm		
Primary absorbance range	131-160 mA		
Duplicate n# for samples and controls (*)	2	CV Authorized (*)	10.0 %
Duplicate n# for calibration (*)	2	CV Authorized (*)	10.0 %
Duplicate n# of summary of raw data	Arithmetic mean		
dilution Factor	1 : 10	(1 = no dilution)	
Dilution Buffer	Sodium chloride		

(*) To be filled by the user

Evaluate and Control Method			
Lag time 10 sec			
Absorption min/max 5000mA			
Delta A per minute			
from	10	Sec to	80
			Sec
Correlation Factor	1.0		Absolute

Pipetting Sequence

Test Procedure Number

Test Procedure Name

BIOPHEN plasminogen

Sampling cycles

Number of Cycle	Arm of transmission	Final Washing	Action of Rotors/transfer	Time Range [sec]	
				Min	Max
1	Arm	normal	No mix. in priority	0	0
2	Reagent Arm	normal	Incubation	180	190
3	Reagent Arm	normal	No mix/measure	0	0

Sampling milieu

Num de transm.	Mean/pipetting goal	Reagent	Speed	Volume (µl)
1	Take air		Slow	20
1	Take the Reagent	Sodium chloride	Slow	50
1	Take the sample		Slow	50
1	In the external cup		Slow	100
2	Take air		Slow	20
2	Take the Reagent	Plg Act	Slow	100
2	In the central cup		Slow	100
3	Take air		Slow	20
3	Take the Reagent	Plg Substrate	Slow	100
3	In the internal cup		Slow	100

Test Definition

General Information				
Former Change		<input type="text"/>		
Name		<input type="text" value="BIOPHEN plasminogen"/>		°Test N°
Result Unit		PC ID t		
Behring	% dN	Conversion factor		Abbreviation
User	% dN	1%dN	1.0	1%dN
PC	% dN	1 %dN	1.0	1%dN
				Place in the list
				Digit number for results
				100,0

Test detailed procedures

Procedure Name test	Evaluation Method
BIOPHEN Plasminogen	Delta A per minute

Calculate Formula of raw data

Formula Type	Minimum	<input type="text"/>
User own Formula	<input type="text"/>	
Unit	mE/min	<input type="text"/>

Detailed test evaluation

		Determination of result		
Usual Values	<input type="text" value="75"/> % dN	at	<input type="text" value="150"/> %dN	Standard Concentrations
				% dN
Test used for Calibration	<input type="text"/>			120**
				100
				75
Evaluation	Regression lin/lin			50
Reference Curve	Measure calibration curve			25
Minimum raw value	<input type="text" value="0.0"/>			
Upper Extrapolation	<input type="text" value="2"/>	x	maxi concentration	
Lower Extrapolation	<input type="text" value="1"/>	x	mini concentration	
Calibrator	calibrator			
Authorized SD respectively to the reference curve	<input type="text" value="0"/>			%
Maximal number repeats	<input type="text" value="0"/>			

** or indicated concentration "C" for the calibrator, to be adjusted.

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