Manufactured By: HYPHEN BioMed

BIOPHEN Plasminogen ON STA-R

Adaptation of BIOPHEN Plasminogen ON STA-R

1. Reconstitution of BIOPHEN plasminogen(Ref A221502)reagents.

Chromogenic determination of plasminogen activity.

	NAME	Reconstitution	Stability	Stabilization in T°
R1	Activator (SK)	2.5 ml of distilled water (*)	1 month at 2-8°C * 7 days at room T° Do not freeze	30 mn before any use (**)
R2	Substrate	2.5 ml of distilled water (*)	1 month at 2-8°C * 7 days at room T° Do not freeze	30 mn before any use (**)

Reagents required but not provided:

Physiological saline (9 g/L NaCl)

Reconstitution:

(*) After reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature.

Stabilization of reagents:

(**) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate board before any use.

Storage of reagents:

Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Foot-note: Do not interchange the reagents from different lots.

NOTE: If necessary use microcups when running the assay, to improve management of the reagents volumes by the automate.

2. Determination of Plasminogen.

Preparation of the calibrators

- Calibration is performed with normal pooled citrated plasma with the assigned value of 100% Plasminogen. The assay includes a standard plasma dilution of 1:20 (managed by the instrument). By definition, this later dilution of the pool represents the 100% Plasminogen activity. The dynamic range is from 0 to 200% Plasminogen. The 200% Plasminogen activity is the 1:10 dilution of the plasma pool.

Or

-Calibration is performed with plasma calibrator commercially available, with a known Plasminogen concentration. The 1:20 dilution corresponds to the indicated Plasminogen concentration, and the 1:10 to twice this concentration.

Tested plasma:

In these conditions, the tested plasmas and controls are used undiluted and assayed at the 1:20 dilution (managed by the instrument).

A calibration curve must be carried out for each new batch of reagents.

Controls:

The control is performed with commercially available control plasmas, titrated for Plasminogen.

NAME	Reconstitution	Stability	Stabilization in T°
Calibration Biophen Plasma Calibrator (ref A222101)	1 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T°	30 minutes on STA board before any use (**)
Quality controls Biophen Normal Control (ref A223201) Biophen Abnormal Control (ref A223301)	1 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T°	30 minutes on STA board before any use (**)

Reconstitution: (*) Following reconstitution of calibrators or controls with distilled water, let them to stabilize for 30 minutes at room temperature.

Storage of reagents: (**) Take care of strictly respecting the 30 minutes temperature stabilization time for *calibrators* and *controls* at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C. Homogeneize before each use.

Foot-note: Do not freeze calibrators or controls.

Foot-note: A calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range.

3. Results:

- *CALIBRATION RANGE (0 to 200 %)*: The calibration curve (working range) is of the Lin (absorbance) – Lin (concentration) type.
- The values obtained for the patients and controls (tested at the standard 1:20 dilution) are directly calculated from the calibration curve.
- The results are expressed as % Plasminogen activity (100 % being equivalent to 1 Unit/ml)

4. Programming the STA-R analyzer:

Click on the icon **set up software** for the **manager program** and create the program according to:

Concentration working range : 0 -200% of Plasminogen

TESTS				Configuration				
Method								
IDENTIFICATION								
Abbreviation		BPLG		Last up date				
Name		BIOPHEN Plasminogen		Method		Colo-Kinetics OD /min		
SAMPLE				DILUENT				
Vol. μ l	Incu.	Dilution	Id.	Name			Stab. h	Continued
100 μ l	0 sec	1/20	PS	PHYSAL			24	
REAGENTS								
	Id.	Name		Stab. h	Vol. μ l	Incu. sec	Prec	Continued
Ra	R1	R1 PLG		168	100	180	<input type="checkbox"/>	
Rb							<input type="checkbox"/>	
Rc							<input type="checkbox"/>	
Rd	R2	R2 PLG		168	100			

Enter data concerning washing

	Id.	Name	Stab. h	Vol. μ l	Incu. sec	Prec.	Vial. ml	Min Vol ml
Diluent	PS	PHYSAL	24				5	1
Ra	R1	R1 PLG	168	100	180	<input type="checkbox"/>	3	0.35
Rb						<input type="checkbox"/>		
Rc						<input type="checkbox"/>		
Rd	R2	R2 PLG	168	100			3	0.35
Washing								
	Ra		Rb		Rc		Rd	
Before	No						No	
After	Special						Intensive	
Wash	12227							
Name	STA-DESORB U							
Stab. H	120							
Vial ml	15							
Vol in. ml	0.9							

Click on the icon **Result software** for the **manager program** and create the program according to:

TESTS		Configuration			
		Result			
METHOD Chro KINETICS OD/mn		RESULT			
Linearity Min.	0.95	Primary Units	%		
First Point	15	Corrector	1.00		
Second Point	90	Determination	<input checked="" type="radio"/> Single <input type="radio"/> Duplicate		
		Precision (%)	5.00	Validation (%)	
		Redilution	1/	<	>
		Condition	1/	>	<
<input checked="" type="checkbox"/> Rd Heating <input type="checkbox"/> X Stirring				Min.	Max.
				0	250

Write calibration data.

TESTS		Configuration			
		Calibration			
MODE :		Linear	Determination		
			<input checked="" type="checkbox"/> Duplicate <input type="checkbox"/> Single		
ASSAYS		Concentrations			
SCALE : <input checked="" type="checkbox"/> Linear <input type="checkbox"/> Log		<input checked="" type="checkbox"/> Linear <input type="checkbox"/> Inverse <input type="checkbox"/> Log			
	Id.	Key	Name	Stab. H	Dilution
Calibrator 1	PS	0	PHYSAL (ZEROCAL)	24	1/1
Calibrator 2	CAL	0	CAL	*	1/10
Calibrator 3	CAL	0	CAL	*	1/20
Calibrator 4	CAL	0	CAL	*	1/40
Calibrator 5	CAL	0	CAL	*	1/80
Ctrl, Niv. 1					
Ctrl, Niv. 2					

*according to the calibrator used

Click on the icon **Printout/Transmission software** for the **manager program** and create the program according to:

TESTS		Configuration				
		Printout /Trans.				
PARAMETERS						
	Unit	Conver Factors.	Print	Transmission Test number	Usual Values Min. Max.	
Main	%		<input checked="" type="checkbox"/>	0	70	120
Aux1	OD/min		<input checked="" type="checkbox"/>			
Aux2			<input type="checkbox"/>			
Aux3			<input type="checkbox"/>			
Printout Limits						
Min.	<input type="text" value="0"/>	Max.	<input type="text" value="250"/>			

Show the last page and enter data concerning the quality control.

TESTS		Configuration				
		Q.C				
LEVEL1						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol. ml
Control	BNC		BN CONT	8*	1*	0.2*
LEVEL 2						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control	BAC		BA CONT	8*	1*	0.2*
LEVEL3						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control						
Period : <input type="text"/> hours <input type="text"/> <input type="checkbox"/> Tests <input type="checkbox"/> Vial						

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