NOTE : reconstitution, preparation, calibration and controls, expression of results: usual cautions and realization according to the specific technical insert D.750.02/BI/1605.

Reagents are reconstituted as per the insert, well equilibrated, and preheated at 37°C

PHASE	CONDITION	OPERATION	PARAMETERS
1 st loading step	YES	Reagent loading (pos 18)	0µl
		Washing after reagent loading (pos 18)	0
		Diluent (pos DIL)	0µL
		Specimen (working dilution : prediluted 1 : 50 in prepared R3 buffer)	70µl
		Reagent (nacelle 2) Reagent R1 Ecarin	70µL
		Washing between specimen	2
		Optical ref. (cycle abs.) Diluent (pos DIL) Reagent (nacelle 2)	80µL 80µL
		Washing (at the end of loading)	3
		Incubation time (T1)	180 s
		2 nd loading step	YES
Specimen	0µL		
Washing between specimen	0		
Reagent (nacelle 3) Reagent R2 substrate	70µL		
Washing after reagent loading	3		
Mix time (T2)	120 s		
3 rd loading step	NO	Activation time (T3)	0 s
		Reagent (pos 17)	0µL
		Washing after reagent loading	0
		Mix time (T4)	0 s
Washing	NO	Activation time (T5)	0 s
		Reagent (pos 1)	0µL
		Reagent (pos 16)	0µL
Levels detection	NO	Cycles number	1
		Washing at the end of loading	1
Acquisition parameters	YES	Ramp	yes
		Inter-ramp interval (T6)	1s
		Delay before acquisition (T7)	5s
		Acquisition time (T8)	125s*
		Canal used	405nm
		Rotation speed	1200 rpm

*ie 120 seconds acquisition; could be reduced to 65 sec, but lower OD values expected.

Example of calibration curve obtained:

	Acquisition 120 sec
% FII	A405
200% (1 :25 dil)	1.58
100% (1 :50 dil)	0.86
50% (1 :100 dil)	0.45
25% (1 :200 dil)	0.23
12.5% (1 :400 dil)	0.12
r^2 (logA405-log conc)	0.999

Note: The stability at 2-8°C and at room temperature (18-25°C) claimed are validated on reconstituted vials kept closed, protected from any evaporation or contamination. Stabilities of reagents on each automate must be controlled, adjusted and validated under each exact laboratory working conditions.

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