## 1. Reconstitution of Hemoclot Thrombin Inhibitors reagents.

- Refer to each specific insert.
- The reagents are equilibrated at RT and then gently mixed and equilibrated at 37°C on the BCS prior to testing.

- # ACK002K presentation :

	NAME	Reconstitution	Stability*	T° Stabilization
R1	Normal Pooled Plasma	1 ml of distilled water*	Refer to the insert*	** 30 mn on board before any use
R2	Human Calcium Thrombin	1 ml of distilled water*		

- # ACK002L presentation :

	NAME	Reconstitution	Stability*	T° Stabilization
R1	Normal Pooled Plasma	2.5 ml of distilled water*	Refer to the insert*	** 30 mn on board before any use
R2	Human Calcium Thrombin	2.5 ml of distilled water*		

**\*The stability data claimed on the insert were obtained on reconstituted vials, kept closed, protected from any contamination or evaporation. Stability must be controlled, and can be adjusted and validated if required, according to the exact use conditions for each laboratory.**

**Reconstitution:** (\*) After reconstitution with distilled water, let the reagent to stabilize for 30 minutes at room temperature (18-25°C).

**Stabilization of reagents:** (\*\*) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate board before any use.

**Storage of reagents:** Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C, and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

If the reagents are kept on the automate board, take care and use reducers to limit as much as possible any evaporation of the reagents.

**Homogenize the reagents before each use.**

Any reagent of biological origin must be handled with all the required cautions, as being potentially infectious.

**Do not interchange the reagents from different lots.**

## 2. Reagents required but not provided:

- **DILUENT:** Owren Koller buffer (eg AAR003A/K) or Physiological Saline buffer (NaCl 9 g/L).

**The same diluent must be used for all the tests performed.**

- **CALIBRATION:** Internal Reference Calibrator titrated for the assayed Thrombin Inhibitor (DTI). Alternatively the following calibrators are available:

**For the assay of Hirudin:**

Plasma Hirudin Standard Low (Ref : ASC020K) using the low range calibration.

OR

Plasma Hirudin Standard High (Ref : ASC020L) using the high range calibration.

**For the assay of Argatroban®:**

Argatroban® Plasma Calibrator (Ref: ASC030K) using the low range calibration.

**For the assay of Dabigatran:**

Dabigatran Plasma Calibrator (Ref: A222801) using the low range calibration.

- **QUALITY CONTROL:** Internal Reference Quality Controls titrated for the assayed Thrombin Inhibitor (DTI). Alternatively the following controls are available:

**For the assay of Hirudin:** Plasma Hirudin Control (Ref : ASC025K)  
(C1 better for low range; C2 better for High range).

**For the assay of Argatroban®:** Argatroban® Control Plasma (Ref: ASC035K)

**For the assay of Dabigatran:** Dabigatran Control Plasma (Ref: A224701)

## 3. Calibrators and Controls:

**Cautions:** Refer to each specific insert for the **reconstitution** volume, the **stability** after reconstitution (indicated provided any contamination or evaporation is avoided. Stability can be adjusted according to the exact use conditions), and the **stabilization before use**.

**Refer to the exact concentrations indicated on the flyer for each lot.**

Nota : For lyophilized calibrators and controls, following reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature. It is recommended to run the calibration curve with a **freshly reconstituted calibrator**. It is necessary to let the reagent temperature to stabilize for at least 30 minutes onto the automate before any use. Take care avoiding any contamination or evaporation of the reagents. Stability can be adjusted according to the exact use conditions. Take care of putting up the specific caps back on the bottles before storing them at 2-8° C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

**Homogenize before each use.**

**Do not freeze calibrators and quality control plasmas.**

**Quality controls must be run regularly, and for each new batch of reagents, after an important maintenance of the instrument, or if measured values are not in compliance with the one expected for the method.**

**Performances may present slight variations according to the assayed DTI and the instrument used. Validate the expected values in the exact laboratory working conditions.**

**Check the validity of the series by including quality control plasmas at different levels in each one.**

**“LOW RANGE Concentration Curve”:**

- **Prepare the appropriate standard curve:**
  - **From about 0 to 2µg/ml Hirudin** (5 levels eg: 0, 0.5, 1.0, 1.5, 2.0 or “C” µg/ml, obtained by mixing appropriate volumes of the two provided “low range” calibrators, #ASC020K, according to the following table:

Hirudin	µg/ml	0	C:4	C:2	3C:4	C
Plasma at “C” (about 2 µg/ml) in µl		0	75	150	225	300
Normal Plasma (about 0µg/ml) in µl		300	225	150	75	0

If required, volumes can be proportionally increased to ensure sufficient volume for a good management by the automate.

- **From about 0 to 2µg/ml Argatroban®** (5 levels of ready to use calibrators, #ASC030K)
- **From about 0.05 to 0.50 µg/ml Dabigatran** (3 levels of ready to use calibrators, #A222801)
- **Refer to the exact concentration indicated on the flyer for each lot.**
- **The standards then have to be assayed at the 1:8 dilution (which is directly managed by the automate).**

**“HIGH RANGE Concentration Curve” (for Hirudin):**

- **Prepare a high range standard curve from about 0 to 5 µg/ml Hirudin** (5 levels eg: 0, 1.25, 2.50, 3.75, 5.0 or “C” µg/ml, obtained by mixing appropriate volumes of the two provided “high range” calibrators, #ASC020L, according to the following table:

Hirudin	µg/ml	0	C:4	C:2	3C:4	C
Plasma at “C” (about 5 µg/ml) in µl		0	75	150	225	300
Normal Plasma (about 0µg/ml) in µl		300	225	150	75	0

If required, volumes can be proportionally increased to ensure sufficient volume for a good management by the automate.

- **Refer to the exact concentration indicated on the flyer for each lot.**
- **The standards then have to be assayed at the 1:20 dilution (which is directly managed by the automate).**

#### 4. Results:

- The calibration curve is of the Lin (Clot time, in seconds), on ordinates - Lin ( $\mu\text{g/ml}$  assayed DTI) type (on abscissae).
- The values obtained for patients and controls are directly calculated from the calibration curve (when the standard dilution is used for the test).
- The results are expressed as  $\mu\text{g/ml}$  of assayed DTI.

The calibration curve is validated when linearity ( $r^2 \geq 0.98$ ), as well as measured control values (at least one in each series, better two levels), are in compliance.

A new calibration curve must be carried out better with each test series, and at least for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).

Clotting times and performances may present variations according to the assayed DTI, the instrument used, and to the clot detection sensitivity adjustment. Validate the expected values in the exact laboratory working conditions. Performances, as well as expected ranges, and values for each new lot of quality controls used, must then be confirmed (and adjusted if necessary) for each specific application in the exact laboratory working conditions and for the specific application.

**NOTE: To avoid any confusion, take care to create and identify appropriate reagents and specific procedure for each assayed DTI and concentration range.**

## 5. Programming the analyzer:

### Creation of reagents files

Return to definition of test. Go into reagents without data of lot. Create the reagent files.

<b>Name:</b> HIRUDIN R1 <b>Abbreviation:</b> HIR 1 clean when the reagent is modified <b>High</b> <input checked="" type="checkbox"/> New calibration with any new lot of reagents <b>Authorized lines</b> from 1 to 14 <b>Stability (*)</b> < 15C    > = 15C <b>Mix by jet</b> Intensity 1 Frequency No agitation <b>Reference numbers</b> add    remove    replace	<b>Name:</b> HIRUDIN R2 <b>Abbreviation:</b> HIR 2 clean when the reagent is modified <b>Special</b> <input checked="" type="checkbox"/> New calibration with any new lot of reagents <b>Authorized lines</b> from 1 to 14 <b>Stability (*)</b> < 15C    > = 15C <b>Mix by jet</b> Intensity 1 Frequency No agitation <b>Reference numbers</b> add    remove    replace
--	---

### Creation of calibration files

Create a table for each point of the calibration curve (eg : set of calibrators) (for the high and low range); according to the examples here below:

<b>Name:</b> HIRUDIN Low1 <b>Abbreviation:</b> HIR Low 1 clean when the reagent is modified <b>Low</b> <input checked="" type="checkbox"/> New calibration with any new lot of reagents <b>Authorized lines</b> from 1 to 14 <b>Stability (*)</b> < 15C    > = 15C <b>Mix by jet</b> Intensity 1 Frequency No agitation <b>Reference numbers</b> add    remove    replace	<b>Name:</b> HIRUDIN High1 <b>Abbreviation:</b> HIR High 1 clean when the reagent is modified <b>Low</b> <input checked="" type="checkbox"/> New calibration with any new lot of reagents <b>Authorized lines</b> from 1 to 14 <b>Stability (*)</b> < 15C    > = 15C <b>Mix by jet</b> Intensity 1 Frequency No agitation <b>Reference numbers</b> add    remove    replace
---	---

(\*) To be filled by the user

### Creation of Control files

Create quality control plasmas for the DTI to be assayed (based on same parameters as calibrators), preferably at least 2 appropriate levels, within the assay range.

(For information, if Hirudin Control Plasmas (#ASC025K) are used, C1 is more representative for the low range, and C2 for the high range).

### Creation of the procedure

Once the various reagents, calibration and controls files are created, the procedure file must be created.

Go to the "Definition tests, procedure" and enter the following parameters.

### **NOTE:**

**According to the tested DTI, and concentration range used for the calibration, two procedures must be created:**

**Procedure for the LOW RANGE (0 to 2 µg/ml Hirudin or 0 to 2 µg/ml Argatroban® or 0.05 to 0.5µg/ml Dabigatran), specific for each assayed inhibitor**

**OR**

**Procedure for the HIGH RANGE (0-5 µg/ml Hirudin)**

**LOW RANGE :**

General Information			
Former Change			
Test Procedure Number			
Test Procedure Name	TI Low		
Terminate measurement after	150 sec	Or	mA (0 = any measurement time)
	0		
Speed	normal		
Mixing	normal		
Wave-Lenght	405 nm		
Primary absorbance range	201-240 mA		
Duplicate n# for samples and controls	*	CV Authorized	<input type="text"/> %
Duplicate n# for calibration	*	CV Authorized	<input type="text"/> %
Duplicate n# of raw data summary	Arithmetic mean		
Dilution Factor	1 : 8	(1 = no dilution)	
Dilution Buffer			

Evaluation and Control Methods			
<b>Fixed Absorption</b>			
Lag Phase	<input type="text" value="7"/>	Sec	Threshold <input type="text" value="300"/> mA
<b>Angle</b>			
Lag	<input type="text" value="7"/>	Sec	Threshold <input type="text" value="300"/> mA
Window size	<input type="text" value="3"/>	Sec	
Minimum angle	<input type="text" value="15"/>	angle	
<b>Decrease</b>			
Lag	<input type="text" value="7"/>	Sec	
Max Absorbance Decrease	<input type="text" value="10"/>	%	
<b>Sinusoide</b>			
Lag Phase	<input type="text" value="7"/>	Sec	Threshold <input type="text" value="300"/> mA
Minimum Amplitude	<input type="text" value="100"/>	mA	Threshold Factor <input type="text" value="2.0"/>
Number of Steps	<input type="text" value="10"/>		Limit <input type="text" value="20"/>
<b>A MaxMin</b>			
Lag Phase	<input type="text" value="7"/>	Seconds	
Minimum/maximum Absorption	<input type="text" value="7200"/>	mA	

**\*to be filled and adjusted by the user. It is recommended to perform the tests in duplicate.**

## Pipetting Sequence

Test Procedure Number

Test Procedure Name

### Sampling cycles

Number of Cycles	Transmission Arm	Final Washing	Rotors/transfert Action	Time Window [sec]	
				Min	Max
1	Arm	normal	No mix./Measure	0	0
2	Reagent Arm	normal	Incubation	220	260
3	Reagent Arm	High	No mix./Measure	0	0
4	Reagent Arm	High	...		

Transm. Number	Mean/pipetting purpose	Reagents	Speed	Volume (µl)
1	1 Take air		Slow	20
	2 Take the sample		Slow	50
	3 In the external cup		Slow	50
2	1 Take air		Slow	20
	2 Take the Reagent	R1 Hirudin	Slow	100
	3 In the central cup		Slow	100
3	1 Take air		Slow	20
	2 Take the Reagent	R2 Hirudin	Slow	100
	3 In the Internal cup		Slow	100
	4** Take the Reagent	SCS Cleaner	Slow	0

**\*\*it can be suppressed if an error is generated on the instrument, or for economical or time gain.**

Test Definition					
<b>General Information</b>					
Former Change					
Name	TI Low			°Test N°	
Result Unit				PC ID t	
Behring	µg/ml	Conversion factor		Abbreviation	TI Low
User	µg/ml	1 µg/ml	1.0	µg/ml	Place in the liste
PC	µg/ml	1 µg/ml	1.0	µg/ml	...
					Digit number for results
					100.0

### Test detailed procedures

Test Procedure Name	Evaluation Method
TI Low	Delta A fix (Fixed absorption)

### Calculate Raw Data Formula

Formula Type	Minimum	
User own's Formula		
Raw data unit	mE/min or <b>SEC</b>	(or sec)

### Detailed test evaluation

<b>Reference range</b>	Result Determination		<b>Standard Concentrations</b>
	0.0* µg/ml	at 2.0 or 0.5) * µg/ml	
<b>Test used for Calibration</b>			
<b>Evaluation</b>	Regression Lin-Lin		
<b>Reference Curve</b>	Measure curve		
<b>Minimum raw value</b>	0.0		
<b>Upper Extrapolation</b>	1.0	x maxi concentration	
<b>Lower Extrapolation</b>	1.0	x mini concentration	
<b>Calibrator</b>	Hirudin calibrator low		
<b>Authorized SD respectively to the reference curve</b>	0	%	
<b>Maximal number repeats</b>	0		

## **HIGH RANGE:**

Suitable to assay samples in the expected range >2 to 5 µg/ml Hirudin (eg in ECC).

Proceed as previously, but with the following changes:

- Well identify a **new test / procedure name** (eg : TI HIGH) each time is it required
- The calibration curve must be in the range 0 to about 5 µg/ml (or “C” µg/ml) for Hirudin.
- In “general information”, **the dilution factor must be programmed at 1:20** (instead of 1:8 previously)
- The reference range is then from about 0 to 5µg/ml

---

**Note :** as an indication, the basic clotting time (CT) obtained for the « 0µg/ml » concentration (normal plasma) is expected near 30 sec ; and near 90 seconds (variable according to the reagents lots and test conditions) for 2µg/ml Hirudin or Argatroban for example, when a 1:8 dilution is used (low range).

The obtained CT remain similar when the High range procedure (1:20 dilution) is used (eg near 30 sec for 0µg/ml and near 90 seconds for 5µg/ml Hirudin)

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