Manufactured by Hyphen BioMed.

Hemoclot Thrombin Inhibitors on Sysmex CA1500

Determination of Hirudin (and other DTIs) using the CA1500 instrument

PROVISORY PROPOSAL

1. Reconstitution of Hemoclot Thrombin Inhibitors reagents

- Refer to each specific insert.
- The reagents are equilibrated at room temperature and then gently mixed and transferred onto the CA1500 to equilibrate to 15°C prior to testing.

• # ACK002K presentation :

	NAME	Reconstitution	Stability*	T° Stabilization
R1	Normal Pooled Plasma	1 ml of distilled water*	Refer to the insert*	** 30 mn on board before any use
R2	Human Calcium Thrombin	1 ml of distilled water*		

• # ACK002L presentation :

	NAME	Reconstitution	Stability*	T° Stabilization
R1	Normal Pooled Plasma	2.5 ml of distilled water*	Refer to the insert*	** 30 mn on board before any use
R2	Human Calcium Thrombin	2.5 ml of distilled water*		

*The stability data claimed on the insert were obtained on reconstituted vials, kept closed, protected from any contamination or evaporation. Stability must be controlled, and can be adjusted and validated if required, according to the exact use conditions for each laboratory.

Reconstitution: (*) After reconstitution with distilled water, let the reagent to stabilize for 30 minutes at room temperature (18-25°C).

Stabilization of reagents: (**) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate board before any use.

Storage of reagents: Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C, and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

If the reagents are kept on the automate board, take care and use reducers to limit as much as possible any evaporation of the reagents.

Homogenize the reagents before each use.

Any reagent of biological origin must be handled with all the required cautions, as being potentially infectious.

Do not interchange the reagents from different lots.

2. Reagents required but not provided:

- **DILUENT:** Owren Koller buffer (eg AAR003A/K) or Physiological Saline buffer (NaCl 9 g/L).

The same diluent must be used for all the tests performed.

- **CALIBRATION:** Internal Reference Calibrator titrated for the assayed Thrombin Inhibitor (DTI). Alternatively the following calibrators are available:

For the assay of Hirudin:

Plasma Hirudin Standard Low (Ref : ASC020K) using the low range calibration.

OR

Plasma Hirudin Standard High (Ref : ASC020L) using the high range calibration.

For the assay of Argatroban®:

Argatroban® Plasma Calibrator (Ref: ASC030K) using the low range calibration.

For the assay of Dabigatran:

Dabigatran Plasma Calibrator (Ref: A222801) using the low range calibration.

- **QUALITY CONTROL:** Internal Reference Quality Controls titrated for the assayed Thrombin Inhibitor (DTI). Alternatively the following controls are available:

For the assay of Hirudin: Plasma Hirudin Control (Ref : ASC025K)
(C1 better for low range; C2 better for High range).

For the assay of Argatroban®: Argatroban® Control Plasma (Ref: ASC035K)

For the assay of Dabigatran: Dabigatran Control Plasma (Ref: A224701)

3. Calibrators and Controls:

Cautions: Refer to each specific insert for the **reconstitution** volume, the **stability** after reconstitution (indicated provided any contamination or evaporation is avoided. Stability can be adjusted according to the exact use conditions), and the **stabilization before use**.

Refer to the exact concentrations indicated on the flyer for each lot.

Nota : For lyophilized calibrators and controls, following reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature. It is recommended to run the calibration curve with a **freshly reconstituted calibrator**. It is necessary to let the reagent temperature to stabilize for at least 30 minutes onto the automate before any use. Take care avoiding any contamination or evaporation of the reagents. Stability can be adjusted according to the exact use conditions. Take care of putting up the specific caps back on the bottles before storing them at 2-8° C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Homogenize before each use.

Do not freeze calibrators and quality control plasmas.

Quality controls must be run regularly, and for each new batch of reagents, after an important maintenance of the instrument, or if measured values are not in compliance with the one expected for the method.

Performances may present slight variations according to the assayed DTI and the instrument used. Validate the expected values in the exact laboratory working conditions.

Check the validity of the series by including quality control plasmas at different levels in each one.

“LOW RANGE Concentration Curve”:

- Prepare the appropriate standard curve:
 - From about 0 to 2µg/ml Hirudin (5 levels eg: 0, 0.5, 1.0, 1.5, 2.0 or “C” µg/ml, obtained by mixing appropriate volumes of the two provided “low range” calibrators, #ASC020K, according to the following table:

Hirudin	µg/ml	0	C:4	C:2	3C:4	C
Plasma at “C” (about 2 µg/ml) in µl		0	75	150	225	300
Normal Plasma (about 0µg/ml) in µl		300	225	150	75	0

If required, volumes can be proportionally increased to ensure sufficient volume for a good management by the automate.

- From about 0 to 2µg/ml Argatroban® (5 levels of ready to use calibrators, #ASC030K)
- From about 0.05 to 0.50 µg/ml Dabigatran (3 levels of ready to use calibrators, #A222801)
- Refer to the exact concentration indicated on the flyer for each lot.
- The standards are then loaded “undiluted” and have to be assayed in duplicate at the 1:8 dilution (which is directly managed by the automate, as defined in the LOW RANGE program used).

Note: Controls and Samples are loaded “undiluted” and tested at the final 1:8 dilution is directly managed by the instrument (defined in the LOW RANGE program : as first step: 30µl sample + 90µl diluent and second step 25µl + 25µl diluent)

“HIGH RANGE Concentration Curve” (for Hirudin):

- Prepare a high range standard curve **from about 0 to 5 µg/ml Hirudin** (5 levels eg: 0, 1.25, 2.50, 3.75, 5.0 or “C” µg/ml, obtained by mixing appropriate volumes of the two provided “high range” calibrators, #ASC020L, according to the following table:

Hirudin	µg/ml	0	C:4	C:2	3C:4	C
Plasma at “C” (about 5 µg/ml) in µl		0	75	150	225	300
Normal Plasma (about 0µg/ml) in µl		300	225	150	75	0

If required, volumes can be proportionally increased to ensure sufficient volume for a good management by the STAR.

- Refer to the exact concentration indicated on the flyer for each lot.
- The standards are then loaded “undiluted” and have to be assayed in duplicate at the 1:20 dilution (which is directly managed by the automate, as defined in the HIGH RANGE program used).

Note: Controls and Samples are loaded “undiluted” and tested at the final 1:20 dilution which is directly managed by the instrument (defined in the HIGH RANGE program, as 1:10 dilution in a first step (12µl sample + 108µl diluent), and subsequent 1:2 dilution).

4. Results:

- The calibration curve is of the Lin (Clot time, in seconds), on ordinates - Lin ($\mu\text{g/ml}$ assayed DTI) type (on abscissae).
- The values obtained for patients and controls are directly calculated from the calibration curve (when the standard dilution is used for the test).
- The results are expressed as $\mu\text{g/ml}$ of assayed DTI.

The calibration curve is validated when linearity ($r^2 \geq 0.98$), as well as measured control values (at least one in each series, better two levels), are in compliance.

A new calibration curve must be carried out better with each test series, and at least for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).

Clotting times and performances may present variations according to the assayed DTI, the instrument used, and to the clot detection sensitivity adjustment. Validate the expected values in the exact laboratory working conditions. Performances, as well as expected ranges, and values for each new lot of quality controls used, must then be confirmed (and adjusted if necessary) for each specific application in the exact laboratory working conditions and for the specific application.

NOTE: To avoid any confusion, take care to create and identify appropriate reagents and specific procedure for each assayed DTI and concentration range.

Note: R1 Plasma pool needs to be programmed as a “Factor deficient plasma “ and R2 Thrombin as “reagent”.

5. PROGRAMMING THE ANALYZER FOR THE LOW RANGE :

This allows measuring concentrations in the range of about 0.1 to 2µg/ml Hirudin or Argatroban®, or from 0.05 to 0.50 µg/ml Dabigatran.

All calibrators, controls and patient plasma are loaded “pure” (undiluted); the 1:8 dilution being directly managed by the the automate.

Click on the window Set up software for the “manager program” and create the program according to: “Clotting for TT” or OTHER

Change the data with the data indicated on the table below

Use the program and create the name for the reagents, Calibration plasmas and control plasmas

System		Ready		
Parameter : TI LOW Para Code (*)				
Sample.Vol		30	µl	
Diluent Vol	Owren or PHYSAL	90	µl	
Rinse	None			
Second Dilution Diluent Vol	Owren or PHYSAL	25µl	25µl	7 8 9
Rinse	None			
Factor Plasma	R1 Plasma 100µl			
Rinse (Pre/Post)	None			
First reagent	R2Thrombin 100µl 60sec			
Push out solution	No	0µl		4 5 6
Rinse (Pre/Post)	None	x0	Cleanl x1	1 2 3
Second reagent	None	0µl	0 sec	0 Enter
Push out solution	No	0µl		
Rinse (Pre/Post)	None	x0	None x0	
Third reagent	None	0µl	0 sec	
Push out solution	No	0µl		
Rinse (Pre/Post)	None	x0	None x0	
Detector	Clot “for T T or other...”			
Sens	Low sens			
Maximum time	260 sec			
Select Test	Set Testname	Special	↑	↓
				quit

(*) User defined.

- **Replication; Duplicate is recommended**
- **Detector settings;** 50% or replicates = 2
- **Group settings;** One.
- **Parameter settings;** Calculated parameter 1 – Linear/Linear Regression

Standard curve

- **Sampler/Holder;** user defined
- **Stds** –All in duplicate.

6. PROGRAMMING THE ANALYZER FOR THE HIGH RANGE :

This calibration allows measuring concentrations in the range of about 2 to 5µg/ml Hirudin.

All calibrators, controls and patient plasma are loaded “pure” (undiluted); the 1:20 dilution being directly managed by the automate.

Click on the window Set up software for the “manager program” and create the program according to: “ Clotting for TT” or other. Change the data with the data indicated on the table below. Use the program and create the name for the reagents, Calibration plasmas and control plasmas

System		Ready		
Parameter : TI HIGH Para Code (*)				
Sample.Vol		12	µl	
Diluent Vol	Owren or PHYSAL	108	µl	
Rinse	None			
Second Dilution		25µl		
Diluent Vol	Owren or Physal	25µl		
Rinse	None			
Factor Plasma	R1 Plasma 100µl			
Rinse (Pre/Post)	None			
First reagent	R2Thrombin 100µl 60sec			
Push out solution	No	0µl		
Rinse (Pre/Post)	None x0	Cleanl x1		
Second reagent	None	0µl	0 sec	
Push out solution	No	0µl		
Rinse (Pre/Post)	None x0	None x0		
Third reagent	None	0µl	0 sec	
Push out solution	No	0µl		
Rinse (Pre/Post)	None x0	None x0		
Detector	Clot “for T T or other...”			
Sens	Low sens			
Maximum time	260 sec			
Select Test	Set Testname	Special	↑	↓
				quit

(*)User defined

- **Replication; Duplicate is recommended**
- **Detector settings;** 50% or replicates =2
- **Group settings;** One.
- **Parameter settings;** Calculated parameter 1 – Linear/Linear Regression

Standard curve

- **Sampler/Holder;** user defined
- **Stds** –All in duplicate.

Fabricant: HYPHEN BioMed

Hemoclot Thrombin Inhibitors sur SYSMEX CA1500

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