Manufactured By: HYPHEN BioMed

HEMOCLOT PS

ON STA-R

Assay of PS clotting activity on STAR

1. Reconstitution of Hemoclot PS (# ACK041K) reagents.

	NAME	Reconstitution	Stability*	T° Stabilization
R1	PS deficient plasma	1 ml of distilled water *	Refer to the insert	** 30 mn on board before any use
R2	Activator	1 ml of distilled water *		

*The stability data claimed on the insert were obtained on reconstituted vials, kept closed, protected from and provided any contamination or evaporation is avoided. Stability must be controlled, and can be adjusted and validated if required, according to the exact use conditions for each laboratory.

Reconstitution: (*) After reconstitution with distilled water, let the reagent to stabilize for at least 15 minutes at room temperature (18-25°C).

Stabilization of reagents: (**) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate board before any use.

Storage of reagents: Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C, and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

If the reagents are kept on the automate board, take care and use reducers to limit as much as possible any evaporation of the reagents.

Homogenize the reagents before each use.

Any reagent of biological origin must be handled with all the required cautions, as being potentially infectious.

Do not interchange the reagents from different lots.

Reagents required but not provided:

- Plasma Calibrator titrated for PS (eg: Biophen Plasma Calibrator #A222101)
- Normal and Abnormal quality control plasmas titrated for PS (ex: BIOPHEN Normal Control Plasma (#A223201) and BIOPHEN Abnormal Control Plasma (#A223301)).
- Distilled water, preferentially sterile.
- Calcium chloride 0.025M (eg: # AAR001A/AAR001K).
- Imidazole buffer (eg: # AAR021A/AAR021K/AAR021L);

Note: using microcups can help the STAR to better manage the reagents volumes.

2. Preparation of the calibration curve and controls/samples.

• Calibration curve:

- Calibration is performed with normal pooled citrated plasma with the assigned value of 100% PS. The assay includes a standard plasma dilution of **1:10 (managed by the automate)**. By definition, this later dilution of the pool represents the **100%** PS activity. If the standard is less than 100%, it is possible to adjust its concentration by using a lower dilution factor.

Or

-Calibration is performed with a commercially available plasma calibrator, with a known PS Concentration (eg **Biophen Plasma Calibrator #A222101, to be titrated, in progress**). The **1:10 (managed by the automate)** corresponds to the indicated PS concentration ("C%").

Nota: If the standard is less than 100%, it is possible to adjust its concentration by using a lower dilution factor.

• Tested plasma and controls:

In these conditions, tested plasmas and controls are loaded "undiluted" and assayed at the 1:10 dilution (managed by the STAR). Samples expected at a normal concentration (>100%, can be assayed at the 1:20 dilution (test dependent).

Controls: The control is performed with commercially available control plasmas, titrated for PS.

Various control plasma are available (**to be titrated, in progress**): **Biophen Normal Control Plasma (#A223201)**, and **Biophen Abnormal Control Plasma (#A223301)**

Nota : For lyophilized calibrators and controls, following reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature. It is recommended to run the calibration curve with a freshly reconstituted calibrator. It is necessary to let the reagent temperature to stabilize for at least 30 minutes onto the automate before any use. Take care avoiding any contamination or evaporation of the reagents. Stability can be adjusted according to the exact use conditions.

Homogenize before each use.

Do not freeze calibrators and quality control plasmas.

Quality controls must be run regularly, and for each new batch of reagents, after an important maintenance of the instrument, or if measured values are not in compliance with the one expected for the method.

Performances may present slight variations according to the instrument used. Validate the expected values in the exact laboratory working conditions.

Check the validity of the series by including quality control plasmas at different levels in each one.

3. Results:

- The calibration curve (working range) is of the Log (CT, in seconds) – Lin (PS concentration, in %) type.
- The values obtained for patients and controls are directly calculated from the calibration curve (when the standard 1:10 dilution is used for the test).
- The results are expressed as % PS activity.

The calibration curve is validated when linearity (r²), as well as measured control values, are in compliance.

A new calibration curve must be carried out for each new batch of reagents, after each important maintenance of the instrument, or when measured values for controls are out of the acceptance range for the method (after checking all other parameters for the system).

Clotting times and performances may present variations according to the instrument used, and to the clot detection sensitivity adjustment. Validate the expected values in the laboratory working conditions. Performances, as well as values for each new lot of quality controls used, must then be confirmed (and adjusted if necessary) in the laboratory working conditions.

4. Programming the STA-R analyser:

Click on the icon **set up software** for the **manager program** and create the program according to:

Allows measuring concentrations in the range of about 10-100% PS (recommended).

TESTS				Configuration			
Method							
IDENTIFICATION							
Abbreviation			HEMPS		Last up date		
Name		HEMPS			Method	CLOT-BASED	
SAMPLE				DILUENT			
Vol. µl	Incu.	Dilution	Id.	Name		Stab. h	Continued
50µl	0 sec	1/10	IM	Imidazole		(*)	
REAGENTS							
	Id.	Name		Stab. h	Vol. µl	Incu. sec	Prec.
Ra	R1	R1 DEF PS		(*)	50	120	<input type="checkbox"/>
Rb	R2	R2 Activator		(*)	50	180	<input checked="" type="checkbox"/>
Rc							<input type="checkbox"/>
Rd	CaCl2	Calcium 0.025M		(*)	100		

Enter data concerning washing

REAGENTS								
	Id.	Name	Stab. h	Vol. µl	Incu. sec	Prec.	Vial. ml	Min.Vol ml
Diluent	IM	Imidazole	(*)				25	5
Ra	R1	R1 DEF PS	(*)	50	120	<input type="checkbox"/>	1	0.2
Rb	R2	R2 Activator	(*)	50	180	<input checked="" type="checkbox"/>	1	0.2
Rc						<input type="checkbox"/>		
Rd	R3	Calcium 0.025M	(*)	100			25	5
Washing								
	Ra		Rb		Rc		Rd	
Before	Non		Non				Special	
After	Normal		Normal				Normal	
Wash							12227	
Name							STA- DESORB U	
Stab. H							120	
Vial ml							15	
Vol in. ml							0.9	

(*) To be filled by the user, according to the instructions of the insert, and to the exact use conditions.

Click on the icon **Result software** for the **manager program** and create the program according to:

TESTS		Configuration	
		Result	
METHOD		RESULT	
Min. Time	35	Primary Unit	%
Max Time	200	Correction factor	1.00
Mean Time	150		
Clot	<u>Normal</u>	Single/Duplicate	
		<input checked="" type="radio"/> Single*	
		<input type="radio"/> Duplicate	
		Precision (%)	10.0*
		Redilution	Condition
		1/	<
		1/	>
<input checked="" type="checkbox"/> Rd Heating		VALIDATION %	
<input checked="" type="checkbox"/> Stirring			
		Min.	Max.
		10*	100*

*Or user adjusted. Testing in duplicate is recommended. For better accuracy, sample redilution is recommended for <10% or >100% PS

Write calibration data.

TESTS		Configuration	
		Calibration	
MODE :		Linear regression	
		Determination	
		<input checked="" type="checkbox"/> Single * Duplicate	
ASSAYS		Concentrations	
SCALE :	<input type="checkbox"/> Linear <input checked="" type="checkbox"/> Log	<input checked="" type="checkbox"/> Linear <input type="checkbox"/> Inverse <input type="checkbox"/> Log	
	Id.	Key	Name
Std 1	CAL		CAL
Std 2	CAL		CAL
Std 3	CAL		CAL
Std 4	CAL		CAL
Std 6	IM		Imidazole
Ctrl, Niv. 1			
Ctrl, Niv. 2			

*or duplicate, user adjusted. **According to the calibrator used (homogenize before each use).

Note: it is recommended to run the assay alone on the machine. the instrument can sometimes have difficulties to manage accurate incubation times; in this case it is recommended to run a limited number of samples at the same time so that the instrument can satisfiably manage the test.

Click on the icon **Printout/Transmission software** for the **manager program** and create the program according to:

TESTS		Configuration				
		Printout/Trans.				
PARAMETERS						
	Units	Factor Convers	Print	Transmission Test numbers	Usual Values Min. Max.	
Main	%		v	0	60(*)	140(*)
Aux1	SEC		v	0		
Aux2						
Aux3						
Printout Limits						
Min.	<input type="text" value="0"/>	Max.	<input type="text" value="200"/>			

(*) to be adjusted by the user.

Show the last page and enter data concerning the quality control.

TESTS		Configuration				
		Q.C				
LEVEL 1						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control	BNC		BN CONT	*	*	*
LEVEL 2						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control	BAC		BA CONT	*	*	*
LEVEL 3						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control				*	*	*
Period : <input type="text"/> hours <input type="text"/> Tests <input type="text"/> Vial						

*According to the controls used (homogenize before each use).

INDICATIVELY: CT (0% PS) = about 50 ± 10sec; CT (100%PS) = about 90 ± 120sec

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