Manufactured By: HYPHEN BioMed

BIOPHEN HEMOCLOT VII-X ON STA

Adaptation of BIOPHEN HEMOCLOT VII-X ON STA

1. Reconstitution of BIOPHEN HEMOCLOT VII-X (Ref ACK051K) reagent.

Determination of factors VII-X activity with a clotting method

NAME	Reconstitution	Stability	Stabilization at room Temperature
Hemoclot VII-X Reagent	1 ml Distilled water (*)	12 hours 2-8°C 8 hours at room T°	30 mn before any use (**)

Reconstitution:

(*) After reconstitution with distilled water, let the factor VII-X reagent to stabilize 15 minutes at room temperature.

Conservation of reagents:

Take care of putting up the specific caps back on the bottles before storing them at 2°-8°C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Stabilization of reagents:

(**) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate before any use.

Preparation of the calibrators and controls

Calibration curve:

- Prepare 1 ml of normal human pooled plasma diluted 1:10 in Owren Koller buffer (by definition the normal citrated plasma pool corresponds to a concentration of 100% of Factors VII + X).

Or

- The calibration curve can also be established with the BIOPHEN Plasma Calibrator (# A222101), using the factor VII-X activity indicated for the lot used.

Name	Reconstitution	Stability	Stabilization in T°
CALIBRATION Biophen Plasma Calibrator (ref A222101)	1 ml of distilled water(*)	24 hours at 2° - 8°C 8 hours at Room T°	30 minutes on STA board before any use (**)

Controls

Nom	Reconstitution	Stability	Stabilization in T°
CONTROLES DE QUALITE <u>Biophen Normal Control</u> (ref <u>A223201</u>) <u>Biophen Abnormal Control</u> (ref <u>A223301</u>)	1 ml of distilled water(*)	24 hours at 2° - 8° C 8 hours at Room T°	30 minutes <i>on STA board</i> before any use (**)

Reconstitution:

(*) After reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature.

Stabilization of reagents:

(**) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate before any use.

Conservation of reagents:

Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Do not freeze calibrators and controls

Nota : A calibration curve must be carried out of each new batch of Hemoclot Factor VII-X reagents

NOTA: Plasmas used for the pool and Thrombin preparation were tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions as being potentially infectious.

Results:

The calibration curve is of the Log (% activity VII-X) - Log (Clot time) type.

The values obtained for the patients and controls are directly calculated from the calibration curve.

The results are expressed in % activity, the 100% activity being the activity in reference normal pool plasma.

5. Programming of the STA analyser

Click on the icon **set up software** for the **manager program** and create the program according to:

Identification									
Abbreviation		Name			Method			Date	
VII-X		FACT VII-X			CLOT-BASED				
Sample				Diluent					
Volume	Incu	Dil.	Id.	Name		vial.	stab	Min.Vol	
50µl	0 sec	1/10	11361	OWREN KOLLER		15 ml	72h	0.9 ml	
Reagents									
Id.	Name	Incu sec	Vol. µl	Vial ml	Stab h	Min. Vol ml	Washing		
							Before	After	
Ra VII-X	FACT VII-X	240	50	1	8	0.2	no	Normal	
Rb									
Rc									
Rd THR	THROMBO		75	5	24	0.5	no	Intensive	
Analysis		Result				Validation (%)			
Min time	10 sec	Primary Units		%		Standard	min	max	
Max time	300 sec	Corrector		1.00			0	150	
Mean time	40 sec	Single/Duplicate		single					
Rd Heating	Oui	Precision		5%					
Stirring	Oui	Redil. Condition							
Clot type	normal	<							
		>							

Enter data concerning washing.

Washing	Ra	Rb	Rc	Rd
Before	No			No
After	Normal			Intensive
Wash				
Name				
Stab. H				
Vial ml				
Vol in. ml				

Show the second page and write calibration data.

CONFIGURATION – Page 2 / 3

Calibration										
Mode graph – linear regression						Scale				
Calibrators				Vial ml	Stab. H	Min. Vol ml	Raw Data Concentration Single/Duplicate	Log Log Single		
Id.	Key	Name	Dil.							
ETVII+X		ETVII+X	1/10	1	8	0.2				
ETVII+X		ETVII+X	1/20	1	8	0.2				
ETVII+X		ETVII+X	1/40	1	8	0.2				
ETVII+X		ETVII+X	1/80	1	8	0.2				
ETVII+X		ETVII+X	1/160	1	8	0.2				
Control										
	Id.	Key	Name	Vial. ml	Stab. H	Min. Vol. ml				
Level 1										
Level 2										

Go to the last page and enter data concerning the quality control.

The STA can manage 3 levels of quality control. For example: the following configuration is possible.

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