Adaptation of BIOPHEN HEMOCLOT VII-X ON STA-R

1. Reconstitution of BIOPHEN HEMOCLOT VII-X (Ref ACK051K) reagent.

Determination of factors VII-X. activity with a clotting method

NAME	Reconstitution	Stability	Stabilization at room Temperature
Hemoclot VII-X Reagent	1 ml Distilled water(*)	12 hours 2-8°C 8 hours at room T°	30 mn before any use (**)

Reconstitution:

(*) After reconstitution with distilled water, let the factor VII-X reagent to stabilize 15 minutes at room temperature.

Conservation of reagents:

Take care of putting up the specific caps back on the bottles before storing them at 2°-8°C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Stabilization of reagents:

(**) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate before any use.

Preparation of the calibrators and controls

Calibration curve:

- Prepare 1 ml of normal human pooled plasma diluted 1:10 in Owren Koller buffer (by definition the normal citrated plasma pool corresponds to a concentration of 100% of Factors VII + X).

Or

- The calibration curve can also be established with the BIOPHEN Plasma Calibrator (# A222101), using the factor VII-X activity indicated for the lot used.

Name	Reconstitution	Stability	Stabilization in T°
CALIBRATION Biophen Plasma Calibrator (ref A222101)	1 ml of distilled water(*)	24 hours at 2° - 8°C 8 hours at Room T°	30 minutes <i>on STA board</i> before any use (**)

Controls

Nom	Reconstitution	Stability	Stabilization in T°
CONTROLES DE QUALITE <u>Biophen Normal Control</u> (ref A223201) <u>Biophen Abnormal Control</u> (ref A223301)	1 ml of distilled water(*)	24 hours at 2° - 8° C 8 hours at Room T°	30 minutes <i>on STA board</i> before any use (**)

Reconstitution:

(*) After reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature.

Stabilization of reagents:

(**) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate before any use.

Conservation of reagents:

Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Do not freeze calibrators and controls

Nota : A calibration curve must be carried out of each new batch of Hemoclot Factor VII-X reagents

NOTA: Plasmas used for the pool and Thrombin preparation were tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions as being potentially infectious.

Results:

The calibration curve is of the Log (% activity VII-X) - Log (Clot time) type.

The values obtained for the patients and controls are directly calculated from the calibration curve.

The results are expressed in % activity, the 100% activity being the activity in reference normal pool plasma.

5. Programming of the STA R analyser

Click on the icon **set up software** for the **manager program** and create the program according to:

TESTS				Configuration			
Method							
IDENTIFICATION							
Abréviation			VII+X		Last up date.		
Name			BIOPHEN VII+X		Method		
			CLOT-BASED				
SAMPLE			DILUENT				
Vol. µl	Incu.	Dilution	Id.	Name	Stab. h	Min. Vol.	
50 µl	0 sec	1/10	11361	OWREN KOLLER	72	0.90	
REAGENTS							
	Id.	Name	Stab. h	Vol. µl	Incu. sec	Prec.	Continued
Ra	VII+X	VII+X	8	50	240	<input type="checkbox"/>	
Rb						<input type="checkbox"/>	
Rc						<input type="checkbox"/>	
THR	THROMBO	THR	24	75			

Enter data concerning washing

	Id.	Nom	Stab. h	Vol. µl	Incu. sec	Prec.	Vial. ml	Min Vol. ml
Diluant	11361	OWREN-KOLLER	72				15	0.9
Ra	VII+X	VII+X	8	50	240	<input type="checkbox"/>	1	0.2
Rb						<input type="checkbox"/>		
Rc						<input type="checkbox"/>		
Rd	THR	THROMBO	24	75			6	0.5
Washing								
	Ra		Rb		Rc		Rd	
Before	No						No	
After	Normal						Intensif	
Wash								
Name								
Stab. H								
Vial ml								

Click on the icon **Result software** for the **manager program** and create the program according to:

TESTS		Configuration	
I. Result			
METHOD: CLOT-BASED		RESULT	
Min. Time	10	Primary Units	%
Max Time	300	Corrector	1.00
Mean Time	30		
		Determination	
Clot	<input checked="" type="radio"/> Normal	<input checked="" type="radio"/> Single	
	<input type="radio"/> weak	<input type="radio"/> Duplicate	
	<input type="radio"/> Mean		
<input checked="" type="checkbox"/> Rd Heating		Precision (%)	5.00
<input checked="" type="checkbox"/> Stirring		Redilution	Condition
		1/	<
		1/	>
		VALIDATION (%)	
		Min.	Max.
		0	150

Write calibration data.

TESTS		Configuration	
II. Calibration			
MODE : Linear		Détermination	
		<input checked="" type="checkbox"/> Duplicate <input type="checkbox"/> Single	
SCALE		Concentrations	
RAW DATA : <input type="checkbox"/> Linear <input checked="" type="checkbox"/> Log		<input type="checkbox"/> Linear <input type="checkbox"/> Inverse <input checked="" type="checkbox"/> Log	
	Id.	Key	Name
Calibrator 1	BPC		BPCALVII+X
Calibrator 2	BPC		BPCALVII+X
Calibrator 3	BPC		BPCALVII+X
Calibrator 4	BPC		BPCALVII+X
Calibrator 5	BPC		BPCALVII+X
	Stab. H	Dilution	
	8	1/10	
	8	1/20	
	8	1/40	
	8	1/80	
	8	1/160	
Ctrl, Niv. 1			
Ctrl, Niv. 2			

Click on the icon **Printout/Transmission software** for the **manager program** and create the program according to:

TESTS							Configuration
III.		Printout/Trans.					
PARAMETERS							
	Unit	Factors. Conver	Print	Transmission Test number	Usual Values Min. Max.		
Main	%		<input checked="" type="checkbox"/>	0	70	150	
Aux1	SEC		<input checked="" type="checkbox"/>	0			
Aux2			<input type="checkbox"/>				
Aux3			<input type="checkbox"/>				
Printout Limits							
Min.	<input type="text" value="0"/>	Max.	<input type="text" value="150"/>				

Show the last page and enter data concerning the quality control.

TESTS							Configuration
							Q.C
LEVEL1							
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml	
Control	BNC		BN CONT	8	1	0.2	
LEVEL 2							
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml	
Control	BAC		BA CONT	8	1	0.2	
LEVEL3							
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml	
Control							
Period : <input type="text"/> hours <input type="text"/> <input type="checkbox"/> Tests <input type="checkbox"/> Vial							

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