**Adaptation of BIOPHEN HEMOCLOT V-L on STA**

**1. Reconstitution of BIOPHEN HEMOCLOT V-L (Ref ACK061K) reagent.  
(Ref ACK061L)**

Determination of factor V-Leiden. activity with a clotting method

<b>NAME</b>	<b>Reconstitution</b>	<b>Stability</b>	<b>Stabilization at the room Temperature</b>
<b>R1</b>	4ml distilled water (*)	8 hours at room T° 24hours at 2-8 °C 1 month at -20°C	30 mn before any use (**)
<b>R2A</b>	1ml distilled water (*)	8 hours at room T° 24hours at 2-8 °C 1 month at -20°C	30 mn before any use (**)
<b>R2B</b>	1ml distilled water (*)	8 hours at room T° 24hours at 2-8 °C 1 month at -20°C	30 mn before any use (**)
<b>R3</b>	4 ml Distilled water (*)	8 hours at room T° 24hours at 2-8 °C 1 month at -20°C	30 mn before any use (**)

**Reconstitution:**

(\*) Following reconstitution with distilled water, let the reagents to stabilize for 30 minutes at room temperature.

**Conservation of reagents:**

Take care of putting up the specific caps back on the bottles before storing them at 2°-8°C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

**Stabilization of reagents:**

(\*\*) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate before any use.

## 2. Determination of Factor V-Leiden activity

NAME	Reconstitution	Stability	Stabilization in T°
Quality controls Biophen Normal Control (ref A223201)	1 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T° Do not freeze	30 minutes <i>on STA board</i> before any use (**)
Biophen ACT PCr Control. (ref A223405)	0.5 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T° Do not freeze	30 minutes <i>on STA board</i> before any use (**)

### Reconstitution:

(\*) Following reconstitution controls with distilled water, let them to stabilize for 30 minutes at room temperature.

### Conservation of reagents:

(\*\*) Following care of strictly respecting the 30 minutes temperature stabilization time for controls at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C. Homogeneize before each use.

**Foot-note: Do not freeze controls.**

## 3. Results:

- The clotting times obtained for patients plasmas and controls with (CT2) or without (CT1) APC, are expressed in seconds, for each patient, the ratio CT2/CT1 between both clotting times must be calculated.
- The ratio of the Clotting Times obtained with or without APC, CT2/CT1, allows measuring the sensitivity of Factor V, in the tested specimen, to the action of Activated Protein C.
- Normal plasma, containing normal Factor V, is sensitive to this action, and yields a ratio:  $CT2/CT1 \geq 2.00$
- Plasmas from patients carrying the R506Q mutation of Factor V, i.e. containing Factor V Leiden, yield a ratio:  $CT2/CT1 \leq 1.80$ .
- Molecular biology allows confirming the diagnosis, and classifying patients as heterozygous or homozygous.

**NOTA: For a same reagent lot and a same tested specimen, CT1 and CT2 may present variations according to the instrument used, and to the clot detection sensitivity adjustment. This might affect the CT2/CT1 ratio. Always validate the cut-off values used for normals and for patients carrying the Factor V mutation, and, if necessary, adjust them.**

#### 4. Programming of the STA. and STA C analyzer:

Click on the icon **set up software** for the **manager program** and create the program according to:

##### A. Create the program for V-L2 (With Activated Prot C (APC))

CONFIGURATION – Page 1 / 3

Identification										
Abbreviation		Name			Method			Date		
VL2		HEMO VL2			CLOT-BASED					
Sample				Diluent						
Volume	Incu	Dil.	Id.	Name			Vial.	stab	Min Vol	
50µl	0 sec	1/10	11361	OWREN KOLLER			15 ml	72h	0.9 ml	
Reagents										
Id.		Name		Incu sec	Vol. µl	Vial ml	Stab h	Min Vol ml	Washing	
Ra R1		R1 HEM VL		0	50	4	8	0.2	Before	after
Rb R2B		R2B HEM VL		120	25	1	8	0.2	no	Special
Rc										
Rd R3		R3 HEM VL			50	4	8	0.2	no	Intensive
Analysis			Result				Validation (sec)			
Min time		40 sec	Primary Units		sec		Standard	min		max
Maxi time		250 sec	Corrector		1.00			10		250
Mean time		70 sec	Single/Duplicate		single					
Rd Heating		Yes	Precision		5%					
Stirring		Yes	Redil. Condition		sec					
Clot type		normal	<							
			>							

Enter data concerning washing.

Washing				
	Ra	Rb	Rc	Rd
Before	No	no		No
After	Normal	Special		Intensive
Wash		12227		
Name		DESORB U		
Stab. H		15		
Vial ml		120		
Vol in. ml		0.90		

Show the second page and write calibration data.

**CONFIGURATION – Page 2 / 3**

Control						
	Id.	Key	Name	Vial. ml	Stab. H	Min Vol ml
<b>Level 1</b>						
<b>Level 2</b>						

Go to the last page and enter data concerning the quality control.

Le STA can manage 3 levels of quality control. For example: the following configuration is possible.

**CONFIGURATION – Page 3 / 3**

Printout / Transmission							
Units	Convers. Factor	Print	Transmission Test number	Usual Values		Sec	
				STANDARD	Min.	Max.	
Main	Sec	yes	0		10	250	
Aux.1							
Aux.2							
Aux.3							
Printout limits		min	10				
		max	250				
Controls							
	Id.	Key	Name	Period H	Vial. ml	Stab. H	Min Vol ml
<b>Level 1</b>	BNC		B N CONT	24	1	8	0.2
<b>Level 2</b>	VL AB		VL AB CONT	24	1	8	0.2
<b>Level 3</b>							

**B. Create the program for V-L1 (without Activator Prot C)**

Identification										
Abbreviation		Name			Method			Date		
VL 1		HEMO VL1			CLOT-BASED					
Sample				Diluent						
Volume	Incu	Dil.	Id.	Name			Vial.	stab	Min Vol	
50µl	0 sec	1/10	11361	OWREN KOLLER			15 ml	72h	0.9 ml	
Reagents										
Id.		Name		Incu sec	Vol. µl	Vial ml	Stab h	Min Vol ml	Washing	
Ra R1		R1 HEM VL		0	50	4	8	0.2	no	Normal
Rb R2B		R2A HEM VL		120	25	1	8	0.2	no	Intensive
Rc										
Rd R3		R3 HEM VL			50	4	8	0.2	no	Intensive
Analysis			Result				Validation (sec)			
Min time		20 sec	Primary Units		sec		Standard	min	max	
Max time		90 sec	Corrector		1.00			20	90	
Mean time		40 sec	Duplicate/single		single					
Rd Heating		Yes	Precision		5%					
Stirring		Yes	Redil. Condition		sec					
Clot type		normal	<							
			>							

Enter data concerning washing.

Washing				
	Ra	Rb	Rc	Rd
Before	No	No		No
After	Normal	Intensive		Intensive
Wash				
Name				
Stab. H				
Vial ml				
Vol in. ml				

Show the second page and write calibration data.

**CONFIGURATION – Page 2 / 3**

Control						
	Id.	Key	Name	Vial. ml	Stab. H	Min Vol ml
<b>Level 1</b>						
<b>Level 2</b>						

Go to the last page and enter data concerning the quality control.

Le STA can manage 3 levels of quality control. For example: the following configuration it is possible.

**CONFIGURATION – Page 3 / 3**

Printout / Transmission							
Units	Convers. Factors	Print	Transmission Test number	Usual Values	Sec		
				<b>STANDARD</b>	<b>Min.</b>	<b>Max.</b>	
Main	Sec	yes	0		10	90	
Aux.1							
Aux.2							
Aux.3							
Printout limits		min	10				
		max	90				
Controls							
	Id.	Key	Name	Period H	Vial. ml	Stab. H	Min Vol ml
<b>Level 1</b>	BNC		B N CONT	24	1	8	0.2
<b>Level 2</b>	VL AB		VL AB CONT	24	1	8	0.2
<b>Level 3</b>							

Come back to the program “Test set up” and create a new test. Choose in the window the option “calculated Results”.  
 You get the following window,

Identification			
Id	Name	Method	Date
VLRATIO	VLRATIO	calculated Results	
Test		Printout /Transmission	
	Id	Name	Unit Ratio
1	VL2	HEMO VL2	Printout Yes
2	VL1	HEMO VL1	Range of transmission 0
Calculated			Printout limit (ratio)
Formula (Res. 1/Res. 2) x K			Minimum 0.10
Coefficient K 1.000			Maximum 5.00

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