Adaptation of BIOPHEN HEMOCLOT V-L on BCS

1. Reconstitution of BIOPHEN HEMOCLOT V-L (Ref ACK061K) reagent.
(Ref ACK061L)

Determination of Factor V-Leiden survival to the APC action using a clotting method

NAME	Reconstitution	Stability	Stabilization at the room Temperature
R1	4ml distilled water (*)	8 hours at room T° 24hours at 2-8 °C 1 month at -20°C	(**) 30 mn before any use
R2A	1ml distilled water (*)	8 hours at room T° 24hours at 2-8 °C 1 month at -20°C	(**) 30 mn before any use
R2B	1ml distilled water (*)	8 hours at room T° 24hours at 2-8 °C 1 month at -20°C	(**) 30 mn before any use
R3	4 ml distilled water (*)	8 hours at room T° 24hours at 2-8 °C 1 month at -20°C	(**) 30 mn before any use

Reconstitution:

(*) Following reconstitution with distilled water let the reagents to stabilize for 30 minutes at room temperature.

Conservation of reagents:

Take care of putting up the specific caps back on the bottles before storing them at 2°-8°C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Stabilization of reagents:

(**) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate before any use.

2. Determination of Factor V-Leiden activity

NAME	Reconstitution	Stability	Stabilization in T°
Quality controls Biophen Normal Control (réf A223201)	1 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T° Do not freeze	30 minutes on <i>BCS board</i> before any use (**)
Biophen ACT PCr Control. (réf A223405)	0.5 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T° Do not freeze	30 minutes on <i>BCS board</i> before any use (**)

Reconstitution:

(*) Following reconstitution of controls with distilled water, let them to stabilize for 30 minutes at room temperature.

Conservation of reagents:

(**) Following care of strictly respecting the 30 minutes temperature stabilization time for controls at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C. Homogenize before each use.

Foot-note: Do not freeze controls.

3. Results:

- The clotting times obtained for patients plasmas and controls with (T2) or without (CT1) APC, are expressed in seconds, for each patient, the ratio VL2/VL1 or T2/T1 between both clotting times must be calculated.
- The ratio of the Clotting Times obtained with or without APC, T2/T1, allows measuring the sensitivity of Factor V, in the tested specimen, to the action of Activated Protein C.
- Normal plasma, containing normal Factor V, is sensitive to this action, and yields a ratio: $T2/T1 \geq 2.00$
- Plasmas from patients carrying the R506Q mutation of Factor V, i.e. containing Factor V Leiden, yield a ratio: $T2/T1 \leq 1.80$.
- Molecular biology allows confirming the diagnosis, and classifying patients as heterozygous or homozygous.

NOTA: For a same reagent lot and a same tested specimen, T1 and T2 may present variations according to the instrument used, and to the clot detection sensitivity adjustment. This might affect the T2/T1 ratio. Always validate the cut-off values used for normals and for patients carrying the Factor V mutation, and, if necessary, adjust them.

5. Programming the BCS analyzer

Click on the icon **set up software** for the **manager program** and create the program according to:

A. Create Reagents

Name: R1 HEMOCLOT VL			
Abbreviation: R1 HVL		clean when the reagent is modified Normal	
<input type="checkbox"/> New calibration with any new lot of reagents			
Authorized lines			
from	1	to	14
Stability *			
< 15C	unlimited	> = 15C	
Mix by jet			
Intensity	1		
Frequency	No agitation		
Reference numbers			
	Number		
	add	remove	replace

Name: R2A HEMOCLOT VL			
Abbreviation: R2A HVL		clean when the reagent is modified Intensive	
<input type="checkbox"/> New calibration with any new lot of reagents			
Authorized lines			
from	1	to	14
Stability *			
< 15C	unlimited	> = 15C	
Mix by jet			
Intensity	1		
Frequency	No agitation		
Reference numbers			
	Number		
	add	remove	replace

Name: R2B HEMOCLOT VL			
Abbreviation: R2B HVL		clean when the reagent is modified Special	
<input type="checkbox"/> New calibration with any new lot of reagents			
Authorized lines			
from	1	to	14
Stability *			
< 15C	unlimited	> = 15C	
Mix by jet			
Intensity	1		
Frequency	No agitation		
Reference numbers			
	Number		
	add	remove	replace

Name: R3 HEMOCLOT VL			
Abbreviation: R3 HVL		clean when the reagent is modified Intensive	
<input type="checkbox"/> New calibration with any new lot of reagents			
Authorized lines			
from	1	to	14
Stability *			
< 15C	unlimited	> = 15C	
Mix by jet			
Intensity	1		
Frequency	No agitation		
Reference numbers			
	Number		
	add	remove	replace

B. Create in the same way control plasmas files.

Name: Biophen Normal Control			
Abbreviation: BNC		clean when the reagent is modified low	
<input checked="" type="checkbox"/> New calibration with any new lot of reagents			
Authorized lines			
from	1	to	14
Stability *			
< 15C		> = 15C	
Mix by jet			
Intensity	1		
Frequency	No agitation		
Reference numbers			
	Number		
	add	remove	replace

Name: Biophen plasma ACT PCr Control			
Abbreviation: V-L AB CONT		clean when the reagent is modified low	
<input checked="" type="checkbox"/> New calibration with any new lot of reagents			
Authorized lines			
from	1	to	14
Stability *			
< 15C		> = 15C	
Mix by jet			
Intensity	1		
Frequency	No agitation		
Reference numbers			
	Number		
	add	remove	replace

* To be filled by the user

C. Create the program for V-L1 (without Ativator Prot C)

General Information

Change Former			
Procedure Test Number			
Procedure Test Name	HEMOCLOT VL1		
Measurement time stopped after	90 sec	Or	mA
	0		
Speed	normal		
Stirring	normal		
Wavelength	405 nm		
Primary absorbance range	131-160 mA		
Duplicate n# for samples and controls	1	CV Authorized	0.0 %
Duplicate n# for calibration	1	CV Authorized	0.0 %
Duplicate n# of summary of raw data	Arithmetic mean		
Dilution Factor	1 : 10	(1 = no dilution)	
Dilution Buffer	OWREN KOLLER		

Evaluate and Control Method

Absorption				
Lag Phase	15	Sec	Threshold	200 mA
Angle				
Lag	15	Sec	Threshold	100 mA
Window size	3	Sec		
Angle minimum	15.0			
Down				
Lag Phase	15	Seconds		
Max decrease OD	10	%		
Sinusoid				
Lag Phase	15	Sec	Threshold base	100 mA
Minimal Amplitude	100	mA	Threshold factor	2.0
Number of tray	10		Threshold	20.0

Pipetting Sequence

Procedure Test Number

Procedure Test Name

HEMOCLOT VL1

Sampling cycles

Number of Cycle	Arm of transmission	Final Washing	Action of Rotors/transfer	Range Time [sec]	
				Min	Max
1	Arm	normal	No mix. in priority	0	0
2	Reagent Arm	normal	No mix/measure	0	0
3	Reagent Arm	normal	Incubation	110	130
4	Reagent Arm	normal	No mix/measure	0	0

Number of transm.		Mean/pipetting goal	Reagent	Speed	Volume (µl)
1	1	Take air		Slow	20
	2	Take the sample		Slow	50
	3	In the external cup		Slow	50
2	1	Take air		Slow	20
	2	Take the Reagent	R1 HVL	Slow	50
	3	In the central cup		Slow	50
3	1	Take air		Slow	20
	2	Take the Reagent	R2 A HVL	Slow	25
	3	In the central cup		Slow	25
4	1	Take air		Slow	20
	2	Take the Reagent	R3 HVL	Slow	50
	3	In the internal cup		Slow	50
5	1	Take the Reagent	SCS Cleaner	Slow	0

Test Definition

General Information				
Former Change		<input type="text"/>		
Name		HEMOCLOT VL1		Test N°
Result Unit				PC ID t
Behring	Sec	Factor Conversion		Abbreviation
User	Sec	1 Sec	1.0	Sec
PC	Sec	1 Sec	1.0	Sec
				Place in the list
				Digit number for results
				HEMOVL1
				100.0

Detailed test procedures

Procedure Name test	Evaluation Method
HEMOCLOT VL1	Fix Absorption

Calculate Formula of raw data

Formula Type	Minimum
User own's Formula	<input type="text"/>
Unit	Sec

Detailed test evaluation

		Determination of result		Standard Concentrations
Validation	<input type="text" value="20"/> Sec	<input type="text" value="90"/> Sec		
Test used for Calibration	<input type="text"/>		Sec	
Evaluation	Raw-mode = Result			
Reference Curve	Without calibration curve			
Minimum raw value	<input type="text"/>	sec		
Upper Extrapolation	<input type="text"/>	x maxi concentration		
Lower Extrapolation	<input type="text"/>	x mini concentration		
Calibrator				
Authorized SD respectively to the reference curve	<input type="text"/>	%		
Maximal number repeats	<input type="text"/>			

D. Create the program for CT2 (With activator prot C)

General Information

Former Change					
Procedure Test Number					
Procedure Test Name	HEMOCLOT VL2				
Measurement time stopped after	250 sec	Or	mA		
	0				
Speed	normal				
Stirring	normal				
Wavelength	405 nm				
Primary absorbance range	131-160 mA				
Duplicate n# for samples and controls	1	CV Authorized	<table border="1"><tr><td>0.0</td><td>%</td></tr></table>	0.0	%
0.0	%				
Duplicate n# for calibration	1	CV Authorized	<table border="1"><tr><td>0.0</td><td>%</td></tr></table>	0.0	%
0.0	%				
Duplicate n# of summary of raw data	Arithmetic mean				
Dilution Factor	1 : 10 (1 = no dilution)				
Dilution Buffer	OWREN KOLLER				

Evaluate and Control Method

Absorption					
Lag Phase	15	Sec	Threshold	200	mA
Angle					
Lag	15	Sec	Threshold	100	mA
Window size	3	Sec			
Angle minimum	15.0				
Down					
Lag Phase	15	Seconds			
Max decrease OD	10	%			
Sinusoid					
Lag Phase	15	Sec	base limit	100	mA
Minimal Amplitude	100	mA	Threshold factor	2.0	
Number of tray	10		Threshold	20.0	

Pipetting Sequence

Procedure Test Number

Procedure Test Name

HEMOCLOT VL2

Sampling cycles

Number of Cycle	Arm of transmission	Final Washing	Action of Rotors/transfer	Range Time [sec]	
				Min	Max
1	Arm	normal	No mix. in priority	0	0
2	Reagent Arm	normal	No mix/measure	0	0
3	Reagent Arm	normal	Incubation	110	130
4	Reagent Arm	normal	No mix/measure	0	0

Num of transm.		Mean/pipetting goal	Reagent	Speed	Volume (µl)
1	1	Take air		Slow	20
	2	Take the sample		Slow	50
	3	In the external cup		Slow	50
2	1	Take air		Slow	20
	2	Take the Reagent	R1 HVL	Slow	50
	3	In the central cup		Slow	50
3	1	Take air		Slow	20
	2	Take the Reagent	R2 B HVL	Slow	25
	3	In the central cup		Slow	25
4	1	Take air		Slow	20
	2	Take the Reagent	R3 HVL	Slow	50
	3	In the internal cup		Slow	50
5	1	Take the Reagent	SCS Cleaner	Slow	0

Test Definition

General Information				
Former Change <input style="width: 150px;" type="text"/>				
Name		<input style="width: 150px;" type="text" value="HEMOCLOT VL2"/>		°Test N° <input style="width: 50px;" type="text"/>
Result Unit		PC ID t		<input style="width: 50px;" type="text"/>
Behring	Sec	Conversion factor		Abbreviation <input style="width: 50px;" type="text" value="HEMOVL2"/>
User	Sec	1 Sec	<input style="width: 30px;" type="text" value="1.0"/>	Sec
PC	Sec	1 Sec	<input style="width: 30px;" type="text" value="1.0"/>	Sec
				Place in the list <input style="width: 50px;" type="text"/>
				Digit number for results <input style="width: 50px;" type="text" value="100.0"/>

Detailed Test procedures

Procedure test Name	Evaluation method
<input style="width: 150px;" type="text" value="HEMOCLOT VL2"/>	<input style="width: 150px;" type="text" value="Fix Absorption"/>
<input style="width: 150px;" type="text"/>	<input style="width: 150px;" type="text"/>
<input style="width: 150px;" type="text"/>	<input style="width: 150px;" type="text"/>
<input style="width: 150px;" type="text"/>	<input style="width: 150px;" type="text"/>

Calculate Formula of raw data

Type Formula	Minimum	<input style="width: 100px;" type="text"/>
User own's Formula	<input style="width: 250px;" type="text"/>	
Unit	Sec	<input style="width: 100px;" type="text"/>

Detailed test evaluation

		Determination of result		
Validation	<input style="width: 50px;" type="text" value="40.0"/> Sec	<input style="width: 50px;" type="text" value="250"/> Sec	Standard Concentrations Sec	
Test used for Calibration	<input style="width: 150px;" type="text"/>			
Evaluation	<input style="width: 150px;" type="text" value="Raw-mode= Result"/>			
Reference Curve	<input style="width: 150px;" type="text" value="Without curve calibration"/>			
Minimum raw value	<input style="width: 50px;" type="text"/>	sec		
Upper Extrapolation	<input style="width: 50px;" type="text"/>	x maxi concentration		
Lower Extrapolation	<input style="width: 50px;" type="text"/>	x mini concentration		
Calibrator				
Authorized SD respectively to the reference curve	<input style="width: 30px;" type="text"/>	%		
Maximal number repeats	<input style="width: 30px;" type="text"/>			

E. Metatest definition

General Information				
Name	HEMOCLOT VL RATIO		Test N°	
Result Unit			ID PC	
Behring	Sec	Factor Conversion	Abbreviation	Ratio HEMO VL
User		1 = 1.0	List sequence	
PC		1 = 1.0	Digit format	10.0

Basic details tests

Test name
HEMO VL2
HEMO VL1

Individual Tests Combination

Formula Type	User formula	
User formula	HVL2/HVL1	
Control value		%

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