Adaptation of BIOPHEN HEMOCLOT V-L on STA-R

1. Reconstitution of BIOPHEN HEMOCLOT V-L (Ref ACK061K) reagent.
(Ref ACK061L)

Determination of factor V-Leiden activity with a clotting method

NAME	Reconstitution	Stability	Stabilization at the room Temperature
R1	4ml distilled water (*)	8 hours at room T° 24hours at 2-8 °C 1 month at -20°C	30 mn before any use (**)
R2A	1ml distilled water (*)	8 hours at room T° 24hours at 2-8 °C 1 month at -20°C	30 mn before any use (**)
R2B	1ml distilled water (*)	8 hours at room T° 24hours at 2-8 °C 1 month at -20°C	30 mn before any use (**)
R3	4 ml Distilled water (*)	8 hours at room T° 24hours at 2-8 °C 1 month at -20°C	30 mn before any use (**)

Reconstitution:

(*) Following reconstitution with distilled water, let the reagents to stabilize for 30 minutes at room temperature.

Conservation of reagents:

Take care of putting up the specific caps back on the bottles before storing them at 2°-8°C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Stabilization of reagents:

(**) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate before any use.

2. Determination of Factor V-Leiden activity

NAME	Reconstitution	Stability	Stabilization in T°
Quality controls Biophen Normal Control (ref A223201)	1 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T° Do not freeze	30 minutes <i>on STA-R board</i> before any use (**)
Biophen ACT PCr Control. (ref A223405)	0.5 ml of distilled water (*)	24 hours at 2-8°C 8 hours at room T° Do not freeze	30 minutes <i>on STA-R board</i> before any use (**)

Reconstitution:

(*) Following reconstitution controls with distilled water, let them to stabilize for 30 minutes at room temperature.

Conservation of reagents:

(**) Following care of strictly respecting the 30 minutes temperature stabilization time for controls at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C. Homogenize before each use.

Foot-note: Do not freeze controls.

3. Results:

- The clotting times obtained for patients plasmas and controls with (CT2) or without (CT1) APC, are expressed in seconds, for each patient, the ratio CT2/CT1 between both clotting times must be calculated.
- The ratio of the Clotting Times obtained with or without APC, CT2/CT1, allows measuring the sensitivity of Factor V, in the tested specimen, to the action of Activated Protein C.
- Normal plasma, containing normal Factor V, is sensitive to this action, and yields a ratio: $CT2/CT1 \geq 2.00$
- Plasmas from patients carrying the R506Q mutation of Factor V, i.e. containing Factor V Leiden, yield a ratio: $CT2/CT1 \leq 1.80$.
- Molecular biology allows confirming the diagnosis, and classifying patients as heterozygous or homozygous.

NOTA: For a same reagent lot and a same tested specimen, CT1 and CT2 may present variations according to the instrument used, and to the clot detection sensitivity adjustment. This might affect the CT2/CT1 ratio. Always validate the cut-off values used for normals and for patients carrying the Factor V mutation, and, if necessary, adjust them.

4. Programming of the STA.-R analyzer:

Click on the icon **set up software** for the **manager program** and create the program for CT2 (VL2) according to:

TESTS				Configuration				
Method								
IDENTIFICATION								
Abbreviation			CT2		Last up date			
Name		HEMO VL2			Method	CLOT-BASED		
SAMPLE				DILUENT				
Vol. μ l	Incu.	Dilution	Id.	Name		Stab. h	Min Vol	
50 μ l	0 sec	1/10	11361	OWREN KOLLER		72	0.9	
REAGENTS								
	Id.	Name		Stab. h	Vol. μ l	Incu. sec	Prec	Continued
Ra	R1	R1 HEM VL		8	50	0	<input type="checkbox"/>	
Rb	R2B	R2B HEM VL		8	25	120	<input type="checkbox"/>	
Rc							<input type="checkbox"/>	
Rd	R3	R3 HEM VL		8	50			

Enter data concerning washing

	Id.	Name		Stab. h	Vol. μ l	Incu. sec	Prec.	Vial. ml	Min Vol ml
Diluent	11361	OWREN KOLLER		72				15	0.9
Ra	R1	R1 HEM VL		8	50	0	<input type="checkbox"/>	4	0.2
Rb	R2B	R2B HEM VL		8	25	120	<input type="checkbox"/>	1	
Rc							<input type="checkbox"/>		
Rd	R3	R3 HEM VL		8	50			4	0.2
Washing									
	Ra		Rb		Rc		Rd		
Before	No		No				No		
After	Normal		Special				Intensive		
Wash			12227						
Name			DESORB U						
Stab. H			120						
Vial ml			15						
Vol in. ml			0.90						

Show the last page and enter data concerning the quality control.

TESTS						Configuration
						Q.C
LEVEL 1						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control	BNC		BN CONT	8	1	0.2
LEVEL 2						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control	VLAB		VLAB CONT	8	1	0.2
LEVEL 3						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control						
Period: <input type="text"/> hours <input type="text"/> <input type="checkbox"/> Tests <input type="checkbox"/> Vial						

2. Configuration CL1 (VL1)

TESTS							Configuration	
Method								
IDENTIFICATION								
Abbreviation		CT1			Last up date			
Name		HEMO VL1			Method		CLOT-BASED	
SAMPLE				DILUENT				
Vol. μ l	Incu.	Dilution	Id.	Name			Stab. h	Min Vol
50 μ l	0 sec	1/10	11361	OWREN KOLLER			72	0.9
REAGENTS								
	Id.	Name		Stab. h	Vol. μ l	Incu. sec	Prec	Continued
Ra	R1	R1 HEM VL		8	50	0	<input type="checkbox"/>	
Rb	R2A	R2A HEM VL		8	25	120	<input type="checkbox"/>	
Rc							<input type="checkbox"/>	
Rd	R3	R3 HEM VL		8	50			

Enter data concerning washing

	Id.	Name	Stab h	Vol. μ l	Incu. sec	Prec.	Vial. ml	Min Vol ml
Diluent	11361	OWREN KOLLER	72				15	0.9
Ra	R1	R1 HEM VL	8	50	0	<input type="checkbox"/>	4	0.2
Rb	R2A	R2A HEM VL	8	25	120	<input type="checkbox"/>	1	
Rc						<input type="checkbox"/>		
Rd	R3	R3 HEM VL	8	50			4	0.2
Washing								
	Ra		Rb		Rc		Rd	
Before	No		No				No	
After	Normal		Intensive				Intensive	
Wash								
Name								
Stab. H								
Vial ml								
Vol in. ml								

Click on the icon **Result software** for the **manager program** and create the program according to:

TESTS		Configuration					
		IV. Result					
METHOD CLOT-BASED		RESULT					
Min. Time	20	Primary units.	Sec				
Max Time	90	Corrector	1.00				
Mean Time	40	Single/Duplicate					
Clot type	<input checked="" type="radio"/> Normal			<input checked="" type="radio"/> Single			
	<input type="radio"/> weak			<input type="radio"/> Duplicate			
	<input type="radio"/> Mean						
		Precision (%)	5.00	VALIDATION			
<input checked="" type="checkbox"/> Rd Heating <input checked="" type="checkbox"/> Stirring		Redilution	Condition			Min.	Max.
		Sec	<				10
		>					

Write calibration data.

TESTS		Configuration	
		V. Calibration	
MODE :	Brut : The results are given in measured units (sec)	Determination	
		<input type="checkbox"/> Duplicate <input checked="" type="checkbox"/> Single	

Click on the icon **Printout/Transmission software** for the **manager program** and create the program according to:

TESTS		Configuration			
VI.		Printout/Trans.			
PARAMETERS					
	Unit	Conver Factor	Print	Transmission Test number	Usual Values Min. Max.
Main	SEC		<input checked="" type="checkbox"/>	0	10 90
Aux1			<input type="checkbox"/>		
Aux2			<input type="checkbox"/>		
Aux3			<input type="checkbox"/>		
Printout Limits					
Min.	10	Max.	90		

Show the last page and enter data concerning the quality control.

TESTS							Configuration
							Q.C
LEVEL 1							
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml	
Control	BNC		BN CONT	8	1	0.2	
LEVEL 2							
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml	
Control	VLAB		VLAB CONT	8	1	0.2	
LEVEL 3							
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml	
Control							
Period : <input type="text"/> hours <input type="text"/> <input type="checkbox"/> Tests <input type="checkbox"/> Vial							

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