Manufactured By: HYPHEN BioMed

HEMOCLOT FV reagent ON STA-R

Adaptation of HEMOCLOT FV ON STA-R

1. Reconstitution of HEMOCLOT FV(Ref ACK071K)reagent.

Determination of factor V activity with a clotting method

NAME	Reconstitution	Stability	Stabilization at room Temperature
Hemoclot FV Reagent	1 ml Distilled water(*)	72 hours 2-8°C 24 hours at room T° (18-25°C)	30 mn before any use (**)

Reconstitution:

(*) After reconstitution with distilled water, let the factor V reagent to stabilize 15 minutes at room temperature.

Storage of reagents:

Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Stabilization of reagents:

(**) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate board before any use.

2. Preparation of the calibrators and controls

Calibration curve:

- Prepare 1 ml of normal human pooled plasma diluted 1:10 in Owren Koller buffer (by definition the normal citrated plasma pool corresponds to a concentration of 100% of Factors V).

Or

- The calibration curve can also be established with the BIOPHEN Plasma Calibrator (# 222101), using the factor V activity indicated for the lot used.

Name	Reconstitution	Stability	Stabilization in T°
CALIBRATION Biophen Plasma Calibrator (ref A222101)	1 ml of distilled water(*)	8 hours at 2° - 8°C 4 hours at Room T°	30 minutes <i>on STA board</i> before any use (**)

Controls

Nom	Reconstitution	Stability	Stabilization in T°
CONTROLES DE QUALITE Biophen Normal Control (ref A223201) Biophen Abnormal Control (ref A223301)	1 ml of distilled water(*)	8 hours at 2° - 8°C 4 hours at Room T°	30 minutes <i>on STA board</i> before any use (**)

Reconstitution:

(*) After reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature.

Stabilization of reagents:

(**) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate before any use.

Storage of reagents:

Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Do not freeze calibrators and controls

Note : A calibration curve must be carried out for each new batch of Hemoclot Factor V reagent, after each important maintenance of the instrument, and if QC values are out of the acceptance range for the method.

NOTA: Plasmas used for the reagent preparation were tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions as being potentially infectious.

3. Results:

The calibration curve is of the Lin (Clot time in seconds) - Log (% activity V) type.

The values obtained for the patients and controls are directly calculated from the calibration curve, when tested at the standard 1:10 dilution in the test.

The results are expressed in % FV activity.

4. Programming of the STA R analyser

Click on the icon **set up software** for the **manager program** and create the program according to:

TESTS				Configuration			
Method							
IDENTIFICATION							
Abréviation		HEMO FV		Last up date.			
Name	HEMOCLLOT FV reagent			Method	CLOT-BASED		
SAMPLE			DILUENT				
Vol. µl	Incu.	Dilution	Id.	Name	Stab. h	Min. Vol.	
50 µl	0 sec	1/10	11361	OWREN KOLLER	72	0.90	
REAGENTS							
	Id.	Name	Stab. h	Vol. µl	Incu. sec	Prec.	Continued
Ra	R1	HEMO FV	24	50	60	<input type="checkbox"/>	
Rb						<input type="checkbox"/>	
Rc						<input type="checkbox"/>	
THR	THR	THROMBO	24	100			

Enter data concerning washing

	Id.	Nom	Stab. h	Vol. µl	Incu. sec	Prec.	Vial. ml	Min Vol. ml
Diluant	11361	OWREN-KOLLER	72				15	0.9
Ra	R1	HEMO FV	24	50	60	<input type="checkbox"/>	1	0.2
Rb						<input type="checkbox"/>		
Rc						<input type="checkbox"/>		
Rd	THR	THROMBO	24	100			5	0.5
Washing								
	Ra		Rb		Rc		Rd	
Before	No						No	
After	intensive						Intensive	
Wash								
Name								
Stab. H								
Vial ml								

Click on the icon **Result software** for the **manager program** and create the program according to:

TESTS		Configuration					
		I. Result					
METHOD: CLOT-BASED		RESULT		VALIDATION (%)			
Min. Time	10	Primary Units	%				
Max Time	300	Corrector	1.00				
Mean Time	35	Determination					
Clot	<input checked="" type="radio"/> Normal	<input checked="" type="radio"/> Single					
	<input type="radio"/> weak					<input type="radio"/> Duplicate	
	<input type="radio"/> Mean						
<input checked="" type="checkbox"/> Rd Heating		Precision (%)	5.00			Min.	Max.
<input checked="" type="checkbox"/> Stirring		Redilution	Condition	0	130		
		1/	<				
		1/	>				

Write calibration data.

TESTS		Configuration			
		II. Calibration			
MODE :		Linear	Détermination		
		<input checked="" type="checkbox"/> Duplicate <input type="checkbox"/> Single			
SCALE		Concentrations			
RAW DATA : <input checked="" type="checkbox"/> Linear <input type="checkbox"/> Log		Linear <input type="checkbox"/> Inverse <input checked="" type="checkbox"/> Log			
	Id.	Key	Name	Stab. H	Dilution
Calibrator 1	BPC		BPCALV	8*	1/10
Calibrator 2	BPC		BPCALV	8*	1/20
Calibrator 3	BPC		BPCALV	8*	1/40
Calibrator 4	BPC		BPCALV	8*	1/80
Calibrator 5	BPC		BPCALV	8*	1/160
Ctrl, Niv. 1					
Ctrl, Niv. 2					

*according to the calibrator used.

Click on the icon **Printout/Transmission software** for the **manager program** and create the program according to:

TESTS		Configuration				
III.		Printout/Trans.				
PARAMETERS						
	Unit	Factors. Conver	Print	Transmission Test number	Usual Values Min. Max.	
Main	%		<input checked="" type="checkbox"/>	0	70	120
Aux1	SEC		<input checked="" type="checkbox"/>	0		
Aux2			<input type="checkbox"/>			
Aux3			<input type="checkbox"/>			
Printout Limits						
Min.	<input type="text" value="1"/>	Max.	<input type="text" value="150"/>			

Show the last page and enter data concerning the quality control.

TESTS		Configuration				
		Q.C				
LEVEL 1						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control	BNC		BN CONT*	8*	1*	0.2*
LEVEL 2						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control	BAC		BA CONT**	8*	1*	0.2*
LEVEL 3						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control						
Period : <input type="text"/> hours <input type="text"/> <input type="checkbox"/> Tests <input type="checkbox"/> Vial						

(*) to be filled by the user according to the controls used

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