Manufactured By: HYPHEN BioMed

FIBRIPHEN

on STA-R

Quantitative determination of Fibrinogen by clotting assay

1. Reagents preparation:

NAME	Reconstitution	Stability	Stabilization at the room Temperature
Fibriphen 1 (ACK571K)	1 ml distilled water*	7 days at room T° (18-25°C)* 14 days at 2-8 °C*	30 mn before any use (**)
Fibriphen 2 (ACK572K)	2 ml distilled water*	7 days at room T° (18-25°C)* 14 days at 2-8 °C*	30 mn before any use (**)
Fibriphen 5 (ACK575K)	5 ml distilled water	7 days at room T° (18-25°C)* 14 days at 2-8 °C*	30 mn before any use (**)

(*) Provided any contamination or evaporation is avoided.

Reconstitution: (*) After reconstitution with distilled water, let the reagent to stabilize 15 minutes at room temperature.

Stabilization of reagents: (**) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate before any use.

Storage of reagents: Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Note : Homogenize the reagent before each use.

NOTE: Any reagent of biological origin must be handled with all the required cautions as being potentially infectious.

2. Preparation of the calibrators and controls

Calibration curve:

- Prepare 1 ml of normal citrated human plasma pool (the 1:20 dilution in Owren Koller buffer, directly managed by the automate, corresponds by definition to a concentration of 3 g/L Fibrinogen).

Or

- The calibration curve can also be established with the BIOPHEN Plasma Calibrator (# A222101), using the Fibrinogen concentration indicated for the lot used.

Name	Reconstitution	Stability	Stabilization in T°
CALIBRATION Biophen Plasma Calibrator (ref A222101)	1 ml of distilled water(*)	See insert	30 minutes <i>on</i> STAR <i>board</i> before any use (**)

Controls

Nom	Reconstitution	Stability	Stabilization in T°
CONTROLES DE QUALITE Biophen Normal Control (ref A223201) Fibrinogen Control Low (ref ASC070K)	1 ml of distilled water(*)	See insert	30 minutes <i>on</i> STAR <i>board</i> before any use (**)

Reconstitution: (*) After reconstitution with distilled water, let the reagent to stabilize 30 minutes at room temperature.

Stabilization of reagents: (**) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate before any use.

Storage of reagents: Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Do not freeze calibrators and controls

Note : A calibration curve must be carried out for each new batch of reagent, after each important maintenance of the instrument, and if QC values are out of the acceptance range for the method.

NOTE: Plasmas used for the reagent and calibrators/controls preparation were tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of biological origin must then be handled with all the required cautions as being potentially infectious.

3. Results:

The calibration curve is of the Log (Clot time in seconds) - Log (g/L Fibrinogen) type.

The values obtained for the patients and controls are directly calculated from the calibration curve, when tested at the standard 1:20 dilution in the test.

The results are expressed in g/L Fibrinogen.

4. Programming of the STA R analyser

Click on the icon **set up software** for the **manager program** and create the program according to:

TESTS							Configuration	
Method								
IDENTIFICATION								
Abréviation			FIBRIPHEN		Last up date.			
Name		FIBRIPHEN			Method		CLOT-BASED	
SAMPLE			DILUENT					
Vol. µl	Incu.	Dilution	Id.	Name	Stab. h	Min. Vol.		
100 µl	240 sec	1/20	OK	OWREN KOLLER	*	*		
REAGENTS								
	Id.	Name	Stab. h	Vol. µl	Incu. sec	Prec.	Continued	
Ra						<input type="checkbox"/>		
Rb						<input type="checkbox"/>		
Rc						<input type="checkbox"/>		
THR	R1	Fibriphen	168	50				

Enter data concerning washing

	Id.	Nom	Stab. h	Vol. µl	Incu. sec	Prec.	Vial. ml	Min Vol. ml
Diluant	OK	OWREN-KOLLER	72				*	*
Ra						<input type="checkbox"/>		
Rb						<input type="checkbox"/>		
Rc						<input type="checkbox"/>		
Rd	R1	Fibriphen	168	50			1 or 2 or 5*	0.2 or 0.5*
Washing								
	Ra		Rb		Rc		Rd	
Before							No	
After							Special	
Wash							12227	
Name							STA-DESORB U	
Stab. H							120	
Vial ml							150	
Vol min ml							0.9	

*according to the reagent used

Click on the icon **Result software** for the **manager program** and create the program according to:

TESTS		Configuration	
		I. Result	
METHOD: CLOT-BASED		RESULT	
Min. Time	3	Primary Units	g/L
Max Time	100	Corrector	1.00
Mean Time	22	Determination	
Clot	<input type="radio"/> Normal	<input checked="" type="radio"/> Single <input type="radio"/> Duplicate	
	<input checked="" type="radio"/> weak		
	<input type="radio"/> Mean		
Rd Heating		Precision (%)	10%
Stirring		Redilution	Condition
		1/	<
		1/40	>
		12	
		VALIDATION (%)	
		Min.	Max.
		0.6	12

Write calibration data.

TESTS		Configuration	
		II. Calibration	
MODE : <input type="text" value="Linear"/>		Détermination	
		<input checked="" type="checkbox"/> Duplicate <input type="checkbox"/> Single	
SCALE		Concentrations	
RAW DATA : X Log Linear		Linear <input type="checkbox"/> Inverse X Log	
	Id.	Key	Name
Calibrator 1	BPC		BPC
Calibrator 2	BPC		BPC
Calibrator 3	BPC		BPC
Calibrator 4	BPC		BPC
			Stab. H
			Dilution
			*
			1/5
			*
			1/10
			*
			1/20
			*
			1/40
Ctrl, Niv. 1			
Ctrl, Niv. 2			

*according to the calibrator used.

Click on the icon **Printout/Transmission software** for the **manager program** and create the program according to:

TESTS		Configuration				
III.		Printout/Trans.				
PARAMETERS						
	Unit	Factors. Conver	Print	Transmission Test number	Usual Values Min. Max.	
Main	g/L		<input checked="" type="checkbox"/>	0	2	4
Aux1	SEC		<input checked="" type="checkbox"/>	0		
Aux2			<input type="checkbox"/>			
Aux3			<input type="checkbox"/>			
Printout Limits						
Min.	<input type="text" value="0.1"/>	Max.	<input type="text" value="15"/>			

Show the last page and enter data concerning the quality control.

TESTS		Configuration				
		Q.C				
LEVEL1						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control	BNC		BN CONT*	8*	1*	0.2*
LEVEL 2						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control	ACF		A CONT Fng*	8*	1*	0.2*
LEVEL3						
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml
Control						
Period : <input type="text"/> hours <input type="text"/> <input type="checkbox"/> Tests <input type="checkbox"/> Vial						

(*) to be filled by the user according to the controls used

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