Assay of FVIII:C or FIX or FXI or FXII by clotting assay on STA-R

APTT based assay, using a deficient plasma for the tested factor as substrate, in which all the other factors are present and in excess. The measured factor, present in the tested plasma, is then the limiting factor.

1. Reconstitution of the reagents:

NAME	Reconstitution	Stability **	Stabilization at the room Temperature
Deficient Plasma : FVIII:C: ADP040A/K Or FIX : ADP050A/K Or FXI: ADP070A/K Or FXII : ADP080A/K	1ml distilled water (*)	8 hours at room T° (18-25°C) 24hours at 2-8 °C	30 mn before any use (***)

Reagents required but not provided:

- Distilled Water;
- Calcium chloride 0.025M (ex AAR001A/AAR001K)
- Owren Koller buffer (ex: # AAR003A/AAR003K).
- Normal pool human plasma or Factor VIII:C or IX (or XI or XII) calibrator (ex: BIOPHEN Plasma Calibrator - # A222101).
- Normal and Abnormal quality control plasmas, titrated for Factor VIII:C or IX (or XI or XII) (ex: BIOPHEN Normal Control Plasma-#A223201 and BIOPHEN Abnormal Control Plasma - #A223301).
- Activated cephalin reagent (ex : CK Prest or Cephén # ACK511K/ACK512K/ACK515K-L).

Reconstitution: (*) Following reconstitution with distilled water, let the reagents to stabilize for 15 minutes at room temperature.

(**) Stability indicated are validated provided any evaporation or contamination of reconstituted reagents are avoided. They must be adjusted according to the actual laboratory working conditions.

Storage of reagents: Take care of putting up the specific caps back on the bottles before storing them at 2°-8° C and of strictly respecting the temperature stabilization time of 30 minutes before using the reagents on the automate.

Stabilization of reagents: (***) It is necessary to let the reagent temperature to stabilize for at least 30 minutes on the automate before any use.

Note : Homogenize the reagents before each use.

NOTE: Any reagent of biological origin must be handled with all the required cautions as being potentially infectious.

2. Calibration and controls

Ex:

NAME	Reconstitution	Stability	Stabilization in T°
Calibration: BIOPHEN Plasma Calibrator (ref A222101)	1 ml of distilled water (*)	*according to the factor, refer insert	30 minutes <i>on STA board</i> before any use (**)
Quality control: Biophen Normal Control (ref 223201)	1 ml of distilled water (*)	*according to the factor, refer insert	30 minutes <i>on STA board</i> before any use (**)
Quality control: Biophen Abnormal Control (ref 223301)	1 ml of distilled water (*)	*according to the factor, refer insert	30 minutes <i>on STA board</i> before any use (**)

Reconstitution: (*) Following reconstitution of calibrators and controls, let them to stabilize for 30 minutes at room temperature.

Storage of reagents:(**) Following care of strictly respecting the 30 minutes temperature stabilization time for controls at room temperature, then the 30 minutes on the automate, particularly if they were stored at + 2°-8°C.

Homogeneize before each use.

Foot-note: Do not freeze calibrators and controls.

Foot-note: Quality controls at different levels must be run regularly in each series, a new calibration curve must be carried out for each new batch of reagents, after an important maintenance of the instrument, or if measured values for controls are not in compliance with the one expected.

NOTE : It is recommended to perform a fresh calibration curve, patient's plasmas and controls successively without interruption for optimal performances.

Samples are loaded undiluted on the instrument; the dilutions are managed by the STAR.

3. Results:

- The calibration curve is of the Log (CT, in sec) - Log (concentration, in %) type.
- The values obtained for the patients and controls, tested at the standard 1:10 dilution, are directly calculated from the calibration curve, and are expressed in %.

NOTE: Clotting times and performances may present slight variations according to the activated cephalin reagent, the reagent lot and the instrument used, and to the clot detection sensitivity adjustment. Validate the expected values in the laboratory working conditions. Performances, as well as target values and acceptance ranges for each new lot of reagent or quality controls used, must then be confirmed (and adjusted if necessary) in the laboratory working conditions.

4. Programming the STA-R analyzer:

Click on the icon **set up software** for the **manager program** and create the program according to:

TESTS						Configuration		
Method								
IDENTIFICATION								
Abbreviation		aPTT DEF		Last up date				
Name		aPTT deficient		Method		CLOT-BASED		
SAMPLE			DILUENT					
Vol. μ l	Incu.	Dilution	Id.	Name			Stab. h	Min Vol
50 μ l	0 sec	1/10	OK	OWREN KOLLER			72	1.5
REAGENTS								
	Id.	Name		Stab. h	Vol. μ l	Incu. sec	Prec	Continued
Ra	DEF	Deficient plasma		8	50	0	<input type="checkbox"/>	
Rb	CK	Activated cephalin reagent		*	50	240	<input type="checkbox"/>	
Rc							<input type="checkbox"/>	
Rd	CaCl2	Calcium		*	50			

Enter data concerning washing

	Id.	Name		Stab. h	Vol. μ l	Incu. sec	Prec.	Vial. ml	Min Vol ml
Diluent	OK	OWREN KOLLER		72				25	1.5
Ra	DEF	Deficient plasma		8	50	0	<input type="checkbox"/>	1	0.2
Rb	CK	Activated cephalin reagent		*	50	240		*	*
Rc							<input type="checkbox"/>		
Rd	CaCl2	Calcium		*	50			*	*
Washing									
	Ra		Rb		Rc		Rd		
Before	No		No				No		
After	Special		Normal				Normal		
Wash	12227								
Name	DESORB U								
Stab. H	120								
Vial ml	15								
Vol in. ml	0.90								

*according to the reagent used

Click on window **Printout/Transmission** and enter following data:

TESTS						Configuration		
					Printout/Trans.			
PARAMETERS								
	Units	Convers Factor	Print	Transmission Test Number	Usual Values			
					Min.	Max.		
Main	%		<input checked="" type="checkbox"/>	0				
Aux1	sec		<input checked="" type="checkbox"/>	0				
Aux2			<input type="checkbox"/>					
Aux3			<input type="checkbox"/>					
Printout limits								
Min.	<input type="text" value="0"/>	Max.	<input type="text" value="150"/>					

Show the last page and enter data concerning the quality control.

TESTS						Configuration	
					Q.C		
LEVEL 1							
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml	
Control	BNC		BN CONT	*	1*	0.2*	
LEVEL 2							
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml	
Control	VLAB		VLAB CONT	*	1*	0.2*	
LEVEL 3							
	Id.	Key	Name	Stab. h	Vial. ml	Min Vol ml	
Control							
Period: <input type="text"/> hours <input type="text"/> <input type="checkbox"/> Tests <input type="checkbox"/> Vial							

* **Must be stirred and stabilized before each use**

* To be filled according to the controls used.

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