# BIOPHEN DTI

## Ref. A220202

Quantitative determination of Direct Thrombin Inhibitors (DTI) with a chromogenic method.

In vitro research use only

Manufactured By: HYPHEN BioMed

Last revision: 29/08/2008

### INTENDED USE:

BIOPHEN DTI kit is a quantitative assay for measuring the anti-IIa activity of hirudin and other Direct Thrombin Inhibitors (DTIs) in human citrated plasma, using a chromogenic method, manual or automated.

### CLINICAL APPLICATIONS:

Preventive or curative indications of Hirudin and other DTIs: quantitative determination of the anti-IIa activity in patients' plasma is an aid for monitoring the therapy and adjusting drug dosage.

### ASSAY PRINCIPLE:

Biophen DTI is a kinetics method based on the inhibition of a constant and in excess amount of Thrombin (FIIa), by the tested hirudin or DTI, and hydrolysis of a Thrombin specific chromogenic substrate, by the Thrombin in excess. pNA is then released from the substrate. The amount of pNA released is then a relation of the residual Thrombin activity. There is an inverse linear relationship between the concentration of hirudin (or DTI) and color development, measured at 405 nm.

[DTI] + [FIIa (excess)] → [FIIa-DTI] + [residual FIIa]  
[FIIa (residual)] + Substrate → Peptide + pNA

### SPECIMEN:

Plasma prepared from citrated anticoagulated blood, where hirudin (or any other direct thrombin inhibitor) activity must be measured.

### REAGENTS:

#### R1: Reagent 1: Thrombin Substrate

Thrombin specific chromogenic substrate (CS-01(81)), lyophilised in presence of stabilizers, of an heparin neutralizing substance, and of a fibrin polymerization inhibitor.  
2 vials (to be reconstituted with 2.5 mL of distilled water).

#### R2: Reagent 2: Thrombin

Purified human thrombin, lyophilised, in presence of stabilizers.  
2 vials (to be reconstituted with 2.5 mL of distilled water).

#### R3: Reagent 3: Tris-BSA Buffer

Tris NaCl reaction Buffer, pH 8.40, containing 1% BSA, and sodium azide as preservative.

Ready to use.  
2 vials of about 25 mL.

#### Warning:

-Purified thrombin was prepared by activation of prothrombin extracted from human plasma, which was tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. Bovine Serum Albumin (BSA) was prepared from bovine plasma, which was tested for the absence of infectious agents, and collected from animals free from BSE. However, no assay may warrant the total absence of infectious agents. Any product of biological origin must then be handled with all the required cautions, as being potentially infectious.  
- Sodium azide (0.9 g/l) may react with lead and copper plumbing to form highly explosive metal azides. Flush with large volumes of water when discarding into a sink.  
- This assay was optimized for minimizing any possible interference for heparinized samples.  
- The thrombin concentration is adjusted for each lot for providing the right reactivity in the assay

### REAGENTS AND MATERIAL REQUIRED BUT NOT PROVIDED:

#### Reagents:

- Distilled water, preferentially sterile.
- Acetic Acid (20%) or Citric Acid (2%) (End point method).
- Plasma Calibrator titrated for hirudin (eg: **Plasma Hirudin Standard Low (# ASC020K)** or **High (# ASC020L)**) Or Reference material for the DTI to be assayed (international or internal, pharmaceutical preparation...).
- Quality Control Plasmas titrated for hirudin (eg: **Plasma Hirudin Control (# ASC025K)** or for the DTI to be assayed.

#### Material:

- Spectrophotometer, photometer or automates for chromogenic assays, with a wave-length set up at 405 nm.
- Stop watch.
- Calibrated pipettes.

### TRACEABILITY TO THE REFERENCE MATERIAL:

Internal reference for Hirudin, or standard pharmacopoeia preparations (when available) or pharmaceutical preparations, diluted into plasma.

### STORAGE CONDITIONS:

Reagents must be stored at 2-8°C, in their original packaging box. They are then stable until the expiration date printed on the box.

**Note:** Stability studies for 3 weeks at 30°C show that the reagents can be shipped at room temperature for a short period without damage.

### PREPARATION AND STABILITY OF REAGENTS:

**Note:** Reconstitution volumes can vary according to the automate used. Refer to each specific instrument adaptation.

#### R1: Reagent 1: Thrombin specific chromogenic substrate

- Reconstitute each vial with 2.5 mL of distilled water. Shake thoroughly until complete dissolution of the content (vortex).
- Let to homogenize for 30 minutes at room temperature (18-25 °C), while shaking the vial from time to time.
- Homogenize the content before each use.

Stability of restored reagent R1, kept in its original vial:

- 72 hours at 2-8°C.
- 24 hours at room temperature (18-25 °C).
- 2 months at -20°C or below.

#### R2: Reagent 2: Thrombin

- Reconstitute each vial with 2.5 mL of distilled water. Shake thoroughly until complete dissolution of the content (vortex).
- Let to homogenize for 30 minutes at room temperature (18-25 °C), while shaking the vial from time to time.
- Homogenize the content before each use.

Stability of restored reagent R2, kept in its original vial:

- 72 hours at 2-8°C.
- 24 hours at room temperature (18-25 °C).
- 2 months at -20°C or below.

#### R3: Reagent 3: Tris-BSA Buffer

Ready to use buffer. Shake before use.

Stability of the buffer, protected from any bacterial contamination:

- In its original vial, until the expiration date printed on the label, at 2-8°C.
- When open, 7 days at 2-8 °C

#### Cautions:

- In order to improve stability, reagents must be closed with their original screw cap following each use (yellow cap for R1, white cap for R2, and white cap for buffer R3).
- Reagents must be handled with care, in order to avoid any contamination during use.
- If the substrate becomes yellow, this indicates the presence of a contaminant. It must be rejected, and a new vial must be used.

#### Note:

- R1 and R2 vials are closed under vacuum. Remove carefully the stopper, in order to avoid any lost of powder when opening the vials.
- According to the automated method used, the reagents can be reconstituted with volumes different from those recommended. In any case, the established reactive ratios (respective reagent concentrations in the reactive milieu) between R1 and R2 must be adhered to.
- Use only reagents from kits with the same lot number. Do not mix reagents from kits with different lots when running the assay. Reagents R1 and R2 are optimized for each lot of kits.

### TESTED SPECIMEN:

Blood (9 vol.) must be collected on 0.109M (or 0.129M) trisodium citrate anticoagulant (1 vol.); plasma supernatant is decanted following a 20 min. centrifugation at 2,500 g; citrated plasma must be tested within 4 hours when stored at room temperature (18-25°C), or it can be frozen at -20°C or below for up to 1 month. Just before use, the plasma must be thawed for 15 min. in a water bath at 37°C.

**Note:** Refer to GEHT or CLSI recommendations for further instructions on specimen collection, handling and storage. Discard any plasma presenting an unusual aspect (haemolysed, lipaemic aspect...).

### TEST PROCEDURE:

Biophen DTI kit is designed for being used with automated kinetics methods but it can also be used for end point manual methods. Adaptations to the various automates are available upon request. The assay is performed at the controlled temperature of 37°C and the colour development is measured at 405 nm.

D.750.02/BI/0202

6560 Gove Court · Mason, OH 45040

Phone: 513.770.1991

Toll Free: 866.783.3797

Fax: 513.573.9241

Email: info@aniara.com

www.aniara.com

#### CALIBRATION:

According to the expected hirudin (or DTI) concentration in tested plasma, the "low range" or the "high range" protocol must be used.

**Note:** Calibration of the assay can be performed with normal plasma supplemented with hirudin or with the Plasma Hirudin Standard Low (# ASC020K) and Plasma Hirudin Standard High (# ASC020L) kits (low or high range). When other direct thrombin inhibitors are used, a calibration curve can be prepared by spiking the inhibitor in normal plasma. Alternatively, inhibition can be expressed as "hirudin equivalent".

#### Low hirudin concentration range:

This protocol is proposed for testing hirudin (or DTI) in citrated plasma, and is used for hirudin concentrations of about 1 µg/mL. It must be adapted to the thrombin inhibitor used.

**Calibration curve:** The Plasma Hirudin Standard kit, low range (# ASC020K), can be used for preparing the calibration curve. If a homemade calibration is used, prepare 1 mL of normal citrated plasma (or a normal plasma pool) containing 2 µg/mL of hirudin, using preferably the hirudin used for patient's treatment. Using this preparation, or with the calibration kit, the calibration curve is obtained as follows:

Hirudin	µg/mL	0	0.5	1	1.5	2
Normal plasma with 2 µg/mL hirudin (µL)		0	25	50	75	100
Normal plasma (µL)		100	75	50	25	0
R3 buffer (µL)		900	900	900	900	900

The calibration plasmas, already diluted 1:10, are obtained.

#### High hirudin concentration range:

This protocol is used for hirudin concentrations in plasma of about 3 to 4 µg/mL.

**Calibration curve:** The Plasma Hirudin Standard kit, high range (# ASC020L), can be used for preparing the calibration curve. If a homemade calibration is used, prepare 1 mL of normal citrated plasma (or a normal plasma pool) containing 5 µg/mL of hirudin, using preferably the hirudin used for patient's treatment. Using this preparation or the calibration kit, the calibration curve is obtained as follows:

Hirudin	µg/mL	0	1.25	2.50	3.75	5
Normal plasma with 5 µg/mL hirudin (µL)		0	25	50	75	100
Normal plasma (µL)		100	75	50	25	0
R3 buffer (µL)		2400	2400	2400	2400	2400

The calibration plasmas, already diluted 1:25, are obtained.

In order to get the full assay performances, the calibration curve must be run within two hours.

#### TESTED SAMPLES AND CONTROLS:

According to the protocol used, samples and controls are assayed at the 1:10 dilution (low range protocol, for expected concentrations ranging from 0 to 2µg/ml), or at the 1:25 dilution (high range protocol, for expected concentrations ranging from 0 to 5µg/ml) in R3 Buffer. The diluted samples must be tested within 2 hours.

#### ASSAY PROTOCOL:

##### Manual Method:

Into the microwell or the plastic test tube, incubated at 37°C, introduce:

Reagents	Microplate	Test Tube
Calibrators, or diluted tested plasmas, or Controls	50µl	200µl
Incubate for 1 min at 37°C, then introduce:		
R1 : Substrate	50µl	200µl
Mix and incubate for 1 min at 37°C, then introduce:		
R2 : Thrombin, preincubated at 37°C	50µl	200µl
Mix and incubate for 2 min at 37°C, exactly :		
Stop the reaction by introducing:		
Citric Acid (20g/L), or 20 % Acetic Acid	100µl	400µl
Mix and measure the Absorbance at 405nm against the sample blank.		

The yellow colour obtained is stable for 2 hours.

The sample blank is obtained by mixing the reagents in the reverse order from that of the test i.e.: Citric Acid (20 g/L), thrombin, substrate, diluted plasma.

Measure the Absorbance at 405 nm (A405). Subtract the sample blank from the A405 obtained for the assay.

##### Kinetics mode:

The assay can be read using a kinetics mode. In this case the change in absorbance is recorded from 10 to 100 seconds following the addition of substrate. There is then no need to subtract the sample blank, or to stop the reaction. The results are obtained using the change in absorbance (ΔA405) for calibrators and tested specimen.

##### Automated methods:

Adaptations to the various analysers are available upon request. The assay is then performed kinetically. The reaction does not require to be stopped and sample blanks are automatically subtracted. Reconstitution volumes can vary according to the automate used. Refer to each specific adaptation and specific cautions for each instrument.

#### Note:

- If higher or lower reactive volumes than those indicated here above are required for the method used, the same respective proportions between reagent concentrations and volumes used, must be adhered to, in order to maintain the assay performances.
- Run a sample blank in presence of highly lipemic, icteric or haemolysed plasmas, or if the plasmas has a "colour" different from the usual one.

#### QUALITY CONTROL:

Using commercially available quality control plasmas, titrated for the assayed DTI, allows validating the calibration curve, as well as the homogeneous reactivity from run to run, when using a same lot of reagents. The calibration curve is acceptable when the concentrations measured for controls are within the acceptance range. Various control plasmas are available: Plasma Hirudin Control (# ASC025K). Each laboratory should verify its own target values and acceptance ranges, for each assayed DTI, in the exact working conditions, for each new lot of controls.


**Note:** A new calibration curve must be carried out for each new lot of reagents, after each important maintenance of the analyzer, or when measured values for the quality controls are out of the acceptance range determined for the method. Each laboratory can establish its own acceptance ranges, according to the instruments and protocols used.

#### RESULTS:

- For the end-point method, use a linear graph paper and plot on abscissa the hirudin (or DTI) concentration (µg/ml) and on ordinates the corresponding absorbances (A405). An inverse linear relationship is obtained between hirudin (or DTI) concentrations and measured A405 : the calibration curve is valid if  $r^2 \geq 0.98$ , and if measured values for controls are in compliance for the method.
- The hirudin (or DTI) concentration in the tested sample (diluted 1:10 for the low range or 1:25 for the high range protocol) is directly obtained on the calibration curve. Results are expressed as µg/ml of hirudin (or DTI).
- When the kinetics mode is used, proceed the same way by plotting the ΔA405 values obtained, instead of A405.
- Using automated methods, the concentrations are directly calculated by the analyser, respectively to the calibration curve, and the sample dilution used.
- When the assay dilution is 1:10 (low range) or 1:25 (high range), the hirudin (or DTI) concentration is directly read on the calibration curve.

#### EXAMPLE OF CALIBRATION CURVE:

The calibration curve below is an example only, using the manual method. Only the calibration curve generated for the series of assays performed must be used for calculating the hirudin concentrations.

#### PERFORMANCE CHARACTERISTICS:

- The dynamic range is from 0 to 2µg/ml (low range) or 0 to 5µg/ml (high range) hirudin.
- The detection threshold for the assay is evaluated on the calibration curve by measuring the "apparent" hirudin concentration, which corresponds to the mean A405 value obtained for a sample free of hirudin (normal plasma) less 3 Standard Deviations (SD). This detection threshold is of about 0.15 µg/ml for the low range, and 0.30 µg/ml for the high range.
- Biophen DTI reagent contains an heparin neutralizing substance : the assay can be run on heparinized samples (no significant interference is observed up to 2 IU/ml UFH or LMWH in plasma).
- For an optimal therapeutic efficacy, the measured DTI concentration must be within the therapeutic range defined by the manufacturer for the corresponding indication.
- The assay is calibrated and optimised for hirudin concentration, expressed in µg/mL. The specific activity for hirudin drugs can vary from product to product or from lot to lot (from < 10,000 ATU/mg to > 15,000 ATU/mg). The curves are constructed respectively to the hirudin concentration. If a calibration by hirudin activity, expressed in ATU\*µL, is needed, or when a different thrombin inhibitor is used, the user must take into account the specific anti-thrombin activity of the preparation used. (\*ATU: Anti-Thrombin Unit).
- The influence of progressive thrombin inhibitors can be neglected because of the short incubation time.
- Example of reproducibility results obtained using the STAR instrument:

	µg/ml hirudin	Intra Assay CV% -N=10	Inter Assay CV% (N=4)
Sample 1	1.1	4.8%	<5%
Sample 2	2.1	1.5%	<5%

- **Correlations:** The Biophen DTI assay shows good correlation with Hemoclot Thrombin Inhibitors kit (CK002K), for hirudin assay, performed on STAR instrument.
- **Limitations of the procedure:**
  - Avoid blood activation, during specimen collection and plasma preparation.
  - Any other anti-thrombin substance (than heparin) in the tested specimen is susceptible to affect the results.
  - In order to get the optimal assay performances, the working instructions must be carefully observed.
  - In our experience, this assay is not appropriate for measuring Argatroban® concentrations in plasma. This DTI exhibits a strong anti-thrombin anticoagulant activity, but it has only a weak inhibition for chromogenic thrombin activity, at the therapeutical doses used.

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