CE **BIOPHEN**
NORMAL CONTROL PLASMA
Ref A223201

Normal Control Plasma for coagulation assay quality control

For in vitro diagnostic use only

Manufactured By: **HYPHEN BioMed**

Last revision: 09/06/2011

ENGLISH

INTENDED USE:

BIOPHEN Normal Control Plasma is a set of 12 vials of normal citrated human plasma proposed for the quality control of some coagulation factors assays. The following table shows the various parameters, which are measured using assays from HYPHEN BioMed or other manufacturers, and according to the package inserts:

Parameter	Method
Factor VII + X	Clotting
ATIII	Chromogenic
Protein C	Chromogenic
Protein C	Clotting
Free Protein S (Ag)	Elisa
Protein S	Clotting
Factor VIII:C	Chromogenic
Prothrombin	Clotting / Chromogenic
Factor V	Clotting
Factor VII	Clotting / Chromogenic
Factor IX	Clotting / Chromogenic
Factor X	Clotting / Chromogenic
Factor XI	Clotting
Factor XII	Clotting
Plasminogen	Chromogenic
Fibrinogen	Clotting
Lupus Anticoagulant	DVVtest® / DVVconfirm®
PT/INR and aPTT	Clotting

DVVtest, DVVconfirm are registered trade marks from American Diagnostica Inc. The BIOPHEN Normal Control Plasma is tested for the absence of Lupus Anticoagulant and can be used as a negative control for this investigation.

This control plasma is also tested for the absence of Activated Protein C resistance (Act PC-r). When the APTT is performed with or without Activated Protein C (APC) the ratio obtained (APTT + APC/APTT) is ≥ 2.00 . When tested with Hemoclot Quanti VL kit, the expected value is $<10\%$ FVL.

REAGENTS SUPPLIED:

12 vials of 1 ml of normal citrated human plasma, lyophilised.

Note: It contains an antibiotic as preservative. Each donor unit used for the preparation of control plasmas is a human plasma, which has been tested with registered methods for the presence of Hepatitis B Surface Antigen, Hepatitis C Antibodies (HVC) and antibodies to HIV 1 and 2 and was found negative. However, no test can completely exclude the presence of infectious agents. Any product of human origin, and more especially plasma, must be considered as being potentially infectious and must be handled with all the required cautions for this kind of material.

STORAGE CONDITIONS:

Unopened reagents, must be stored at 2–8 °C, in their original packaging box. They are then stable until the expiration date printed on the label.

Note: The stability studies at 30°C show that the reagents can be shipped at room temperature without damage.

PREPARATION AND STABILITY OF REAGENTS:

Reconstitute each vial with exactly 1 mL of distilled water. Shake thoroughly until complete dissolution of the content (vortex). Incubate at room temperature (18–25°C) for 30 min, while shaking the vial from time to time.

Homogenise the content before each use.

Stability of Normal control plasma, kept in its original vial:

- For AT, PC, VII-X, Plasminogen, FII, FVII, FIX, FX, FXI, FXII, Fibrinogen and aPCr (%FVL): **8 hours** at room temperature (18–25°C), or **24 hours at 2–8°C**.
- For FVIII:C, FV, PS: **4 hours** at room temperature (18–25°C), or **8 hours at 2–8°C**. Do not freeze.

Cautions: In order to improve stability, reagents must be closed with their original screw caps following each use. Reagents must be handled with care, in order to avoid any contamination during use. It is recommended to homogenize each vial before use, in order to have a good reproducibility, all the time.

APPLICATIONS:

It can be used as a normal control plasma for some coagulation assays. For each lot, the concentrations and acceptance ranges are reported for each parameter tested on the flyer provided within the kit. The various parameters reported are accurately determined by reference to the corresponding NIBSC (National Institute

for Biological Standards and Controls, UK) International Standards, when available, or against an internal reference standard. The following table shows the usual ranges expected for the BIOPHEN Normal Control Plasma.

Parameter	Method	Acceptance range
Factor VII + X	Clotting	> 75%
ATIII	Chromogenic	> 80%
Protein C	Chromogenic	> 80%
Protein C	Clotting	>60%
Free Protein S (Ag)	Elisa	>60%
Protein S	Clotting	>50%
Factor VIII:C	Chromogenic	> 60%
Lupus Anticoagulant	DVVtest/DVVconfirm	≤ 1.2
Prothrombin	Clotting / Chromogenic	> 75%
Factor V	Clotting	> 75%
Factor VII	Clotting / Chromogenic	> 75%
Factor IX	Clotting / Chromogenic	> 60%
Factor X	Clotting / Chromogenic	> 75%
Factor XI	Clotting	> 70%
Factor XII	Clotting	> 70%
Plasminogen	Chromogenic	> 75%
Fibrinogen	Clotting	>2 g/L
Act PC-r	Ratio APTT \pm APC	≥ 2.00
ActPCr (% FVL)	Hemoclot Quanti VL	<10%
PT (INR)	Clotting	< 16 sec (INR <1.30)
aPTT	Clotting	37 \pm 5 sec

PERFORMANCE CHARACTERISTICS:

The following obtained values are provided as an example only.

Parameter	[C]	Acceptance range
Factor VII + X	100%	90-110
ATIII	97%	87-107
Protein C (chromogenic)	87%	77-97
Protéine C (coagulante)	82%	66-98
Protéine S libre	75%	64-86
Protéine S (coagulante)	71%	57-85
Factor VIII:C	73%	63-83
DVVtest/DVVconfirm ratio	1.05	-
Prothrombin	91%	81-101
Factor V	87%	77-97
Factor VII	102%	92-112
Factor IX	88%	78-98
Factor X	88%	78-98
Facteur XI	72%	61-83
Facteur XII	79%	67-91
Plasminogen	96%	86-106
Fibrinogen	2.67	2.17 - 3.17
Act PC-r (ratio)	2.56	-
PT	13.6	11.6-15.6
INR	1.06	0.96-1.16
aPTT	34	29-39

For each parameter, the concentration and acceptance range may present variations from lot to lot, but it is exactly measured for each lot and reported on the flyer provided within the kit. When the BIOPHEN Normal Control Plasma is used as quality control plasma for the assay of coagulation factors, the values obtained must be within the acceptance ranges reported for the lot used, in order to validate the test series. Should the value be out of these ranges, the results for the corresponding series must be considered as invalid. It is then recommended to rerun the series and to check all the assay parameters.

Note: The results obtained for VII-X assay can slightly vary according to the thromboplastin reagent type and lot, and the instrument used. The target value and acceptance range for VII-X parameter must consequently be confirmed and adjusted, if necessary, for each new lot of control, in the laboratory working conditions. PT and aPTT results must be checked and adjusted if necessary in the same way.

CAUTIONS:

- Like all lyophilised plasmas, the plasmas from the BIOPHEN Normal Control Plasma are more or less cloudy after reconstitution. This is due essentially to the lipids that, after lyophilisation, become less soluble and can form a light deposit.
- If necessary, let each vial 10 minutes at room temperature and shake gently before use in order to homogenise the content.

Reagents must be handled with care, in order to avoid any contamination or activation during use. Any plasma containing a coagulum or contamination must be rejected.

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