Manufactured By: HYPHEN BioMed

BIOPHEN CS-11(22) Factor Xa Chromogenic Substrate Vial of 25 mg / flacon de 25 mg

Ref. A229015

STORE AT
2-8 °C

For research use only

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CONSERVER A
2-8 °C

Origin:

Synthetic generic Chromogenic substrate, highly purified and stabilized.

Specificity:

Recommended substrate for Factor Xa.

Presentation:

Vial containing 25 mg of BIOPHEN CS-11(22), lyophilised in presence of Mannitol as a bulking agent.

Reconstitution:

According to the research protocol used, the BIOPHEN CS-11(22) chromogenic substrate can be restored with variable volumes of distilled water ; for example 10 mL can be used for a substrate concentration of 2.5 mg/mL, or 20 mL for a substrate concentration of 1.25 mg/mL. Shake thoroughly until complete dissolution (vortex). Let to stabilize for 30 min. at room temperature.

Peptide sequence:

Mixture (50%-50%) of Bz-Ile-Glu (γOCH₃)-Gly-Arg-pNa (form 1) and Bz-Ile-Glu (γOH)Gly-Arg-pNa (form 2)

Purity grade:

> 95%.

Molecular weight:

711.8 (form 1) and 697.7 (form 2).

Free pNA content

< 0.05%.

Use:

For in vitro use only.
All research studies and protocols where a source of chromogenic substrate for Factor Xa is required.

Respective reactivities:

FXa	Thrombin	Plasmin	Kallicrein	aPC
100	1	2	3	0

Assay conditions must be duly established for rendering the assay conditions totally specific for Factor Xa, when this substrate is used.

Stability:

- Restored :**
- 7 days at room temperature
 - 3 months at 2-8 °C
 - **Do not freeze.**

Lyophilized : Until the expiration date printed on the vial.

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