Manufactured By: HYPHEN BioMed

# Bentonite suspension

## Suspension de Bentonite

(Vial of 50 ml / Flacon de 50 ml)

Ref. AAR024A

STORE AT  
2-8 °C

For research use only

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CONSERVER A  
2-8 °C

Lot:

Expiration:

#### Composition:

Buffered bentonite suspension, ready to use; contains 0.9 g/L Sodium Azide as preservative.

**\*CAUTION:** Sodium azide (NaN<sub>3</sub>), may react with lead and copper plumbing to form highly explosive metal azides. In order to avoid this risk, flush with large volumes of water when discarding into a sink.

Ready to use.

#### Presentation:

50 ml vial, ready to use.

#### Use:

Removal of Fibrinogen from plasma for testing FPA. Mix thoroughly before use. When opened, it can be used for 4 weeks, stored at 2-8 °C, and provided that any contamination is avoided.

Ready to use

For in vitro use only.

#### Treatment of plasma:

Before measurement of FPA, cross reactive Fibrinogen must be removed from plasma by bentonite adsorption. Mix thoroughly the bentonite suspension in order to make it homogeneous. To 1 ml of the anticoagulated plasma (special anticoagulant for FPA testing), add 0.5 ml of bentonite suspension. Mix and agitate for 10 min. using an end-over-end agitator. Centrifuge for 20 min. at 2,500 g and collect 1 ml supernatant.

Proceed to a new bentonite adsorption by adding again 0.5 ml of suspension to the 1 ml of supernatant, in a similar manner. The bentonite treated plasma is then fibrinogen free. It must be used within:

- 24 hours at room temperature or at 2-8 °C.
- 1 month when stored frozen at -20 °C or below.

The bentonite treated plasma is two-fold diluted, and the measured FPA concentration must be multiplied by 2.

#### Stability:

In its original vial, unopened and kept at 2-8 °C, until the expiration date printed on the label.

When open and protected from any contamination, this solution is stable for 4 weeks at 2-8 °C.

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