Manufactured By: HYPHEN BioMed

Bovine Thrombin

Thrombine bovine

1000 NIH/ Vial / 1000 NIH/flacon

Ref. ABE102C

STORE AT
2-8°C

For research use only

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CONSERVER A
2-8°C

Origin:

Prepared from a prothrombin complex concentrate extracted from bovine citrated plasma, and fully activated. Thrombin is then purified using an ion exchange chromatography and stabilised. It is a highly purified preparation, mainly in the α form, with a specific activity of more than 2,000 NIH/mg.

Presentation:

Thrombin activity is reported in harmonized NIH or IU or USP units by reference to the WHO/NIBSC International Standard for Human Thrombin (01/580) (refer to "A reunification of the US ("NIH") and International Unit into a single standard for Thrombin", Longstaff et al, Thromb Haemost 2005, 93:261-6)). NIH is a clotting unit for thrombin activity. Stabilised and lyophilised bovine Thrombin. Vial containing approximately 1000 NIH of thrombin.

Reconstitution:

Each vial must be restored with 2 ml distilled water, in order to get a concentration of about 500 NIH/ml. Dilute in the adequate buffer, if necessary.

Excipients: Content per vial: Poly-Ethylene-Glycol 6000 (PEG 6000) (40mg), Sodium Chloride (12mg), Glycine (40mg), stabilizers.

Purity: Bovine Thrombin has one major band of 35,000 daltons on SDS-PAGE.

Activity:

- Clotting time after reconstitution with 2ml and dilution at 1:75:
 - on purified bovine fibrinogen (5 mg/ml) : **8.5" ± 2.5"**
 - on purified human fibrinogen (5 mg/ml) : **11.5" ± 3"**
 - Chromogenic activity on CS-01(38) thrombin substrate: **A405 : > 1.20/min.NIH.ml**
 - Clotting activity determined in NIH (or IU).
 - Chromogenic activity on specific substrates (expressed in nkats):
- When tested in the optimized conditions (0.05M Tris buffer at pH8.40, containing 0.30M NaCl), the following specific activities are observed on the various IIa substrates (expressed in **nkats/µg**):

CS-01(38)	CS-01(81)
4,4 ± 1,0	4,6 ± 1,0

The exact activity in nkats/µg with CS-01(38) is reported on the analysis certificate for each lot.

Viral safety:

The bovine plasma used for Bovine Thrombin preparation was tested with registered methods and found negative for bovine infectious diseases, notably for the bovine spongiform encephalopathy. However, no assay may warrant the total absence of infectious agents. Any product of bovine origin must then be handled with all the required cautions, as being potentially infectious.

Use: For *in vitro* use only. All research studies and protocols where a source of highly purified bovine Thrombin is necessary.

Stability:

- Restored :** - 2 days at room temperature (18-25 °C).
- 7 days at 2-8 °C.
- Lyophilized :** Until the expiration date printed on the vial.

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