Manufactured By: HYPHEN BioMed

Bovine Thrombin
Thrombine bovine
500 NIH/Vial or Set of 10x500 NIH
500 NIH/flacon ou coffret de 10x500 NIH

Ref. ABE 501C / ABE501K

STORE AT
2-8 °C

For research use only

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CONSERVER A
2-8 °C

Origin:

Prepared from a prothrombin complex concentrate extracted from bovine citrated plasma, and fully activated. Thrombin is then purified using an ion exchange chromatography and stabilised. It is a highly purified preparation, mainly in the α form, with a specific activity of more than 2,000 NIH/mg.

Presentation:

Thrombin activity is reported in harmonized NIH or IU or USP units by reference to the WHO/NIBSC International Standard for Human Thrombin (01/580) (refer to "A reunification of the US ("NIH") and International Unit into a single standard for Thrombin", Longstaff et al, Thromb Haemost 2005, 93:261-6)). NIH is a clotting unit for thrombin activity. Stabilised and lyophilised bovine Thrombin.

Vial containing approximately 500 NIH of thrombin.

ABE501C: 500 NIH/vial; ABE501K: Set of 10 vials containing 500 NIH/vial.

Reconstitution:

Each vial must be restored with 5 ml of distilled water, in order to get a concentration of about 100 NIH/ml. Dilute in the adequate buffer, if necessary.

Excipients: Content per vial: Poly-Ethylene-Glycol 6000 (PEG 6000) (50mg), Sodium Chloride (30 mg), Calcium gluconate (110mg), Glycine (50mg), polybren, BSA (10mg), Stabilizers.

Purity: Bovine Thrombin has one major band of 35,000 daltons on SDS-PAGE.

Characteristics:

1. Clotting time on normal plasmas tested 1:20 with Owren Koller buffer (100 μ L diluted plasma and 100 μ L bovine Thrombin):

6 to 20 Sec.

2. Chromogenic activity on thrombin substrate:

A405 : > 1.20/min.NIH.ml

3. Clotting activity determined in NIH (or IU).
4. Chromogenic activity on specific substrates (expressed in nkats):

When tested in the optimized conditions (0.05M Tris buffer at pH8.40, containing 0.30M NaCl), the following specific activities are observed on the various IIa substrates (expressed in **nkats/ μ g**):

CS-01(38)	CS-01(81)
4.4 \pm 1.0	4.6 \pm 1.0

The exact activity in nkats/ μ g with CS-01(38) is reported on the analysis certificate for each lot.

Viral safety:

The bovine plasma used for Bovine Thrombin preparation was tested with registered methods and found negative for bovine infectious diseases, notably for the bovine spongiform encephalopathy. However, no assay may warrant the total absence of infectious agents. Any product of bovine origin must then be handled with all the required cautions, as being potentially infectious.

Use:

For in vitro use only. All research studies and protocols where a source of highly purified bovine Thrombin is necessary, such as fibrinogen testing with the Clauss method.

Stability:

- Restored: - 2 days at room temperature (18-25 °C)
- 7 days at 2-8 °C.

Lyophilized: Until the expiration date printed on the vial.

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