Manufactured by Hyphen BioMed.

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**CONSERVER A**  
**2-8 °C**

**Origin:**

Highly purified from citrated bovine plasma

**Presentation:**

**ABP 102A:** 1 vial containing approximately 100 µg of highly purified bovine Factor X.

**ABP102K** Set of 6 vials containing approximately 100 µg of highly purified bovine Factor X.

**Reconstitution:**

Each vial must be restored with 1 ml distilled water. The solution can be diluted to the desired concentration when required.

**Excipients:**

- Glycine and Hepses
- 9 mg sodium chloride/vial

**Purity:**

One major band of 55,000 daltons on SDS-PAGE.

**Activity:**

Tested for its clotting activity (factor X deficient plasma) and its amidolytic activity, following activation with RVV.

Activity of about 80 U/mg (1 Unit "U" is the amount in 1 ml normal bovine plasma).

**Viral safety:**

The plasma used for Factor X purification was tested with registered methods and found negative for bovine infectious diseases, notably for the bovine spongiform encephalopathy. However, no assay may warrant the total absence of infectious agents. Any product of bovine origin must then be handled with all the required cautions, as being potentially infectious.

**Use:**

For in vitro use only.  
All research studies and protocols where a source of highly purified bovine Factor X is necessary.

**Stability:**

- Restored :**
- 24 hours at room temperature (18-25 °C)
  - 72 hours at 2-8 °C
  - 6 months at -20 °C or below.

**Lyophilized :** Until the expiration date printed on the vial.

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