HEMOCLOT Protein S Ref ACK041K

Measurement of Protein S activity using a clotting method

Manufactured By: HYPHEN BioMed

For *in vitro* diagnostic use only

Last revision: 09/05/2011

INTENDED USE:

The HEMOCLOT Protein S kit is an *in vitro* clotting assay for the quantitative determination of Protein S (PS) in human citrated plasma, with a manual or automated method.


SUMMARY:

In presence of Activated Protein C (APC, for which Protein S is a cofactor), phospholipids and Calcium, Protein S inactivates Factor Va and Factor VIIIa (the activated forms of Factors V and VIII). In physiological conditions, Thrombin binds to Thrombomodulin exposed on endothelial cells and which is abundant in capillaries, and loses its procoagulant activity, whilst it acquires an anticoagulant function capable of activating Protein C to APC. APC then combines with the Free form of Protein S and this complex binds onto phospholipids surface (i.e. activated platelets) in the presence of calcium. Clotting pathways are then stopped by inactivation of Factors Va and VIIIa. This PC/PS inhibition pathway has a decreased activity if the free form of Protein S is deficient or abnormal. Congenital or acquired Protein S deficiencies are a risk factor for venous thrombosis.

ASSAY PRINCIPLE:

The HEMOCLOT Protein S method is an APTT like clotting assay, but triggered by Factor IXa in presence of phospholipids, Calcium, and of a constant an in excess amount of Activated Protein C (APC).

In the first step, the diluted assayed plasma is mixed with Protein S deficient plasma (R1). Then, the Activator Reagent (R2), in a constant and optimised concentration, is added. Clotting is initiated by the addition of Calcium (Ca²⁺). Clotting time is then recorded. Protein S being the limiting factor, there is a direct relationship, between the Protein S concentration and the corresponding clotting time.

SPECIMEN:

Citrated human plasma.

REAGENTS:

HEMOCLOT Protein S contains the reagents required for 3 series of 20 tests with manual or automated methods.

R1: Reagent 1: Protein S Deficient Plasma.

Protein S Deficient plasma, immuno-depleted, lyophilized in the presence of a heparin neutralizing agent: 3 vials (to be reconstituted with 1 ml of distilled water).

R2: Reagent 2: Activator reagent

Lyophilized material containing human Factor IXa, human APC, and phospholipids, in an optimized concentration. 3 vials (to be reconstituted with 1 ml of distilled water).

Precaution and warnings:

- Human plasma used for the preparation of the reagents was tested and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.
- All the required cautions must be respected in order to avoid any risk of ingestion or accidental introduction of reagents in body. In case of skin contact, wash extensively with water. In case of contact with a wound, address to the appropriate medical service, and indicate the biological origin and the nature of the product.
- Avoid contact with skin and eyes. Do not empty into drains. Wear suitable protective clothing.

REAGENTS AND MATERIAL REQUIRED BUT NOT PROVIDED:

Reagents:

- Distilled water, preferentially sterile.
- Calcium chloride 0.025M (eg: # AAR001A/AAR001K).
- Imidazole buffer (eg: # AAR021A/AAR021K/AAR021L); The same buffer must be used for all tests performed.
- Plasma Calibrator titrated for PS (eg: BIOPHEN Plasma Calibrator #A222101) or international or internal reference material, or reference frozen citrated normal plasma pool.
- Normal and Abnormal Quality Control Plasmas titrated for PS (eg: BIOPHEN Normal Control Plasma # A223201, and BIOPHEN Abnormal Control Plasma #A223301).

Material:

- Electromagnetic water bath or semi automatic or automatic clotting instruments.
- Calibrated pipettes and stop watch for manual methods.

STORAGE CONDITIONS:

The reagents must be stored at 2-8°C, in their original packaging box. They are then stable until the expiration date printed on the box.

Note: Stability studies for 3 weeks at 30°C show that the reagents can be shipped at room temperature for a short period without damage.

PREPARATION AND STABILITY OF REAGENTS:

Note: Refer to each specific instrument adaptation.

Preparation:

R1: Reagent 1: PS Deficient Plasma

Reconstitute each vial with **exactly 1 ml** of distilled water. Shake thoroughly until complete dissolution of the content (vortex). Incubate at room temperature (18-25°C) for 15 min, while shaking the vial from time to time. Homogenize the content before each use.

R2: Reagent 2: Activator reagent

Reconstitute each vial with **exactly 1 ml** of distilled water. Shake thoroughly until complete dissolution of the content (vortex). Incubate at room temperature (18-25°C) for 15 min, while shaking the vial from time to time. Homogenize the content before each use.

Stability:

Stability of restored reagents R1 and R2, provided that any contamination or evaporation is avoided, kept in their original vial:

- 24 hours at 2-8°C.
- 8 hours at room temperature (18-25 °C).
- 1 month at -20°C or below (before use, thaw in a water bath at 37°C until complete thawing and homogenize).

Cautions:

- In order to improve stability, reagents must be closed with their original screw cap following each use.
- Reagents must be handled with care, in order to avoid any contamination during use.
- Incubating the reconstituted vials at RT allows stabilizing the reagents, and obtaining a homogeneous reactivity.
- Take care to limit as much as possible any evaporation of the reagents during use, eg. by using chimneys.
- R1, R2 vials are closed under vacuum. Remove carefully the stopper, in order to avoid any lost of powder when opening the vials.
- For any automated method used, in any case, the established reactive ratios (respective reagent concentrations in the reactive milieu) between R1, R2 must be adhered to.
- Use only reagents from kits with the same lot number. Do not mix reagents from kits with different lots when running the assay. **Reagents R1, R2 are optimized for each lot of kits.**

PREPARATION OF PLASMA (SPECIMEN COLLECTION, HANDLING AND STORAGE):

Blood (9 volumes) must be collected on 0.109 M citrate anticoagulant (1 volume), with great care, in a silicon glass or a plastic tube. Sampling must be performed through a net venipuncture, avoiding any blood activation. Within 2 hours, blood must be centrifuged at 2,500 g for 20 min, and plasma decanted into a plastic tube, using a plastic pipette.

Storage of plasma:

- Up to 4 hours at Room Temperature (18-25°C).
- Up to 1 month frozen at -20°C or below (before use, thaw for 15 min. in a water bath at 37°C).

Refer to GEHT or NCCLS/CLSI recommendations for further instructions on specimen collection, handling and storage. Discard any sample presenting an unusual aspect.

TEST PROCEDURE:

HEMOCLOT Protein S kit is a clotting method, manual or automated. Adaptations to automates are available upon request. The assay is performed at 37°C, and the clotting time, triggered by the addition of calcium, is measured.

CALIBRATION:

- Calibration is performed with a normal pooled citrated plasma (made with plasmas from at least 30 normal individuals, males or females, aged between 18 and 55 years, and free of any medication or disease), with the assigned value of **100 % PS**. The assay includes a standard plasma dilution of **1:10**. By definition, this latter dilution of the pool represents the **100 % PS** activity. The dynamic range is from **0 to 100 % PS**. The **100 % PS** activity is then the **1:10** dilution of the **plasma pool** in imidazole buffer.
- Or the kit can also be calibrated with a commercially available plasma calibrator (eg **BIOPHEN Plasma Calibrator #A222101**), with a known PS concentration (C): the **1:10** dilution corresponds to the indicated "C" concentration. The **100%** concentration is obtained (in the assay conditions) by using the following dilution factor: **10 x C :100**.

Prepare 3 ml of the **1:10** dilution of the normal plasma pool, or of a (**10 x C :100**) dilution of the PS titrated calibrator, in imidazole buffer (eg AR021). This corresponds to **100%** PS. The calibration curve can then be obtained by preparing serial dilutions in buffer, as follows:

PS concentration	0%	25%	50%	75%	100%
Dilution	1:1	1:40	1:20	3:40	1:10
Plasma Pool 1:10	0 mL	0.250 mL	0.500 mL	0.750mL	1 mL
Imidazole buffer	1 mL	0.750 mL	0.500 mL	0.250mL	0 mL

In order to get the full assay performances, the calibration curve must be prepared just before running the assay.

D.750.02/CK/041K

7768 Service Center Drive • West Chester OH 45069
Phone: 513.770.1991 Toll Free: 866.783.3797

Fax: 513.573.9241 Email: info@aniara.com

www.aniara.com

TEST PLASMA:

Dilute the tested samples (and controls) 1:10 using the imidazole buffer.

For expected concentrations higher than 100%, exact values can be obtained by assaying the plasma at the 1:20 dilution, and results need then to be multiplied by 2.

Caution: to ensure optimal performances of the assay, perform all tests (calibration, samples, controls) extemporaneously and successively without interruption, and the same buffer.

ASSAY PROTOCOL:

Manual Method:

Into the **plastic** test tube, incubated at 37°C, introduce:

Assay	Manual Method
Calibrators; or 1:10 diluted plasma or control	50µl
R1 reagent preincubated at 37°C	50µl
Mix and incubate for 1 minute at 37°C	
R2 reagent preincubated at 37°C	50µl
Mix and incubate for 3 minutes at 37°C	
CaCl2 0.025M (preincubated at 37°C, and stirred)	100µl
Record Clotting Times, in seconds	CT

Automated methods:

Adaptations to various analysers are available upon request. Refer to each specific adaptation and specific cautions for each instrument.

Note: Isotonic Imidazole buffer (ref AR021) is recommended for diluting calibrators, controls and assayed plasmas for the measurement of Protein S with Hemoclot Protein S, as it corresponds to the easiest and most universal working conditions. However, the Owren Koller type buffer (ref AR003) can also be used for diluting plasmas. The same type of buffer must be used throughout the assay, for diluting calibrators, controls and plasmas.

QUALITY CONTROL:

Using commercially available quality control plasmas, titrated for PS, allows validating the calibration curve, as well as the homogeneous reactivity from run to run, when using a same lot of reagents. The calibration curve is acceptable when the concentrations measured for controls are within the acceptance range. Various control plasmas are available: **BIOPHEN Normal Control Plasma (#A223201)** and **BIOPHEN Abnormal Control Plasma (#A223301)**. Each laboratory should verify its own target value and acceptance range, in the exact working conditions, for each new lot of controls.


Note: A new calibration curve must be carried out for each new lot of reagents, after each important maintenance of the analyzer, or when measured values for the quality controls are out of the acceptance range determined for the method. - Each laboratory can establish its own acceptance ranges, according to the instruments and protocols used. - Include at least one quality control (at different levels) in each test series.

RESULTS:

- For the manual method, use a linear graph paper and plot on abscissae the PS concentration (%) and on ordinates the clotting times (sec). Using a linear graph paper, the assay is linear in the range 10- 100% PS.
- Calculate the "r²" value. Calibration is acceptable if: r² ≥ 0.98, and if measured values for controls are in compliance.
- The Protein S concentration in the tested sample (when assayed at the 1:10 dilution) is directly obtained from the calibration curve. Results are expressed as %, or U/ml (a normal plasma pool contains by definition 1 U/ml of Protein S activity, which in the assay yields a concentration of 100 %).
- Using automated methods, the PS concentrations are directly calculated by the analyser, respectively to the calibration curve, and the sample dilution used.
- When the assay dilution is 1:10, the PS concentration is directly read on the calibration curve. When predilutions are used, multiply the measured PS concentration by the predilution factor in order to get the concentration in the tested specimen (eg x2 for a sample tested at the 1:20 dilution).

EXAMPLE OF CALIBRATION CURVE:

The calibration curve here below, obtained using KC10 instrument, is indicated as an example only. Only the calibration curve generated for the series of measures performed must be used.

PERFORMANCES AND CHARACTERISTICS:

- Dynamic range: 10-100% PS (or depending on instrument adaptation).
- The detection threshold (measured by the "apparent" Protein S concentration, corresponding to the mean CT obtained for a PS deficient sample plus 3 standard deviations (SD)) is of about 10%.
- Specificity: PS deficient plasma is measured < 5%.
- Reproducibility: Inter assay CV on obtained clotting times are <5% (KC10, STAR).

Limitations of the procedure:

- Blood activation, during specimen collection and plasma preparation, may lead to wrong results.
- The assay can be performed in patients treated with heparin (up to 1 IU/ml), or under dicumarol therapy (PS activity is then decreased). In presence of an abnormally prolonged CT, confirming the diagnosis with another method and/or another sample collected is recommended.
- Special caution should be taken for patients with known high levels of FVIII:C, or with LA or mutated FV (R506Q mutation, FVLeiden), and results confirmed with another method (eg immunological method for Free Protein S) and considering the clinical context.
- Aprotinin inhibits Activated Protein C. The "apparent" Protein S activity could be decreased in patients treated with aprotinin. Using a different method for testing protein S (i.e. immunoassay) is then recommended.
- For a better accuracy, samples measured ≤10% can be tested at the 1:5 dilution, and obtained results divided by 2; for samples expected >100%, the 1:20 dilution can be used and obtained results multiplied by 2. If a different dilution factor from the standard 1:10 is used, the concentration must be corrected by the complementary dilution factor, i.e. x2 for 1:20, or x0.5 for 1:5).
- To ensure optimal performances of the assay, perform all assays (calibration, samples, controls) extemporaneously and successively without interruption. Working instructions must be carefully observed.
- For a same reagent lot and a same tested specimen, the Clotting Time may present variations according to the instrument used (especially between mechanical and optical clot detection), and to the clot detection sensitivity adjustment. Performances, as well as the normal range, and target values and acceptance ranges for each new lot of quality controls used, must then be confirmed (and adjusted if necessary) in the exact laboratory working conditions.

USUAL VALUES:

The 100 % Protein S concentration (or 1 Unit/ml) is the amount in a normal human citrated plasma pool, obtained by pooling plasmas from healthy males or females aged from 18 to 55 years, and out of any medication or disease. The Protein S concentration in healthy adults is usually between 60 and 140% (variable according to gender, slightly lower in females). The Protein S concentration is decreased in neonates (Premature: range 14-38%; New born: range 12-60%; 5 days old: 36%; 30 days old: 60% and up to 90 days old: 90%).

PATHOLOGICAL VARIATIONS:

Protein S activity ≤ 60% indicates the presence of a deficiency, which must be confirmed by another measurement, or on another sample collected from the patient. A deficiency of Protein S activity may be congenital or acquired. It can also develop during pregnancy or in some diseases (transitory), including disseminated intravascular coagulation (DIC) and HIV infection. At the beginning of an inflammation process, Protein S activity can also be decreased. PS concentration is decreased during dicumarol treatment.

BIOCHEMISTRY:

Protein S is a vitamin K dependent human protein, synthesized in the liver, of about 69 kDa. PS is usually present in plasma at a concentration of approximately 25 µg/ml. Plasma Protein S has two presentations: about 60% (15 µg/ml) is bound to the complement component C4b β-chain while the remaining 40% is free (10 µg/ml). Only the free Protein S form presents anticoagulant activity as Activated Protein C (APC) cofactor. A factor V mutant (Factor V Leiden, carrying a single amino acid substitution on position 506 (R506Q)) is resistant to the inactivation function of the PS-APC complex.

CLINICAL INFORMATION:

Congenital deficiency:

There are three types of hereditary protein S deficiencies:

- Type I - decreased PS activity: decreased total and free PS concentrations (quantitative)
- Type II - decreased PS activity: normal free and total PS levels (qualitative).
- Type III - decreased PS activity: decreased free PS and normal total PS levels.

Acquired deficiency:

PS deficiency can also be acquired during:
- vitamin K deficiency or treatment with dicumarol/warfarin. Vitamin K dependent coagulation factors (II, VII, IX, X) are also decreased, and this generates bleeding tendency rather than thrombosis.
- systemic sex hormone therapy and pregnancy, liver disease, and certain chronic infections (eg AIDS).
PS deficiency is the underlying cause in a small proportion of cases with disseminated intravascular coagulation (DIC), deep vein thrombosis (DVT) and pulmonary embolism (PE).
Rare deficiencies due to an acquired or transitory auto-antibody to Protein S have been reported (eg during chicken pox disease in some children).

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