Manufactured By: HYPHEN BioMed

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**CONSERVER A**  
**2-8 °C**

**Origin:**

Citrated human plasma, deficient for Protein S, Immuno-depleted.

**Presentation:**

1 mL vial, lyophilized.

**Characteristics:**

This plasma is deficient for Protein S (<1%), whereas all the other coagulation factors are within about the normal range (> 50%).

**Excipients:**

- Glycine, Stabilizers.

**Viral safety:**

Plasmas used for Protein S Deficient Plasma preparation were tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

**Use:**

For *in vitro* use only.  
Lyophilized, human citrated plasma, deficient for Protein S, for any protocol or research study where a source of human Protein S deficient plasma is required.

**Reagent Preparation:**

- Restore the vial with 1mL distilled water
- Leave for 15 min at room temperature
- Mix gently until complete dissolution of the content.

When opened and protected from any contamination, the reconstituted plasma is stable for:

- 8 hours at room temperature (18-25 °C)
- 24 hours at 2-8 °C
- 2 months, frozen in its original vial, or in a plastic tube (before use, thaw in a water bath at 37 °C, for at least 15 min).

**Reagent stability:**

In its original vial, and before any use, when stored at 2-8 °C, the reagent is stable until the expiration date printed on the label.

- Restored:**
- 8 hours at room temperature
  - 24 hours at 2-8 °C
  - 2 months frozen at -20 °C or below

**Lyophilized:** Until the expiration date printed on the vial.

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