vWF Deficient Plasma

ADP150A / #ADP150K

Deficient plasma for the assay of vWF

Manufactured By: HYPHEN BioMed

For in vitro research use only

Recherche in vitro exclusivement

Last revision / Dernière version: 07/04/2008

ENGLISH

Intended use:

For in vitro use only.

Lyophilized, human citrated plasma, deficient for vWF, for any protocol or research study where a source of human vWF deficient plasma is required.

Reagents:

1 vial of 1 ml (#ADP150A) or 6 vials of 1 ml (#ADP150K) of citrated human plasma, deficient for vWF, immuno-depleted, lyophilized in the presence of glycine and stabilizers.

This plasma is deficient for vWF (<1%), whereas all the other coagulation factors are within about the normal range (> 50%).

Reagent preparation and stability:

In the original package, and before any use, when stored at 2-8°C, the reagent is stable until the expiration date printed on the kit.

• Reagent Preparation:

Restore the vial with 1 ml of distilled water; mix gently until complete dissolution of the content (vortex), let for 15 min. at room temperature(18-25°C); homogenize before each use.

• Reagent stability following reconstitution:

When opened and protected from any contamination, the reconstituted plasma is stable for:

- 8 hours at room temperature (18-25°C)
- 24 hours at 2-8°C
- 2 months, frozen at -20°C or below, in its original vial, or in a plastic tube (before use, thaw in a water bath at 37°C, for at least 15 min).

Note: Plasmas used for the Deficient Plasma preparation were tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

Note: Stability studies at 30°C show that the reagent can be shipped at room temperature without damage.

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