## tPA/PAI-1 Deficient Plasma

# ADP210A / # ADP210K

Deficient plasma for the assay of tPA or PAI-1

Manufactured By: HYPHEN BioMed

For in vitro research use only

Recherche in vitro exclusivement

Last revision / Dernière version: 09/03/2009

### ENGLISH

#### Intended use:

For in vitro use only.

Lyophilized, human citrated plasma, tPA and PAI-1 deficient, proposed for any protocol or research study where a source of human tPA/PAI-1 deficient plasma is required.

#### Reagents:

1 vial of 1 ml (#ADP210A) or 6 vials of 1 ml (#ADP210K) of citrated human plasma, tPA and PAI-1 deficient, immuno-depleted, lyophilized in the presence of glycine and stabilizers.

This plasma is deficient for tPA and PAI-1 (<1ng/ml tPA and PAI-1 antigen, <1ng/ml PAI-1 activity), whereas all the coagulation factors are within the normal range (> 50%).

#### Reagent preparation and stability:

In the original package, and before any use, when stored at 2-8°C, the reagent is stable until the expiration date printed on the kit.

##### • Reagent Preparation:

Restore the vial with 1 ml of distilled water; mix gently until complete dissolution of the content (vortex), let for 15 min. at room temperature (18-25°C); homogenize before each use.

##### • Reagent stability following reconstitution:

When opened and protected from any contamination, the reconstituted plasma is stable for:

- 8 hours at room temperature (18-25°C)
- 24 hours at 2-8°C
- 2 months, frozen at -20°C or below, in its original vial, or in a plastic tube (before use, thaw in a water bath at 37°C, for at least 15 min).

**Note:** Plasmas used for the Deficient Plasma preparation were tested with registered methods and found negative for HIV antibodies, HBs Ag and HVC antibodies. However, no assay may warrant the total absence of infectious agents. Any product of human origin must then be handled with all the required cautions, as being potentially infectious.

**Note:** Stability studies at 30°C show that the reagent can be shipped at room temperature without damage.

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